



363-364-365 Aluminum Folding Screen Parts List




PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1		216600	363 SCREEN FRAME ASSEMBLED	1
1		216610	364 SCREEN FRAME ASSEMBLED	1
1		216620	365 SCREEN FRAME ASSEMBLED	1
2		211500	LEG SET 363 COMPLETE W/HARDWARE	1
2		211510	LEG SET 364 COMPLETE W/HARDWARE	1
2		211520	LEG SET 365 COMPLETE W/HARDWARE	1
2A	A,C	216630	MID LEG (363,364,365)	SEE NOTES
2B	C	566512	SCREW-MACH 10-24 X 1 PAN PHILLIP M/S ZINC	1 per 2 legs
2B	C	566203	NUT-10-24 HEX NYLON LOCK NUT ZINC	1 per 2 legs
2C	B,C	216635	WING LEG (363,364,365)	SEE NOTES
2D	C	566516	SCREW-MACH 10-24 X 2 PAN PHILLIP M/S ZINC	1 per 2 legs
2D	C	566203	NUT-10-24 HEX NYLON LOCK NUT ZINC	1 per 2 legs
2E	C	700800	PLASTIC ENDCAP 3/4" WHITE	1 per leg

NOTES

Note	Description
A	Qty (2) for 363 model & Qty (4) for 364-365 models
B	Qty (4) for 363-364 models & Qty (6) for 365 model
C	Included in the above listed "LEG SETS"

363-364-365 Aluminum Folding Screen Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1A	A	299900	363,364,365 VINYL PANEL STANDARD WHITE	SEE NOTES
1A	A,B	2999XX	363,364,365 SURE CHEK PANEL	SEE NOTES
1B	C	216625	SCREEN ROD (363, 364, 365)	SEE NOTES
1C	C	704500	ROD & ROD SPRING RETAINER(SET)	SEE NOTES
NOTES				
Note	Description			
A	363= (3) Panels, 364= (4) Panels, 365= (5) Panels			
B	XX- For SURECHEK use color code for specific part number i.e. 93 - Mint (299993)			
C	Qty. (2) per panel			