

Treatment Table 8500



ISO 9001-2008 Certified

Customize your treatment table!

Our tables offer versatility, comfort and strength. Build onto our basic models using shelves, cabinets, and drawers to meet your specific needs.



8500



8500-SH



8500-C1



8500-D1

KEY FEATURES

- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

TREATMENT TABLE 8500

- 1 Heavy gauge expanded vinyl & medium density foam for superior comfort.
- 2 Standard paper holder.
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 5 OPTIONAL - Fully laminated cabinet and/or drawer



SPECIFICATIONS

	8500
Overall Height*	30" (76cm) STANDARD*
Overall Width	28" (71cm)
Overall Length	72" (183cm)
Weight**	74 lbs. (33.6kg)**
Weight Capacity	400 lbs. (181kg)

*Special height legs available in 19", 25", & 36"
 **Weight varies depending on table options ordered.
 Add 25 lbs. for 8500-SH

OPTIONS



TB-133 **TB**



*Special Height Legs



Shelf **SH**



Cabinet **C1**
 Double Cabinet **CC**
 Cabinet & Drawer **CD**



Drawer **D1**
 Double Drawer **DD**



Matching Pillow
MP00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS