

## Classic 500 Series, Straight Line Treatment Table with Shelving



### MODEL 1035

- No upcharge for adjustable backrest
- Hardwood legs, laminate shelf, storage compartments and adjustable backrest
- Choose from 2 top widths: 27" & 30"
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Scratch resistant topcoat finish
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- Available in Natural finish only
- **500 lbs. load capacity under normal use (181.43 kg)**

### SPECIFICATIONS


1035-27	72" L x 27" W x 31" H
1035-30	72" L x 30" W x 31" H
1035-38	78" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



20–small Pillow



034– Paper Dispenser & Cutter Combo




040– Paper Cutter



030– Paper Dispenser



088–Flat Top



009–Nose Hole



006–Special Length  
48" or 60" lengths available on 24", 27" & 30" widths  
78" length is only available as item number: 1035-38

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black