



## Vital Sign Monitor Wall Mounted Bracket Series

ADIME912-WM-01 ADIME912-WM-03

ADIME912-WM-06 ADIME912-WM-08



**ADIME912-WM-01**  
Spot Monitor  
Wall Mount Bracket  
with Welch Allyn Bracket



**ADIME912-WM-03**  
Spot Monitor  
Wall Mount Bracket  
for Edan Monitor



**ADIME912-WM-06**  
Spot Monitor  
Wall Mount Bracket  
with Edan M3 Bracket



**ADIME912-WM-08**  
Spot Monitor  
Wall Mount Bracket  
with Omni Bracket

Our Spot Monitor Wall Mount Brackets provide secure, stable, and ergonomic placement for patient monitoring devices in clinical environments where space efficiency and workflow organization matter. This collection includes four configurations to match leading monitor platforms: Welch Allyn Connex Spot Monitor, Edan "M" and "iM" series, Edan M3 Vital Signs Monitors, and Viva Diagnostics Omni Spot units. Each bracket is built from durable aluminum and ABS to withstand daily clinical use and features smooth vertical maneuverability so staff can raise or lower the monitor to the optimal working height. Once positioned, the bracket locks firmly into place, ensuring stability during readings and interaction.

Every bracket is model-matched for proper fit and arrives with a rear plastic accessory basket to keep cuffs, probes, and frequently used components within immediate reach. The Welch Allyn model supports Connex Spot and 4400 monitors; the Edan Monitor bracket fits M and iM series as well as X-series; the Edan M3 bracket fits M3/M3A/M3B/iM3/iM3S; and the Omni bracket is compatible with Viva Omni Spot unit.

## FEATURES

**Compatibility:** Each bracket is built to match the mounting interface of its respective monitor series to ensure secure and stable attachment: Welch Allyn Bracket: Connex Spot Monitor & 4400 series; Edan Monitor Bracket: M, iM, and X series models; Edan M3 Bracket: M3 / M3A / M3B / iM3 / iM3S; Omni Bracket: Viva Diagnostics Omni Spot.

**Adjustable:** Monitor height can be smoothly raised or lowered as needed for patient evaluation.

**Support Frame:** Aluminum and plastic construction provides a strong, impact-resistant structure suitable for ongoing clinical environments.

**Accessory Basket:** Integrated plastic basket keeps cuffs, probes, and small accessories immediately accessible where needed.

**Compact:** Maintains clear wall space and organizes monitoring areas without introducing clutter.

**Stability:** Once locked, the bracket prevents unwanted shifting while clinicians perform readings and documentation.

**Streamlined:** Structural layout supports neat routing of power and data cables to reduce tangling and snag points.

**Sold Separately:** Bracket system is supplied independently for use with facilities' existing monitoring equipment.

# FEATURES



## Adjustable

Monitor height can be smoothly raised or lowered as needed for patient evaluation.



## Stable

Once locked, the bracket prevents unwanted shifting while clinicians perform readings.



## Streamlined

Structural layout supports neat routing of power and data cables to reduce tangling and snag points.



## Support Frame

Aluminum and plastic construction provides a strong, impact-resistant structure.



## Accessory Basket

Integrated included plastic basket keeps cuffs, probes, and small accessories immediately accessible.



## Compact

Clears floor space to minimize clutter and prevent tangling.



*\*Monitors: Sold separately*

Spot Monitor Wall Mount Brackets

SKU:	ADIME912-WM-01	ADIME912-WM-03	ADIME912-WM-06	ADIME912-WM-06
Item Name:	Spot Monitor Wall Mount Bracket with Welch Allyn Bracket	Spot Monitor Wall Mount Bracket for Edan Monitor	Spot Monitor Wall Mount Bracket with Edan M3 Bracket	Spot Monitor Wall Mount Bracket with Omni Bracket
Overall Dimensions:	13" L x 12.2" W x 7.8" D	13" L x 12.2" W x 7.8" D	13" L x 12.2" W x 7.8" D	13" L x 12.2" W x 7.8" D
Total Weight:	7.6 lbs.	9 lbs.	8.3 lbs.	8.2 lbs.
Materials:	Aluminum & Plastic	Aluminum & Plastic	Aluminum & Plastic	Aluminum & Plastic
Height Adjustment:	5" to 13"	5" to 13"	5" to 13"	5" to 13"
Compatibility:	<ul style="list-style-type: none"> <li>• Welch Allyn Connex Spot Monitor</li> <li>• Welch Allyn 4400 Monitor</li> <li>• All standard computer display VESA mounts (no work surface)</li> </ul>	<ul style="list-style-type: none"> <li>• M Series (M80 / M70 / M60 / M50 — legacy models)</li> <li>• iM Series (iM80 / iM70 / iM60 / iM50)</li> <li>• X Series (X12 / X10 / X8i / X15 / iX12 / iX10)</li> </ul>	<ul style="list-style-type: none"> <li>• Edan Vital Signs Monitors: M3 / M3A / M3B / iM3 / iM3S</li> </ul>	<ul style="list-style-type: none"> <li>• Viva Diagnostics Omni Spot Vital Signs Monitor</li> </ul>

# COMPATIBILITY

\*MONITORS: SOLD SEPARATELY



**ADIME912-WM-01**  
Spot Monitor  
**Wall Mount Bracket**  
with Welch Allyn Bracket

Compatibility
Welch Allyn Connex Spot Monitor
Welch Allyn 4400 Monitor
All standard computer display VESA mounts (no work surface)



**ADIME912-WM-03**  
Spot Monitor  
**Wall Mount Bracket**  
for Edan Monitor

Compatibility
M Series (M80 / M70 / M60 / M50 - legacy models)
iM Series (iM80 / iM70 / iM60 / iM50)
X Series (X12 / X10 / X8i / X15 / iX12 / iX10)



**ADIME912-WM-06**  
Spot Monitor  
**Wall Mount Bracket**  
with Edan M3 Bracket

Compatibility
Edan Vital Signs Monitors: M3 / M3A / M3B / iM3/ iM3S



**ADIME912-WM-08**  
Spot Monitor  
**Wall Mount Bracket**  
with Omni Bracket

Compatibility
Viva Diagnostics Omni Spot Vital Signs Monitor