

## Upholstered Apron Couch with Double Drawer Storage


### MODEL 3713

- Easy-clean all laminate construction
- Adjustable pillow wedge headrest
- Cabinet style base with laminate & steel drawers (Available with optional -055 locks)
- Built-in paper dispenser
- 2" of firm foam padding
- Heavyweight, knitbacked upholstery
- 300 lbs. load capacity under normal use

**NOTE:** Paper cutter option, 040, not available for model 3713-10



California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



**3713-16 Adjustable pillow wedge headrest**

72"L x 27" W x 20"H



**3713-15 w/Flat Foam Adjustable Headrest**


72"L x 27" W x 20"H



**3713-10w/Non-Adjustable Pillow Wedge**

72"L x 27" W x 20"H

### Options



040 Paper Holder




055 Drawer Locks

### Laminate Colors



Maple



Dark Cherry



Gray



Slate Gray

### Standard Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black