

## Bonkers for Bunnies Kid Couch with Sliding Doors

### MODEL 3770-B

- Fun couch especially designed for elementary schools, pre-school and daycare facilities
- Exclusive Clinton designs
- Graphics embedded in the laminate on both doors
- Slate Gray laminate end panels with Wedgewood upholstery
- Easy-clean all laminate construction
- Various headrest options
- Cabinet style base with sliding doors on one side
- Built-in paper dispenser
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- 300 lbs. load capacity under normal use

**NOTE:** Paper cutter option, 040, not available for model 3770-10-B

**NOTE:** If adding 055 locks, a quantity of 2 must be purchased




California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



040 Paper Holder



055 Door Locks



3770-16-B Adjustable pillow wedge headrest

72"L x 27" W x 20"H



3770-15-B w/Flat Foam Adjustable Headrest

72"L x 27" W x 20"H



3770-10-B w/Non-Adjustable Pillow Wedge

72"L x 27" W x 20"H

### Graphics



### Upholstery Color



Wedgewood