

KD Panel Leg Series Couch with Shelf

MODEL 3920-27

- All laminate construction
- Sleek, easy-clean, panel leg design
- Top is easily replaceable
- Full shelf storage
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 300 lbs. load capacity under normal use (181.43 kg)



SPECIFICATIONS

3920-27

72" L x 27" W x 21.5" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.


2



Options



034– Paper Dispenser & Holder Combo



040– Paper Holder



030– Paper Dispenser



020– Small Pillow



060– Large Pillow


Laminate Colors



Gray



Maple



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black