

## The CardioChek® PA Analyzer

The CardioChek PA is a cost-efficient, user-friendly diagnostic point-of-care analyzer, which accurately determines multiple patient values.

## **Features & Benefits**



Results in less than 2 minutes

# PORTABLE

Handheld, battery-operated

## PATIENT-FRIENDLY

Small sample (fingerstick) 15-40 uL

## QUALITY

CLIA-waived, FDA-cleared, CE-marked; CRMLN-certified for CHOL and HDL test strips; Meets National Cholesterol Education Program (NCEP) Guidelines

# HIGH RETURN ON INVESTMENT (ROI)

Rapid throughput and low waste, with minimal start-up cost



Customize your diagnostic offerings with a variety of PTS Panels test strips for use with the CardioChek \* PA analyzer. Our test strips offer an easy, efficient method for point-of-care diagnostic testing and Health and Wellness Screening.

**Multi-Analyte** 



### LIPID Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Triglycerides Calculates: LDL Cholesterol, Total Cholesterol/HDL Ratio



#### CHOL+HDL+GLU Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Glucose Calculates: Total Cholesterol/HDL Cholesterol Ratio



#### Panel Test Strip

Measures: Total Cholesterol & Glucose



#### CHOL+HDL Panel Test Strip

Measures: Total Cholesterol & HDL Cholesterol Calculates: Total Cholesterol/HDL Cholesterol Ratio



#### TOTAL CHOLESTEROL

Test Strip

HDL CHOLESTEROL

Measures: Total Cholesterol

Test Strip Measures: HDL Cholesterol



#### TRIGLYCERIDES

Test Strip

Measures: Triglycerides



Test Strip

Measures: Glucose

Other test strips available include PTS Panels LDL Cholesterol (direct, not calculated), Creatinine (special order only), and Ketone (special order only).

