



## Adscope® 608 Convertible

AFD technology and convertible design make this the most versatile of our clinician scopes.

### The 608 series features

- ADC's adjustable frequency design (AFD) technology provides the acoustic response of a traditional bell and diaphragm in a convenient two-sided chestpiece
- Extra-large bell (1 1/4" diameter) for unsurpassed low-frequency response
- Ultra-sensitive diaphragm (1 1/4" diameter) for amplification and crisp high-frequency response
- AFD technology further enhances acoustic response: use light pressure for low frequencies, firm pressure for high frequencies
- Non-chill bell and diaphragm retaining rims provide patient comfort
- Clinician headset boasts reinforcing yoke molded into flexible 22" PVC binaurals fixed at a 15° angle for greater comfort
- Adsoft™ Plus silicone eartips for the ultimate in wearing comfort and acoustic seal
- Accessory kit with spare Adsoft™ Plus eartips and scope ID tag
- Weighs 5.8 oz (164.4 g)
- Overall length 31"
- Assembled, inspected, and packaged in the USA
- Lifetime warranty



AFD Technology Diaphragm



Convert pediatric diaphragm to adult bell in seconds

### Ordering Information

Item	Description
<b>608**</b>	Adscope® 608, Convertible Clinician Scope, (Specify Color)

### Color Choices for 608\*\* Only

<b>BK</b>	<b>BD</b>	<b>DG</b>	<b>FG</b>	<b>FL</b>
<b>FRB</b>	<b>FS</b>	<b>FV</b>	<b>G</b>	<b>LV</b>
<b>MCB</b>	<b>MCA</b>	<b>MG</b>	<b>MRS</b>	<b>PBCA</b>
<b>N</b>	<b>O</b>	<b>P</b>	<b>R</b>	<b>RB</b>
<b>TL</b>	<b>TQ</b>	<b>V</b>		

**NOTE:** To ensure correct color when ordering, add color suffix to model number (e.g., 608**DG**)

