

# Home Health Care Product Guide



# THE NEW AMERICAN STANDARD

ADC medical instruments and accessories have been carefully designed to meet the rigorous demands of the most discriminating healthcare professionals throughout the world.

In use on six continents in over twenty different countries, ADC products enjoy a reputation for advanced design combined with skilled craftsmanship.

ADC products are manufactured from the finest materials with an uncompromising concern for quality. Exhaustive inspection of every item by our factory-trained quality control staff ensures strict compliance with ADC standards for performance and durability. And the commitment to manufacturing excellence is backed by the industry's most extensive warranties.

“State-of-the-art” distribution combined with a knowledgeable, courteous staff guarantees that the service you receive will be as outstanding as the products we deliver.

So compare the quality, performance, and most importantly, the value of ADC products to any others you can buy.

You'll discover that there's really no comparison at all!

## CONTENTS

Blood Pressure .....	3-7
Stethoscopes .....	8-9
Thermometry .....	10-11
Instruments .....	12-13
Caseware .....	14-15
Retail Displays .....	16



# Blood Pressure



## PROSPHYG™ 760 SERIES

Our most popular aneroid, the 760 is the model of choice for demanding institutions.

### The 760 series features:

- ≡ 300mmHg manometer, certified to  $\pm 3$ mmHg
- ≡ Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing
- ≡ Dip molded latex inflation bladder & bulb withstands heavy use
- ≡ Chrome-plated air release valve
- ≡ Leatherette carrying case with durable nylon zipper
- ≡ 20 year calibration warranty
- ≡ Assembled, inspected, and packaged in the U.S.A. from foreign components



760 Prosplyg™ Series



775 Prosplyg™ Series

## PROSPHYG™ 775 SERIES

An economical aneroid sphygmomanometer for today's cost conscious healthcare market.

### The 775 series features:

- ≡ 300mmHg manometer, certified to  $\pm 3$ mmHg
- ≡ Navy blue nylon cuff with range markings to facilitate selection of correct cuff size
- ≡ Latex inflation bladder and bulb
- ≡ Chrome-plated deflation valve
- ≡ Complete with leatherette zippered carrying case and operating instructions
- ≡ 10 year calibration warranty
- ≡ Assembled, inspected, and packaged in the U.S.A. from foreign components



**768-641** Pros Combo™ II Series

## 768-641 PRO'S COMBO II POCKET ANEROID KIT

A professional kit that combines our best selling aneroid, with an optional stethoscope in a color coordinated premium nylon case.

### The 768-641 series features:


- ≡ 300mmHg manometer, certified to  $\pm 3$ mmHg
- ≡ Adcuff™ calibrated nylon cuff with ADC's proprietary Size Guide™ marking system prevents mis-cuffing
- ≡ Latex inflation bladder and bulb
- ≡ Chrome-plated deflation valve
- ≡ Extra large, color-matched, dual pocket carry case, measures 10" x 6" x 3"
- ≡ Color matched sprague scope
- ≡ Assembled, inspected, and packaged in the U.S.A. from foreign components

## ADTOONS™ 760 SERIES

Same as our Standard 760 Series, except cuff & zippered case are available in pattern matched Adtoon™ fabric prints.



**760** Adtoon™  
Prophyg™ Series



**768-641** Adtoons™ Pro's Combo™ Series

## ADTOONS™ PRO'S COMBO II POCKET ANEROID KIT

Same as our Best Selling 768-641 Series, except cuff & zippered case are available in pattern matched Adtoon™ fabric prints and color coordinated sprague.



**780**  
Prosplyg™ Series

## PROSPHYG™ 780 SERIES

An easy-to-use instrument, including everything needed for home blood pressure measurement.

### The 780 series features:

- ≡ 300mmHg manometer, certified to  $\pm 3$ mmHg
- ≡ Navy blue nylon D-Ring self-adjusting cuff for one-handed cuff application
- ≡ Latex inflation bladder and bulb
- ≡ Chrome-plated deflation valve
- ≡ Diaphragm only blood pressure scope with black PVC tubing (#660BK)
- ≡ Complete with leatherette zippered carrying case
- ≡ 10 year calibration warranty
- ≡ Assembled, inspected, and packaged in the U.S.A. from foreign components

## PROSPHYG™ 790 SERIES

Same as above, except diaphragm chestpiece attached to cuff for simplified use.



**790**  
Prosplyg™ Series



**6005**  
Manual BP Kit

## MANUAL BP KIT 6005

The most feature rich manual kit on the market.

### The 6005 series features:

- ≡ Dial features bold white numerals on black face for easy reading
- ≡ Self-adjusting D-Ring Adcuff™ for Home™ for one-handed application
- ≡ Size Guide™ marking system ensures optimal accuracy
- ≡ Latex-free
- ≡ Auto-deflation valve
- ≡ Accurate to  $\pm 3$ mmHg of reading
- ≡ Attached nurse scope
- ≡ Zippered carrying case and operating instructions

## ADVANTAGE™ 6012 SEMI-AUTOMATIC BP MONITOR

Compact, and easy-to-use

### The 6012 features:

- ≡ Measures systolic, diastolic pressure, and pulse
- ≡ Manual inflation
- ≡ Self adjusting latex-free D-Bar cuff
- ≡ Automatic deflation valve
- ≡ Auto-off



6012 Advantage™ BP



6013 Advantage™ BP

## ADVANTAGE™ 6013 AUTOMATIC BP MONITOR

Smart Logic™ Technology for the ultimate in simplicity.

### The 6013 features:

- ≡ Measures systolic, diastolic pressure, and pulse
- ≡ Smart Logic™ Technology determines ideal inflation level
- ≡ Automatic inflation and deflation control
- ≡ Self adjusting latex-free D-Bar cuff
- ≡ Recalls last 14 readings with date and time stamp
- ≡ Auto-off
- ≡ Optional AC adapter (sold separately)

## ADVANTAGE™ 6014 ADVANCED BP MONITOR

Feature packed with average mode for truly consistent readings.

### The 6014 features:

- ≡ Measures systolic, diastolic pressure, and pulse
- ≡ Smart Logic™ Technology determines ideal inflation level
- ≡ Automatic inflation and deflation control
- ≡ Average mode automatically averages 3 readings for more meaningful results
- ≡ Self-adjusting latex-free D-Bar cuff
- ≡ Advanced 2 zone memory with date and time stamp recalls last 30 readings for two different people
- ≡ Auto-off
- ≡ Optional thermal printer (sold separately)
- ≡ Optional AC adapter (sold separately)



6014 Advantage™ BP



6014P - Optional Advantage™ Thermal Printer

Refer to full catalog for all model options

## ADVANTAGE™ 6015 DIGITAL WRIST BP MONITOR

Our most popular home care model, the 6015 combines quality construction with ease of use.

### The 6015 features:

- ≡ Measures systolic, diastolic pressure, and pulse
- ≡ Smart Logic™ Technology determines ideal inflation level
- ≡ Latex-free wrist cuff fits 5" to 7" wrists
- ≡ Recalls last 90 readings
- ≡ Auto-off
- ≡ Compact storage case



6015 Advantage™ BP



6016 Advantage™ BP

## ADVANTAGE™ 6016 DIGITAL WRIST BP MONITOR

Feature packed with average mode for truly consistent readings.

### The 6016 features:

- ≡ Measures systolic, diastolic pressure, and pulse
- ≡ New 8 bit microprocessor and advanced algorithm for more accurate readings
- ≡ Smart Logic™ Technology determines ideal inflation level
- ≡ Average mode automatically averages 3 readings for more meaningful results
- ≡ Latex-free wrist cuff fits 5" to 7" wrists
- ≡ Advanced 2 zone memory with date and time stamp recalls last 30 readings for two different people
- ≡ Auto-off
- ≡ Compact storage case
- ≡ Optional thermal printer (sold separately)



6014P - Optional Advantage™ Thermal Printer

## HOME BLOOD PRESSURE PRODUCT FEATURE CHART

Item	Type	Inflation	Smart Logic	Deflation	Scope	Latex	Pulse	Memory	Battery	Available Printer	Available A/C Adapter	Warranty
780	Aneroid	Manual	N/A	Manual	Yes	Latex	N/A	N/A	N/A	N/A	N/A	10 year
790	Aneroid	Manual	N/A	Manual	Attached	Latex	N/A	N/A	N/A	N/A	N/A	10 year
6005	Aneroid	Manual	N/A	Auto valve	Attached	Latex-free	N/A	N/A	N/A	N/A	N/A	5 year
6012	Dig., Arm	Manual	N/A	Auto valve	N/A	Latex-free	Yes	1	4 "AA"	No	No	5 year
6013	Dig., Arm	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	14	4 "AA"	No	Yes	5 year
6013SA	Dig., Arm	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	14	4 "AA"	No	Yes	5 year
6013X	Dig., Arm	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	14	4 "AA"	No	Yes	5 year
6014	Dig., Arm	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	2x30	4 "AA"	Yes	Yes	5 year
6015	Dig., Wrist	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	90	2 "AAA"	No	No	5 year
6016	Dig., Wrist	Automatic	Yes	Elec. Auto	N/A	Latex-free	Yes	2x30	2 "AAA"	Yes	No	5 year

# Stethoscopes



## PROSCOPES™

ADC Proscopes™ combine lightweight comfort, sensitive acoustics, and handsome styling, with truly affordable pricing.

### The Proscope™ series features:

- ≡ Choice of diaphragm only (660 series) or dualhead (670 series)
- ≡ Proprietary diaphragm enhances acoustic response by 50%
- ≡ Flexible 22" PVC tubing
- ≡ Inspected and packaged in the U.S.A.
- ≡ One Year Warranty

### Proscope™ 660 series features:

- ≡ Non-chill diaphragm only nurse style aluminum chestpiece (1¼" dia.)
- ≡ Weighs just 3 oz. / Overall length 32"

### Proscope™ 670 series features:

- ≡ Combination diaphragm/bell aluminum chestpiece for wide frequency response
- ≡ Weighs just 3.5 oz. / Overall length 32.5"



660 / 670 Proscope™ Series



647 Sprague Series

## SPRAGUE RAPPAPORT TYPE SCOPE

Our exclusive design features an improved valve mechanism that virtually eliminates acoustic leakage and a chestpiece key to facilitate after market service.

### The Sprague features:

- ≡ Threaded chestpiece drum permits the use of five inter-changeable chestpiece fittings - an adult and pediatric diaphragm, and three bells (1¼", 1", ¾") - for complete diagnostic capability
- ≡ Proprietary valve mechanism eliminates acoustic leakage
- ≡ 9 piece accessory kit
- ≡ Inspected and packaged in the U.S.A.

### Proscope™ 640 features:

- ≡ Traditional 16" tubing length in black only
- ≡ One Year Warranty

### Adscope™ 641 features:

- ≡ Convenient 22" tubing length in choice of colors
- ≡ Five Year Warranty

### Adscope™ 647 Sprague I features:

- ≡ Latex-free extra thick "Y" tubing in choice of colors
- ≡ Five Year Warranty



640 / 641 Sprague Series

Refer to full catalog for all model options



## ADSCOPE™ 609

The ADSCOPE-lite™ is the only affordably priced stethoscope to combine superior acoustic response with an ultra lightweight design.

### The 609 Series features:

- ≡ Combination adult chestpiece precisely machined to exacting tolerances from lightweight aluminum
- ≡ Non-chill bell and diaphragm retaining rim for patient comfort
- ≡ Custom molded 22" PVC tubing
- ≡ Adsoft™ PVC eartips and spare diaphragm included
- ≡ Weighs just 4.5 oz. / Overall length 31"
- ≡ Five Year Limited Warranty



609 Adscope™ Series



603 Adscope™ Series

## ADSCOPE™ 603

One of our most popular professional acoustic scopes, the 603 combines outstanding acoustic performance with rugged durability.

### The 603 features:

- ≡ Combination adult size chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- ≡ Stainless steel binaurals fixed at 15° angle
- ≡ Adsoft™ deluxe PVC eartips, the ultimate in comfort and acoustic seal
- ≡ Spare diaphragm & scope ID tag included
- ≡ Weighs 6.5 oz. / Overall length 31"
- ≡ Limited Lifetime Warranty

## ADSCOPE™ 602

One of the finest acoustic scopes available, the 602 Cardiology model combines precision crafted performance with rugged durability.

### The 602 Series features:

- ≡ Combination chestpiece precisely machined to exacting tolerances from surgical stainless steel
- ≡ Non-chill bell and diaphragm retaining rim for patient comfort
- ≡ Large internal bore stainless steel binaurals are fixed at 15°
- ≡ Custom molded bi-lumen 16" PVC tubing stays soft and flexible
- ≡ Adsoft™ deluxe PVC eartips in small and large sizes for ultimate in wearing comfort and acoustic seal
- ≡ Spare diaphragm & scope ID tag included
- ≡ Weighs 7.5 oz. / Overall length 25.5"
- ≡ Limited Lifetime Warranty on metal parts



602 Adscope™ Series

# Thermometry



## ADTEMP™ II 413 THERMOMETER

Our most popular model.

### ADTEMP™ II Thermometer features:

- ☞ Compact size - under 5" long, weighs just 11 grams
- ☞ Measures in 60 seconds
- ☞ Audible "Beep" when measurement is complete
- ☞ Auto-off function conserves battery life
- ☞ Range 90°F-109.9°F ±.2°F / 32°C - 43.9°C ±.1°C
- ☞ Replaceable 1.55v (LR41) battery provides 1,500 measurements
- ☞ Includes 5 sheaths & clear plastic case to keep instrument sanitary



413 Adtemp™ II Series



414 Adtemp™ III Series

## ADTEMP™ III FLEX 414 THERMOMETER

Flexible probe for patient safety.

### ADTEMP™ III Flex Thermometer features:

- ☞ Same as above 413 except additional flexible probe feature for added patient safety



415 Adtemp™ IV

## ADTEMP™ IV 415 DIGITAL THERMOMETER

World's fastest peak hold pen-type thermometer.

### The ADTEMP™ IV thermometer features:

- ≡ 18 second oral, 10 second rectal temperature measurement
- ≡ Auto-off function conserves battery life
- ≡ Large LCD for easy reading
- ≡ Range 90°F-109.9°F ±.2°F or 32°C - 43.9°C ±.1°C depending upon scale selection
- ≡ Replaceable 1.55v (LR41) type battery provides up to 1,500 measurements
- ≡ Includes 5 probe sheaths & clear 2 piece plastic storage case

## ADTEMP™ V 418 SUPER FAST DIGITAL THERMOMETER

The fastest pen-type digital thermometer on the market using predictive technology.

### ADTEMP™ V Super Fast Flex thermometer features:

- ≡ 10 second oral, 8 second rectal measurement using proprietary predictive technology
- ≡ Auto-off function conserves battery life
- ≡ Range 90°F-109.9°F ±.2°F or 32°C - 43.9°C ±.1°C depending upon scale selection
- ≡ Replaceable 1.55v (LR41) type battery provides 1,500 measurements
- ≡ Integral carry case
- ≡ Includes 5 probe sheaths



418 Adtemp™ V



421 Adtemp™ Tympanic

## ADTEMP™ 421 DIGITAL EAR THERMOMETER

The most compact Tympanic thermometer on the market.

### Thermometer features:

- ≡ Measures core temperature instantly from the ear reliably and accurately
- ≡ Range 89.6°F-108°F ±.4°F or 32°C-42.2°C ±.2°C depending on scale selection
- ≡ Probe cover FREE design is less costly and more convenient to use
- ≡ Large 6mm back lit display is easy to read
- ≡ Compact (5" Long, 1.9 oz.) ergonomic design fits comfortably in operator's hand and patient's ear
- ≡ Fever alarm at 99.5°F (37.5°C)
- ≡ 12 reading memory
- ≡ Complete with Lithium battery, storage/travel case, and operating instructions

# Instruments



## MEDICUT™ SHEARS

Available in two sizes, ADC's Medicut™ Shears are the world's best selling brand of medical grade utility shears.

### Medicut™ Shears feature:

- ⌘ Razor sharp, super-hardened surgical stainless steel blades with one serrated edge to cut through the toughest materials
- ⌘ Contoured polypropylene handles with large ring provide maximum control and comfort
- ⌘ Fully auto-clavable to 143°C (290°F)
- ⌘ Available in traditional 7¼", or smaller 5½" length
- ⌘ Inspected and packaged in the U.S.A.
- ⌘ Two Year Quality Performance Guarantee

### Medicut™ Shears feature:

- ⌘ Traditional 7¼" length for maximum leverage and strength
- ⌘ Stainless steel safety bandage tip

### MiniMedicut™ Shears feature:

- ⌘ Smaller 5½" length is ideal for cutting bandages, tape, and clothing
- ⌘ Polypropylene safety bandage tip color coordinated to the handle



320 Medicut™ Series  
321 Mini Medicut™ Series

## SURGICAL INSTRUMENTS

Forged from corrosion resistant surgical stainless steel, these popular floor grade instruments are built to last.

### 301 Lister Bandage Scissors feature:

- ⌘ Traditional 5½" length for maximum leverage and strength
- ⌘ Stainless steel safety bandage tip

### 310 Kelly Forceps (Straight)

### 323 Listerette Scissors feature:

- ⌘ Traditional 5½" length for maximum leverage and strength
- ⌘ Stainless steel safety bandage tip
- ⌘ Contoured polypropylene handles with large ring provide maximum control and comfort



301 Lister  
Bandage Scissors



310 Kelly Forceps  
(Straight)



323 Listerette Scissors

Refer to full catalog for all model options



351 Adlite™

## ADLITE™

### The 351 series features:

- ≡ Sealed concentric beam, disposable diagnostic penlight
- ≡ Activated by depressing pocket clip
- ≡ Available with or without pupil gauge
- ≡ Measures 5" long, ½" diameter
- ≡ Sold in 6 pack

## REUSABLE PENLIGHTS

Two distinctly different models to suit your personal preferences.

### The Metalite™ features (#352):


- ≡ Brushed aluminum housing
- ≡ Twist/lock activated
- ≡ Measures 5¼" long by ½" diameter
- ≡ Two "AAA" batteries included

### The Adlite II™ features (#354):

- ≡ Plastic housing for lightweight comfort
- ≡ Pocket clip activated
- ≡ Measures 5" long by ½" diameter
- ≡ Two "AAA" batteries included



352 Metalite™  
354 Adlite™



394 EKG Ruler  
395 EKG Caliper

## EKG RULER & CALIPER

Two of the most popular electrocardiograph accessories.

### The EKG Ruler features:

- ≡ Seconds scale for PR, QRS, Cycle, Length, Coupling, Escape, and Parasystolic Intervals
- ≡ Measures just 1 1/16" x 6 3/8" small enough to fit in a medical coat pocket
- ≡ Complete with protective sleeve

### The EKG Caliper features:

- ≡ Hardened steel arms and points for accurate measurement of tracing intervals
- ≡ Point protector has built-in 40mm metric scale
- ≡ Measure 4 3/4"



## POCKET PAL II™ ORGANIZER

An indispensable uniform accessory that combines fashion and function.

### The 216 series features:

- ≡ Washable 420D white nylon pocket organizer sewn with nylon thread
- ≡ 5 slide-in and one covered pocket in front, one full-width slide-in pocket in rear
- ≡ Handy key ring
- ≡ Velcro® snap closures in choice of 8 vibrant colors
- ≡ Measures 6¼" x 4¾" x ¼", Weighs 1 oz.

### Pocket Pal II™ Kit includes:

- ≡ Pocket Pal II™ (#216)
- ≡ 5½" Lister Bandage Scissor (#301)
- ≡ 3 color pen (#392)
- ≡ Adlite™ Disposable Penlight (#350)

### Pocket Pal III™ Kit includes:

- ≡ Pocket Pal II™ (#216)
- ≡ MiniMedicut™ Shear (#321)
- ≡ 3 color pen (#392)
- ≡ Adlite™ Plus (#356)



- 216** Pocket Pal™ II Organizer Series
- 116** Pocket Pal™ II Kit Series
- 117** Pocket Pal™ III Kit Series (not shown)



**116-641** Nurse Combo™ Series

## NURSE COMBO

Combines the Pocket Pal II™ Kit with color coordinated Adscope™ 641 Sprague scope.

### The Nurse Combo features:

- ≡ The Pocket Pal II™ Kit complete with three instruments and choice of 9 Velcro® tag colors
- ≡ Adscope™ 641 Sprague scope with 22" tubing color coordinated to the Velcro® tags and a 9 piece accessory kit
- ≡ Assembled, inspected, and packaged in the U.S.A.



**1026** Adtoon™ Nurse Carry Tote

## NURSE CARRY TOTE

A convenient carry case for all your professional gear that is as colorful as it is practical.

### The Adtoon™ Nurse Carry Tote features:

- ☞ 210 denier nylon with heavy liner for rugged durability
- ☞ Choice of 4 Adtoon™ Prints
- ☞ Zipper outside pocket
- ☞ 2 mesh inner pockets
- ☞ Padded shoulder strap
- ☞ Measures 12" x 10½" x 5"
- ☞ Weighs 9.6 oz.
- ☞ One year warranty

## NYLON MEDICAL BAG

A heavy duty, padded medical nylon bag to protect your most valuable medical equipment.

### The Medical Bag features:

- ☞ A multi-pocketed, heavy-duty nylon medical bag
- ☞ Double zipper opening for main compartment, 10½" x 9½" x 4¼"
- ☞ Two interior zippered pockets 6½" x 4"
- ☞ Includes an adjustable, removable webbed shoulder strap, and two contoured handles for easy carrying
- ☞ 11 interior gusseted pockets of varying sizes
- ☞ Back file folder pocket 13¾" x 8"
- ☞ Front zippered pocket measures 6½" x 9" x 4"
- ☞ Measures 13½" x 9" x 5½" (closed) / Weighs 1 lb
- ☞ One year warranty

*(Note: Items shown in caseware photo are not included)*



**1024** Medical Bag



**880** Adtoon™ Carry Case (not shown)  
**888** Adtoon™ Premium Carry Case

## ADTOON™ CARRY CASES

- ☞ Choose from standard or premium zippered carry cases in choice of pattern
- ☞ 210 denier nylon construction
- ☞ Rugged self repairing nylon zipper with oversized pull
- ☞ Premium case has one large zipper pocket, end strap, and handy key ring.

## RETAIL DISPLAYS

We've taken some of our best selling products and made them available in attention grabbing display packaging that highlights our color coordinated program of diagnostic products, medical instruments, and accessories.

### Counter Display

- ≡ Two tier counter unit measures 31" tall
- ≡ 19" turning radius
- ≡ Uses 2.5 sq. ft. of counter space
- ≡ 12 hooks in total
- ≡ Weighs 3.5 lbs.

### Floor Display

- ≡ Four tier floor unit measures 65" tall
- ≡ 21" turning radius
- ≡ 18 hooks in total
- ≡ Weighs 7.1 lbs.

### Planograms

- ≡ Three and four foot planograms for slat walls or pegboards are available to configure the right display for you through our customer service department.



Counter Display



Floor Display



55 Commerce Drive  
Hauppauge, New York 11788  
tel: 1-800-232-2670  
fax: 631-273-9659

[www.adctoday.com](http://www.adctoday.com)



HCCAT118  
Printed in U.S.A.