

ADCUFF™ SPU

Single-Patient-Use Blood Pressure Cuffs



Better DESIGNED

Our line of single-patient-use blood pressure cuffs has been completely redesigned with a focus on precision, comfort, and value.

Curved Corners
for greater patient comfort

Range Markings
for accurate sizing,
minimizes risk of
miscuffing


**Strong Hook-and-Loop
Adhesive**
withstands repeated
open/close cycles



Index Line
for accurate sizing,
minimizes risk of miscuffing

Cuff ID Panel
provides vital cuff information
at a glance

Size Icons and Color Coding
for quick visualization
of selected cuff



Adult - 11
27.5 - 36.5 cm

1 Tube	REF: 8450-11A-1HP
2 Tube	REF: 8450-11A-2SC

Printed Model Numbers
simplify stock control
and reordering

**Properly Positioned
Artery Mark**
for optimized arterial compression
(you'd be surprised how many aren't)

**Standardized connections mean fewer model options,
improved stock control, and simplified inventory management.**

**All in a flexible, non-proprietary solution that's compatible
with most major monitor brands.**

- Thirteen sizes for neonatal, child, and adult patients
- One-piece construction with a soft and flexible finish
- 100% individually leak tested
- Economical enough for single-patient application
- Durable enough for repeated use

SPU

Cuff Reference Guide





ADC's single-patient-use cuffs are available in 13 clinical sizes. Cuffs are color coded and available in both one- and two-tube configurations.

Adcuff SPU Blood Pressure Cuffs

- Economical enough for single-patient application
- Durable enough for repeated use
- 100% leak tested to assure dependability
- Five clinical cuff sizes for neonatal patients
- Eight clinical cuff sizes for child and adult patients
- The industry's most precise cuff markings
- Cuff design minimizes miscuffing errors
- Rounded corners and soft, woven construction
- One-tube configuration with bayonet-type (HP) connector (part 8974)
- Two-tube configuration with screw-type (Dinamap) connector (part 8971)
- Neonatal cuffs with Luer connector (part 8975)
- Standardized connectors for improved inventory management
- Add low-cost adapters directly to monitor for a long-lasting, one-time conversion
- Compatible with most brands of monitors

Adcuff™ SPU Cuffs and Accessories

	Connector	Infant	Child	Sm. Adult	Adult	Adult Long	Lg. Adult	Lg. Adult Long	Thigh
Range (cm)	N/A	9-14.8	13.8-21.5	20.5-28.5	27.5-36.5	27.5-40	35.5-46	35.5-50	45-56.5
Color	N/A	Orange	Green	Royal Blue	Navy	Navy	Burgundy	Burgundy	Brown
One-Tube	None	8450-71-1	8450-9C-1	8450-10SA-1	8450-11A-1	8450-11AL-1	8450-12X-1	8450-12XL-1	8450-13T-1
One-Tube	Bayonet (HP) 	8450-71-1HP	8450-9C-1HP	8450-10SA-1HP	8450-11A-1HP	8450-11AL-1HP	8450-12X-1HP	8450-12XL-1HP	8450-13T-1HP
Two-Tube	None	8450-71-2	8450-9C-2	8450-10SA-2	8450-11A-2	8450-11AL-2	8450-12X-2	8450-12XL-2	8450-13T-2
Two-Tube	Screw (Dinamap) 	8450-71-2SC	8450-9C-2SC	8450-10SA-2SC	8450-11A-2SC	8450-11AL-2SC	8450-12X-2SC	8450-12XL-2SC	8450-13T-2SC
Inflation System*	None	8650-71	8650-9C	8650-10SA	8650-11A	8650-11AL	8650-12X	8650-12XL	8650-13T

* Includes 875PLF bulb and valve

Adcuff™ Neonatal SPU Cuffs and Accessories

	Connector	Neonatal 1	Neonatal 2	Neonatal 3	Neonatal 4	Neonatal 5
Range (cm)	N/A	3-6	4-8	6-11	7-13	8-15
Color	N/A	Orange	Light Blue	Green	Navy	Burgundy
One-Tube	Luer Slip	8450-1N-1L	8450-2N-1L	8450-3N-1L	8450-4N-1L	8450-5N-1L
Two-Tube	Luer Slip	8450-1N-2L	8450-2N-2L	8450-3N-2L	8450-4N-2L	8450-5N-2L




Simply STANDARDIZED

ADC's Adcuff SPU standardized connectivity solution is the simplest yet. While some competitors lock you into a platform, or offer hundreds of complicated and confusing cuff configurations that drive up costs, we start with two choices:

Choose a **one-tube cuff** with a bayonet-type (HP) connector, or a **two-tube cuff** with a screw-type (Dinamap) connector.

Then, if needed, use a straight, X, or Y adapter to connect your cuff of choice to virtually any monitor on the market.



It's a low-cost, long-lasting, one-time conversion.

For 1-Tube Monitors

with Male Marquette Connectors

1-Tube Cuff



HPF-S-MQF
Bayonet (HP) Female
to Marquette Female

2-Tube Cuff



SCF2-Y-MQF
Screw (Dinamap) Female
to Marquette Female

with Female Screw Connectors

1-Tube Cuff



HPF-S-SCM
Bayonet (HP) Female
to Screw (Dinamap) Male

2-Tube Cuff



SCF2-Y-HPM
Screw (Dinamap) Female
to Bayonet (HP) Male

For 2-Tube Monitors

with Male Marquette Connectors

2-Tube Cuff



SCF-S2-MQF
Screw (Dinamap) Female
to Marquette Female

with Female Bayonet Connectors

2-Tube Cuff



SCF-S2-HPM
Screw (Dinamap) Female
to Bayonet (HP) Male

For Custom Configurations

Don't see what you need? ADC offers open-ended adapters and will custom assemble for your facility – **always without charge.**

For 1-Tube Monitors

1-Tube Cuff



HPF-S-N
Bayonet (HP) Female
to customizable open end

2-Tube Cuff



SCF2-Y-N
Screw (Dinamap) Female
to customizable open end

For 2-Tube Monitors

2-Tube Cuff



SCF-S2-N
Screw (Dinamap) Female
to customizable open end



Diagnostic instruments of extraordinary design, exceptional quality, and superior performance – all with the industry’s most comprehensive warranties.

In use on six continents in more than 60 countries, ADC products enjoy a reputation for performance, durability, and value. For more than 30 years, healthcare professionals and institutions have relied on ADC to deliver the right equipment at the right price, with a user experience that’s second to none.

Distributed by:



American Diagnostic Corporation
 55 Commerce Drive
 Hauppauge, NY 11788

tel: 631-273-9600
 toll free: 1-800-232-2670
 fax: 631-273-9659

www.adctoday.com
info@adctoday.com

