

*e-sphyg*3™

Digital Blood Pressure Monitor

- **Accurate**
- **Affordable**
- **Clinical Grade**



ADC's all-new e-sphyg 3 captures blood pressure measurements with the touch of a button. Designed for busy professionals who value both accuracy and affordability, the e-sphyg 3 is the feature-packed, clinical-grade solution you've been waiting for.

Three clinically validated modes:

1x **Single measurement mode** with pulse rate and mean arterial pressure (MAP) displayed in as little as 30 seconds.

3x **Customizable, three-reading average mode** with individual measurement display.

Manual mode for auscultation with a stethoscope.



M Last-reading memory


Rechargeable high-capacity battery STANDARD



The e-sphyg³ is perfect for:



Physician practices



Walk-in clinics



Nursing homes

Desktop design with built-in handle
Or choose one of these optional platforms
to fit your practice workflow.

Oscillometric measurements capture:

SYS Systolic blood pressure
(60 to 255 mmHg)

DIA Diastolic blood pressure
(30 to 200 mmHg)

HR Pulse rate
(40 to 200 BPM)

MAP Mean arterial pressure
(35 to 240 mmHg)

IB Irregular
Heartbeat

For adult and pediatric patients

INCLUDES: Three-cuff kit
(Small Adult, Adult, and
Large Adult) for accurate
readings on virtually
all patients.

Choose our traditional
Adcuff design or our new
Adcuff+ one-piece cuffs.

Wall Mount
Custom bracket for
easy mounting to
the wall or on an
ADC Adstation system.


*Includes cuff
storage basket.*



Mobile Stand

Our rugged, easy-to-assemble
mobile stand features five
rolling casters (two locking) for
room-to-room portability.

Includes cuff storage basket.



Desk Caddy

A modern, compact
design that helps keep
cuffs and cords
organized.



Dentist offices



Birth centers



Long-term care facilities

e-sphyg3™

Digital Blood Pressure Monitor

The clinical-grade solution you've been waiting for.

Technical Comparisons

FEATURE	e-sphyg3	Vital Signs Monitors	Home BP Devices	Manual BP Devices
Automated Measurement	Yes	Yes	Yes	No
Clinical Grade NIBP	Yes	Yes	No	N/A
3 Reading Average Mode	Yes	Some	Some	N/A
Automatic and Manual Modes	Yes	Some	No	No
High-Capacity Battery - STANDARD	Yes	No	No	N/A
Includes 3 cuffs	Yes	No	No	No
Multiple Platforms	Yes	Yes	No	Yes
Irregular Heartbeat	Yes	Some	Some	No



Ordering Information

9003K-MCC e-sphyg 3 with two-piece cuffs

9003K-MCC1 e-sphyg 3 with one-piece cuffs

Both configurations include: e-sphyg 3 monitor; built-in, rechargeable, high-capacity battery; 6-foot power cord with four adapters; three Adcuff blood pressure cuffs (small adult, adult, large adult); and 4.5-foot tubing with connector.

Optional Platforms

9003D e-sphyg 3 desk caddy

9003W e-sphyg 3 wall mount bracket (compatible with ADC Adstation systems) and cuff basket

9003M e-sphyg 3 mobile stand and cuff basket

Replacement Parts

9003ADAPT AC adapter set for esphyg 3

9003BAT High-capacity battery for esphyg 3

9003CK Adcuff kit (small adult, adult, large adult) with e-sphyg 3 connectors

9003CT Extension tubing/connectors for esphyg 3



55 Commerce Drive
Hauppauge, NY 11788

Toll Free: 1-800-ADC-2670
(1-800 232-2670)

Email: info@adctoday.com

Web: www.adctoday.com