



Lucas Recovery Couch Series with Adjustable Back

ADI996-11-BLU ADI996-12-BLU



ADI996-11-BLU

Lucas Recovery Couch
with Adjustable Back



ADI996-12-BLU

Lucas Plus Recovery Couch
with Adjustable Back

Features

Adjustable Backrest: Six positioning angles including 27°, 36°, 43°, 50°, 55°, and 60° to accommodate patient comfort and recovery needs.

Professional Appearance: Antimicrobial Midnight Express upholstery offers a sleek look.

Durable: Reinforced wooden frame built to handle heavy daily use and support up to 420 lbs.

Hygienic: Antimicrobial upholstery is easy to clean and helps maintain a sanitary recovery environment.

Comfortable Support: 2-inch-high density foam cushioning provides lasting comfort for patients.

Optional: Pair with our compatible 27" paper roll dispenser and cutter (sold separately) to streamline hygiene between patients.

Outfitted: Rubber feet on the bottom of the wooden frame and feet provides added stability while protecting both the recovery couch and surrounding floor surfaces from scuffs and wear.

Lucas Recovery Couch

With Adjustable Back & Antimicrobial Upholstery



420 lbs.
weight
capacity

Upholstery Color



Midnight Express

Stability




Wooden Feet



High-Quality Plywood

Upholstery



2-inch Thick Foam



Antimicrobial Leather

Adjustability



Adjustable Back Rest

Add-On



Paper Roll Dispenser



Paper Roll Cutter

Specifications

Weight	60.8 lbs.
Dimensions	72" L x 27" W x 20" H
Weight Capacity	420 lbs.
Materials	PVC Leather, Foam, Solid Wood, Rubber, and Particleboard
Thickness of Foam Padding	2"
Adjustable Positions	27°, 36°, 43°, 50°, 55° & 60°
Upholstery	Antimicrobial
Color	Midnight Express
Feet Dimensions	17" H x 2.4" W x 2.4" D

Lucas Plus Recovery Couch

With Adjustable Back & Antimicrobial Upholstery



420 lbs.
weight
capacity

Upholstery Color



Midnight Express

Stability




Wooden Frame



*High-Quality Plywood
and Galvanized Steel*

Upholstery




2-inch Thick Foam



Antimicrobial Leather

Adjustability



Adjustable Back Rest

Add-On

ADI996-11-PA



Paper Roll Dispenser

ADI903-20-PC



Paper Roll Cutter

Specifications

Weight	136 lbs.
Dimensions	72" L x 27" W x 21.5" H
Weight Capacity	420 lbs.
Materials	Particleboard, Foam, PVC, Plywood, Laminate, and Rubber
Thickness of Foam Padding	2"
Adjustable Positions	27°, 36°, 43°, 50°, 55° & 60°
Upholstery	Antimicrobial
Color	Midnight Express
Storage Space Dimensions	63" L x 24" W x 9.5" H