



Xro Series Power Hi-Lo Exam Table Open Base, Flat Top Changing Table

ADI903-19-BLK



The Xro Series Power Hi-Lo Open Base Changing Table with Flat Top provides a reliable, comfortable, and efficient solution for clinical and community care environments. Designed to support a wide range of examinations and procedures, this foot power table offers smooth height adjustability from 17 to 33 inches, allowing healthcare professionals to position patients safely and ergonomically. Its flat-top surface and black medical-grade upholstery ensure easy sanitation, durability, and patient comfort, making it ideal for general exams, therapy, and assisted care applications. Built to meet the demands of high-use medical environments, the table lifts up to 440 lbs. and supports a working capacity of 550 lbs., powered by a single robust motor for quiet and consistent performance.

Built for performance and adaptability, the Xro Series table operates quietly on a standard 110V plug-in power source. Its open base allows easy access for staff and mobility aids, including Hoyer Lifts, enhancing safety and workflow efficiency in patient handling. Optional features, including guardrails, height risers, a paper roll holder, paper cutter, and hand control, allow customization for patient stability, increased table clearance, and improved workflow, making this table a dependable choice for medical offices, clinics, and community wellness facilities.

Features

Convenient: Convenient hands-free foot control or optional hand control allows for effortless height adjustment and seamless operation.

Functional: The smooth height adjustment allows for multiple positioning options for patient comfort and clinician access, safe and secure support for a variety of procedures, easy patient transfers, mobility with lockable wheels, and hygienic, easy-to-clean surfaces for efficient clinical use.

Upholstery: The durable, easy-to-clean black vinyl upholstery is ideal for sterile, high-traffic environments.

Design: The open base supports safe patient transfers and provides full access for caregivers or practitioners.

Optional Guardrails: Powder-coated steel guardrails (39" L) enhance safety during transfers and procedures.

Optional Height Risers: Steel plate risers increase table elevation, allowing higher ground clearance for Hoyer Lift compatibility.

Hygienic: An optional paper roll holder and dispenser provides quick, hygienic, and efficient patient preparation between appointments.

Application: Suitable for use in medical offices, outpatient clinics, physical therapy centers, chiropractic practices, women's health centers, urgent care facilities, hospitals and rehabilitation environments.

Robust: Powered by one motor, the table supports a weight capacity of 550 lbs., making it reliable for a wide range of patients.

Safe: The design incorporates sleek, rounded frame legs to minimize injury risks with 4 fixed knobs and 4 lockable casters for effortless mobility and secure stabilization during use.

Durable: The upholstery is crafted from antimicrobial PVC vinyl, which resists stains, mildew, oil, and wear, ensuring longevity and easy maintenance.

Features



**Black
PVC Leather Vinyl**

Flat-Top Surface

Up to
550 lbs.
weight
capacity



Adjustable Height

Hi-Lo function, operated via a hands-free foot control or optional hand control allows



Open-Base

Medical-Grade Locking Casters
Equipped with four casters for effortless mobility



Foot Pedal

Add-ons



Hand Controller

#903-20HC

Adjusts the table height effortlessly



Guard Rails

#ADI903-19-GR-BLK

Enhance safety during transfers and procedures



Extra Height Riser

#ADI903-19-EH

Increases table elevation, allowing higher ground clearance for Hoyer Lift compatibility



Paper Roll Dispenser

#903-20RH

Provides hygienic, and efficient patient preparation between appointments



Paper Cutter

#903-20PC

Ensures fast liner changes

Specifications

Dimensions	17" – 33" H x 27" W x 71" L
Height Range	17" – 33"
Materials	PVC Leather, Metal, Plastic, Plywood, Stainless Steel
Weight	142 lbs.
Upholstery Color	Black
Lifting Weight Capacity	440 lbs.
Working Weight Capacity	550 lbs.



Compliant

Meets ADA standards, ensuring accessibility for all patients.



Add-ons

Guardrail	
Weight	26 lbs.
Dimensions	1.2" H x 7" W x 39" L
Material	Powder-Coated Steel

Extra Height Riser	
Weight	15 lbs.
Dimensions	8" H x 2.3" W x 21.5" L
Material	Steel Plate
Weight Capacity	154 lbs. per Height Riser

Paper Roll Dispenser	
Weight	4.4 lbs.
Dimensions	6.9" H x 2" W x 25" L
Material	Powder-Coated Metal