



## Xro Open Base Series Power Hi-Lo Exam Tables



ADI903-20   ADI903-21   ADI903-22   ADI903-23   ADI903-24



### **ADI903-20-SERIES**

**Xro Series  
Power Hi-Lo Exam Table**  
Open Base With Adjustable Back &  
Antimicrobial Upholstery



### **ADI903-21-SERIES**

**Xro Open Base Series  
Power Hi-Lo Exam Table**  
With 3-Section Top  
(Adjustable Back, Leg-Rest &  
Antimicrobial Upholstery)



### **ADI903-22-SERIES**

**Xro Open Base Series  
Power Hi-Lo Exam Table**  
With 3-Section Top  
(Adjustable Back, Leg-Rest, Stirrups  
& Antimicrobial Upholstery)



### **ADI903-23-SERIES**

**Xro Shrouded Base Series  
Power Hi-Lo Exam Table**  
With 3-Section Top  
(Adjustable Back, Leg-Rest &  
Antimicrobial Upholstery)



### **ADI903-24-SERIES**

**Xro Shrouded Base Series  
Power Hi-Lo Exam Table**  
With 3-Section Top  
(Adjustable Back, Leg-Rest, Stirrups  
& Antimicrobial Upholstery)

## Xro Series Power Hi-Lo Exam Table

Open Base With Adjustable Back &  
Antimicrobial Upholstery

- Electric height adjustment from 17" to 33.5" inches
- Pneumatic adjustable backrest reclines from 0° to 65° with the support of a reliable gas strut
- Open-base scissor-lift design maximizes floor clearance and delivers stable, vertical movement
- Weight capacity up to 600 lbs.
- Powder-coated heavy metal frame
- 2" high-density foam and antimicrobial PVC vinyl upholstery
- Hands-free foot control allows for effortless height adjustment
- Optional handheld remote control



### Upholstery Colors

	<b>Midnight Express</b> #ADI903-20-BLU		<b>Black</b> #ADI903-20-BLK
---	---	---	--------------------------------

## Features



*Scissor-lift design*



*Hands-free height adjustment*  
With the use of the foot pedal



*Adjustable backrest*



*Medical-grade wheels*  
With locking mechanisms

## Accessories



**Hand Controller**  
#903-20HC



**Paper Roll Holder**  
#903-20RH



**Paper Cutter**  
#903-20PC

### Height

Maximum  
33.5"

Minimum  
17"



## Specifications

<b>Materials</b>	Steel Frame, Sponge, Leather, ABS Caster
<b>Overall Dimensions</b>	72" L x 29" W x 17-33.5" H
<b>Lifting Capacity</b>	400 lbs.
<b>Working Capacity</b>	600 lbs.
<b>Upholstery</b>	Antimicrobial PVC Vinyl
<b>Padding thickness</b>	2"
<b>Motors</b>	1 Dewert motor

## Xro Open Base Series Power Hi-Lo Exam Table

With 3-Section Top (Adjustable Back,  
Leg-Rest & Antimicrobial Upholstery)

- Electric height adjustment from 17" to 33.5" inches
- Pneumatic adjustable backrest reclines from 0° to 65° with the support of a reliable gas strut
- Manual tilting leg-rest section that adjusts from 0° to 80° upright
- Open-base scissor-lift design maximizes floor clearance and delivers stable, vertical movement
- Weight capacity up to 600 lbs.
- Powder-coated heavy metal frame
- 2" high-density foam and antimicrobial PVC vinyl upholstery
- Hands-free foot control allows for effortless height adjustment



### Upholstery Colors

	<b>Midnight Express</b> #ADI903-21-BLU		<b>Black</b> #ADI903-21-BLK
---	---	---	--------------------------------

## Features



*Adjustable leg rest*  
With removable cushion



*Hands-free height adjustment*  
With the use of the foot pedal



*Adjustable backrest*



*Medical-grade wheels*  
With locking mechanisms

## Accessories



**Hand Controller**  
#903-20HC



**Paper Roll Holder**  
#903-20RH



**Paper Cutter**  
#903-20PC

## Specifications

<b>Materials</b>	Steel Frame, Sponge, Leather, ABS Caster
<b>Overall Dimensions</b>	75.5" L x 29" W x 17"-33.5" H
<b>Lifting Capacity</b>	400 lbs.
<b>Working Capacity</b>	600 lbs.
<b>Upholstery</b>	Antimicrobial PVC Vinyl
<b>Padding thickness</b>	2"
<b>Motors</b>	1 Dewert motor





## Xro Open Base Series Power Hi-Lo Exam Table

With 3-Section Top (Adjustable Back, Leg-Rest, Stirrups & Antimicrobial Upholstery)

- Electric height adjustment from 17" to 33.5" inches
- Pneumatic adjustable backrest reclines from 0° to 65° with the support of a reliable gas strut
- Manual tilting leg-rest section that adjusts from 0° to 80° upright
- Open-base scissor-lift design maximizes floor clearance and delivers stable, vertical movement
- Weight capacity up to 600 lbs.
- Stirrups provide lower leg support and rotates from 0° to 115°
- Hands-free foot control allows for effortless height adjustment
- Medical-grade wheels with locking mechanisms



### Upholstery Colors

	<b>Midnight Express</b> #ADI903-22-BLU		<b>Black</b> #ADI903-22-BLK
---	---	---	--------------------------------

## Features



*Small Size Leg Cushion*  
#903-20SC



*Hands-free height adjustment*  
With the use of the foot pedal



*Adjustable backrest*



*Stirrups*  
#903-20ST

## Accessories



*Hand Controller*  
#903-20HC



*Paper Roll Holder*  
#903-20RH



*Paper Cutter*  
#903-20PC

## Specifications



<b>Materials</b>	Steel Frame, Sponge, Leather, ABS Caster
<b>Overall Dimensions</b>	75.5" L × 29" W × 17"-33.5" H
<b>Lifting Capacity</b>	400 lbs.
<b>Working Capacity</b>	600 lbs.
<b>Upholstery</b>	Antimicrobial PVC Vinyl
<b>Padding thickness</b>	2"
<b>Motors</b>	1 Dewert motor

## Xro Shrouded Base Series Power Hi-Lo Exam Table

With 3-Section Top (Adjustable Back,  
Leg-Rest & Antimicrobial Upholstery)

- Electric height adjustment from 21" to 32" inches
- Pneumatic adjustable backrest reclines from 0° to 65° with the support of a reliable gas strut
- Manual tilting leg rest section that adjusts from 0° to 80° upright
- Weight capacity up to 600 lbs.
- Powder-coated heavy metal frame
- Easy-clean ABS plastic shroud enclosed base
- 2" high-density foam and antimicrobial PVC vinyl upholstery
- Hands-free foot control allows for effortless height adjustment



### Upholstery Colors

**Midnight Express**  
#ADI903-23-BLU

**Black**  
#ADI903-23-BLK

## Features



*Shrouded base*



*Hands-free height adjustment*  
With the use of the foot pedal



*Adjustable backrest*



*Adjustable leg rest*  
With removable cushion

## Accessories



*Paper Roll Holder*  
#903-20RH



*Paper Cutter*  
#903-20PC

### Height

Maximum  
32"

Minimum  
21"



## Specifications

<b>Materials</b>	Steel Frame, Sponge, Leather, ABS Caster
<b>Overall Dimensions</b>	75.5" L × 29" W × 21"-32" H
<b>Lifting Capacity</b>	400 lbs.
<b>Working Capacity</b>	600 lbs.
<b>Upholstery</b>	Antimicrobial PVC Vinyl
<b>Padding thickness</b>	2"
<b>Motors</b>	1 Dewert motor

## Xro Shrouded Base Series Power Hi-Lo Exam Table

With 3-Section Top (Adjustable Back, Leg-Rest, Stirrups & Antimicrobial Upholstery)

- Electric height adjustment from 21" to 32" inches
- Pneumatic adjustable backrest reclines from 0° to 65° with the support of a reliable gas strut
- Manual tilting leg rest section that adjusts from 0° to 80° upright
- Weight capacity up to 600 lbs.
- Powder-coated heavy metal frame
- Easy-clean ABS plastic shroud enclosed base
- 2" high-density foam and antimicrobial PVC vinyl upholstery
- Stirrups provide lower leg support and rotates from 0° to 115°



### Upholstery Colors

	<b>Midnight Express</b> #ADI903-24-BLU		<b>Black</b> #ADI903-24-BLK
---	---	---	--------------------------------

## Features



*Small Size Leg Cushion*  
#903-20SC



*Hands-free height adjustment*



*Adjustable backrest*



*Stirrups*  
#903-20ST

## Accessories



*Paper Roll Holder*  
#903-20RH



*Paper Cutter*  
#903-20PC

### Height

Maximum  
32"

Minimum  
21"



## Specifications

<b>Materials</b>	Steel Frame, Sponge, Leather, ABS Caster
<b>Overall Dimensions</b>	75.5" L × 29" W × 21"-32" H
<b>Lifting Capacity</b>	400 lbs.
<b>Working Capacity</b>	600 lbs.
<b>Upholstery</b>	Antimicrobial PVC Vinyl
<b>Padding thickness</b>	2"
<b>Motors</b>	1 Dewert motor