



# HEIGHT ADJUSTABLE EXAM STOOL

- ADI910-01-BLK   ADI910-02-BLK   ADI910-03-BLK
- ADI910-01-BLU   ADI910-02-BLU   ADI910-03-BLU

Step into unparalleled comfort with our line of Height Adjustable Exam Stools, a definitive choice for those seeking an exceptional seating solution. We've meticulously engineered our stool to blend functionality with captivating style to redefine healthcare seating. The body of the stool is constructed from durable stainless steel to ensure your investment will last for many years. Each of our stool bases is outfitted with premium, easy-to-clean antimicrobial upholstery that not only assures luxurious comfort and appeal but creates a hygienic environment by minimizing the risk of cross-contamination. For optimal plushness, our stool base has 4" foam that delivers consistent relief.

Additionally, we've equipped them with 2" caster wheels that supply straightforward maneuverability. You can seamlessly glide across your floors without leaving any marks behind. In order to deliver an adaptable stool, we've equipped our Height Adjustable Exam Stools with the ability to raise or lower to your desired height. This will enable you to maintain proper posture and help cater your seating position depending on the medical procedures or examinations being conducted. Though this may seem like an average piece of medical furniture, it's deceptively crucial. It will provide long-term comfort, hold up to daily wear and tear, and enhance your office's overall aesthetic.



ADI910-01-BLU



ADI910-02-BLK



ADI910-03-BLU



## FEATURES



### UPHOLSTERY COLOR OPTIONS



Black



Midnight  
Express



*Iron & Aluminum  
Material*



*Smooth Gliding  
Caster Wheels*

## SPECIFICATIONS

Pneumatic Height Adjustable Exam Stool	
Materials	Stainless Steel, Foam, Antimicrobial PU Leather and Plywood
Overall Dimensions	18.8" – 26" H x 14.3" W x 14" D
Weight	12.1 lb.
Cushion Thickness	4.1"
Seat Diameter	14"
Wheels	4
Weight Capacity	250 lb.
Height Adjustment	Pneumatic Gas Lift Seat Height: 18.8" - 26"
Wheels Diameter	2"
Color Options	Black and Midnight Express

## FEATURES



### UPHOLSTERY COLOR OPTIONS




Black



Midnight  
Express



*Smooth Gliding  
Caster Wheels*



*Circular  
Foot Ring*

## SPECIFICATIONS

Pneumatic Height Adjustable Exam Stool With Backrest	
Materials	Stainless Steel, Foam, Antimicrobial PU Leather and Plywood
Overall Dimensions	28" - 32.6" H x 16.1 D
Weight	16.5 lb.
Cushion Thickness	4.1"
Seat Diameter	16.1"
Wheels	5
Weight Capacity	250 lb.
Height Adjustment	Pneumatic Gas Lift Seat Height: 18.8" - 23.6"
Wheels Diameter	2"
Color Options	Black and Midnight Express
Backrest	8.2" H x 14.9" W

## FEATURES



### UPHOLSTERY COLOR OPTIONS



Black



Midnight  
Express



*Contoured  
Backrest*



*Pneumatic Lift*

## SPECIFICATIONS

Height Adjustable Exam Stool With Circular Foot Ring	
Materials	Stainless Steel, Foam, Antimicrobial PU Leather and Plywood
Overall Dimensions	19" H x 23" W x 16.1" D
Weight	12.1 lb.
Cushion Thickness	4.1"
Seat Diameter	16.1"
Wheels	4
Weight Capacity	250 lb.
Height Adjustment	Seat Height: 19.1"-24.6"
Wheels Diameter	2"
Color Options	Black and Midnight Express
Circular Foot Ring	9.8"