



Height Adjustable Exam Stool



ADI910-01-BLK ADI910-02-BLK ADI910-03-BLK
ADI910-01-BLU ADI910-02-BLU ADI910-03-BLU

Step into unparalleled comfort with our line of Height Adjustable Exam Stools, a definitive choice for those seeking an exceptional seating solution. We've meticulously engineered our stool to blend functionality with captivating style to redefine healthcare seating. The body of the stool is constructed from durable stainless steel to ensure your investment will last for many years. Each of our stool bases is outfitted with premium, easy-to-clean antimicrobial upholstery that not only assures luxurious comfort and appeal but creates a hygienic environment by minimizing the risk of cross-contamination. For optimal plushness, our stool base has 4" foam that delivers consistent relief.

Additionally, we've equipped them with 2" caster wheels that supply straightforward maneuverability. You can seamlessly glide across your floors without leaving any marks behind. In order to deliver an adaptable stool, we've equipped our Height Adjustable Exam Stools with the ability to raise or lower to your desired height. This will enable you to maintain proper posture and help cater your seating position depending on the medical procedures or examinations being conducted. Though this may seem like an average piece of medical furniture, it's deceptively crucial. It will provide long-term comfort, hold up to daily wear and tear, and enhance your office's overall aesthetic.



ADI910-01-BLU



ADI910-02-BLK



ADI910-03-BLU





250 lbs.
Weight
Capacity

Upholstery Color Options



Black



Midnight
Express

Features



Steel



Smooth Gliding
Caster Wheels


Specifications

ADI910-01-BLK ADI910-01-BLU		Pneumatic Height Adjustable Lux Exam Stool
Materials	Steel, Foam, Antimicrobial PVC Upholstery	
Overall Dimensions	18.3" – 23.6" H x 22" D	
Weight	12.5 lbs.	
Cushion Thickness	4"	
Seat Diameter	16"	
Qty of Wheels	5	
Wheels Diameter	2"	
Height Adjustment	Seat Height: 18.3" - 23.6"	
Weight Capacity	250 lbs.	
Color Options	Black and Midnight Express	



300 lbs.
Weight
Capacity

Upholstery Color Options



Black



Midnight
Express

Features



Smooth Gliding
Caster Wheels



Circular
Foot Ring

Specifications

ADI910-02-BLK ADI910-02-BLU		Height Adjustable Exam Stool With Circular Foot Ring
Materials	Steel, Foam, Antimicrobial PVC Upholstery	
Overall Dimensions	18.5"- 25" H x 14.9" D	
Weight	12 lbs.	
Cushion Thickness	4"	
Seat Diameter	14.17"	
Qty of Wheels	4	
Wheels Diameter	2"	
Height Adjustment	Pneumatic Gas Lift Seat Height: 18.5" - 25"	
Weight Capacity	300 lbs.	
Color Options	Black and Midnight Express	



250 lbs.
Weight
Capacity

Upholstery Color Options



Black



Midnight
Express

Features



Contoured
Backrest



Pneumatic Lift

Specifications

ADI910-03-BLK ADI910-03-BLU		Pneumatic Height Adjustable Exam Stool With Backrest
Materials	Steel, Foam, Antimicrobial PVC Upholstery	
Overall Dimensions	18.3" - 23.6" H x 16.2" D	
Weight	16 lbs.	
Cushion Thickness	4"	
Seat Diameter	16.1"	
Qty of Wheels	5	
Wheels Diameter	2"	
Height Adjustment	Seat Height: 18.3"-23.6"	
Weight Capacity	250 lbs.	
Color Options	Black and Midnight Express	
Backrest	8.2" H x 14.9" W	