



# TREATMENT TABLE

with Shelf and Adjustable Back



ADI996-02-BLK



Treatment Table with Shelf and Adjustable Back, Black

Our foldable Treatment Table with Shelf and Adjustable Back is an indispensable asset for any medical facility, engineered to meet the rigorous demands of modern healthcare environments. Its stainless-steel frame offers exceptional durability and easy cleaning, which is ideal for high-traffic medical environments. The adjustable backrest, with multiple recline angles ranging from 0° to 67°, ensures patient comfort during various treatments. When not in use, the foldable leg design and less than 75 lb. body allow for efficient space-saving storage. Our antimicrobial mattress delivers an extra layer of hygiene, which is crucial for maintaining a sterile treatment environment. The integrated paper roll dispenser reassures a sanitary space for your patients. The outfitted shelf will help provide convenience for medical supplies, as well as spare exam paper rolls.

This table is practical and patient-friendly, built to endure the rigors of daily clinical use. With rubber feet to prevent floor scratches, it offers stability and protection in any setting. Measuring 31.5 inches in height, 27.5 inches in width, and 74.8 inches in length, the table delivers ample surface area for various treatments without compromising maneuverability. Whether used in hospitals, rehabilitation centers, or private practices, this table ensures superior performance with its high-quality materials and meticulous design—an essential tool for professionals who demand reliability and efficiency.

## FEATURES

**Sturdy:** Constructed from stainless steel, this table provides exceptional durability and is easy to clean, making it ideal for high-traffic medical environments.

**Adjustable:** The backrest adjusts to six angles—0°, 32.5°, 45°, 54.5°, 62°, and 67°—ensuring versatile patient positioning and enhanced comfort.

**Space-Saving:** Foldable legs allow for easy storage, maximizing available space when the table is not in use.

**Storage:** Includes a particle board shelf beneath the table for quickly organizing and accessing essential medical supplies.

**Sanitation:** Integrated paper roll holder simplifies maintaining a clean and hygienic surface between patients.

**Sterile:** Features a mattress treated with antimicrobial agents, ensuring a hygienic surface for patient care.

**Protective:** The rubber feet help prevent scratches and damage to the floors, ensuring stability and protection during use.

# SPECIFICATIONS

ADI996-02-BLK	
Weight	73.63 lbs.
Dimensions	75.6" L x 27.6" W x 31.5" H
Weight Capacity	440 lbs
Backrest Adjustable	0°, 32.5°, 45°, 54.5°, 62°, and 67°
Color	Black Upholstery
Mattress Thickness	4"
Materials	201 Stainless Steel, Foam, PU Leather, and Particle Board



## FEATURES



### SPACE-SAVING

FOLDABLE LEGS ALLOW FOR EASY STORAGE, MAXIMIZING AVAILABLE SPACE WHEN THE TABLE IS NOT IN USE.



### STORAGE

INCLUDES A PARTICLE BOARD SHELF BENEATH THE TABLE FOR QUICKLY ORGANIZING AND ACCESSING ESSENTIAL MEDICAL SUPPLIES.



### ADJUSTABLE

THE BACKREST ADJUSTS TO SIX ANGLES—0°, 32.5°, 45°, 54.5°, 62°, AND 67°—ENSURING VERSATILE PATIENT POSITIONING AND ENHANCED COMFORT.



### SANITATION


INTEGRATED PAPER ROLL HOLDER SIMPLIFIES MAINTAINING A CLEAN AND HYGIENIC SURFACE BETWEEN PATIENTS.



*Adjustable Backrest*



*Integrated Paper Roll Holder*



*Foldable Legs*



*Shelf*