

AdirMed 996-01-BLK

Adjustable Treatment Table

Product Highlights

- **440 lbs** Weight Capacity & Incredibly Light-Weight Construction
- **Adjustable** Backrest with Multiple Positioning
- **2" Thick** Foam Padding for Patient Comfort
- **Seamless** Reinforced Upholstery & Safe Rounded Corners
- **Paper** Dispenser with Easy Access
- **2-Year** Warranty

Physical Specifications

Dimensions: 74.8"L x 27.5"W x 31.4"H

Weight: 60 lbs.

Accessories: Paper Dispenser

Board Material: Steel

Upholstery Material: Knit-Backed Upholstery



AdirMed 996-03-MA

Wood Exam Table with Full Shelf

Product Highlights

- **400 lbs** Weight Capacity
- **Aesthetically-Pleasing** and Feature-Rich Exam Table
- **Designed** for Optimal Patient Comfort and Ergonomics
- **Constructed** from Acacia Wood for durability
- **32-inch** Working Height with Backrest
- **Smooth** and Easy-to-Clean Vinyl Surface
- **2-Year** Warranty

Physical Specifications

Dimensions: 74" D x 28" W x 32" H

Weight: 115 lbs.

Accessories: Paper Dispenser

Board Material: Acacia Wood

Upholstery Material: Vinyl Surface

