



# Upholstered Therapy Rehab Mat Table

ADI996-05-BLK



Transform your rehabilitation experience with our oversized Upholstered Therapy Rehab Mat Table. This specialized table is meticulously designed. We've crafted our table with a 2-inch-thick mattress to ensure your patients have a plush, relaxing surface ideal for extended use. The antimicrobial upholstery provides a hygienic environment and safeguards against potential contaminants, prioritizing cleanliness and safety. Additionally, it can handle being regularly sanitized in between patients. With a lengthy size of 72" L x 53" W, this therapy table is optimized for versatility and functionality, accommodating a range of body types and therapy requirements. The sturdy 6-leg design with center supports contributes to its strength and allows it to handle daily rehabilitation sessions. Its thoughtfully engineered leg design will ensure robust support, granting stability and durability for various therapy routines. The wooden legs and chic black mattress enhance its aesthetic appeal, seamlessly blending into any clinic or rehabilitation space. Elevate your space today with our premium Upholstered Therapy Rehab Mat Table, offering comfort and reliability for a holistic rehabilitation experience.

## FEATURES

**Sanitation:** Our antimicrobial mattress is designed to provide a safe and hygienic environment for patients by reducing the risk of cross-contamination.

**Sturdy:** This table is engineered with a sturdy, solid wood 6-leg design and center supports, providing exceptional strength and durability.

**Comfort:** The 2-inch-thick mattress provides optimal extended patient comfort.

**Aesthetic:** The chic black mattress and wooden legs enhance the table's visual appeal while seamlessly blending into any clinic or rehabilitation space.

**Robust:** The generous size will give you and your patient ample space to conduct therapy sessions.

**Threshold:** Its 500-pound weight capacity accommodates individuals of varying sizes while providing support during therapeutic sessions.

# SPECIFICATIONS

ADI996-05-BLK	
Weight	92.5 lb.
Dimensions	72" L x 53" W x 18.1" H
Padding	2" thick high-density foam
Color of Upholstery	Black
Upholstery	Antimicrobial
Weight Capacity	500 lb.

## FEATURES



### CAPACITY

ENGINEERED TO ACCOMMODATE UP TO 500 LB.



### AESTHETIC

CHIC MATTRESS AND WOODEN LEGS ENHANCE ITS VISUAL APPEAL



### COMFORT

2-INCH-THICK MATTRESS PROVIDES OPTIMAL PATIENT COMFORT



### ROBUST

AMPLE SPACE TO CONDUCT THERAPY SESSIONS



*Anti-Microbial Mattress*



*Solid Wood Legs*



*Two Center Supports*