



Wall-Mounted Folding Therapy Table



SKU: ADI996-06-BLK



Maximize space without compromising functionality with our **Wall-Mounted Folding Therapy Table**. Built for stability, this table has a three-front-leg design that locks securely in place, providing a solid foundation for treatments. The generously sized surface, layered with 2-inch foam padding, offers optimal patient support and comfort. The antimicrobial black PU leather upholstery promotes hygiene in high-traffic medical environments. The expansive surface area provides ample room for a therapist to sit or lean alongside the patient, facilitating techniques that require close interaction, such as assisted range-of-motion exercises, soft tissue manipulation, or postural training. Ideal for physical therapy clinics, rehabilitation centers, hospitals, sports medicine facilities, and wellness studios, this space-saving solution seamlessly integrates into any professional setting. When not in use, the foldable design allows for vertical storage, minimizing floor space usage—perfect for small treatment rooms, multi-purpose spaces, and home therapy setups.

Features

Stable: The three-front-leg design locks securely, ensuring a firm, balanced surface.

Compact: The wall-mounted foldable design allows easy storage in compact treatment areas.

Comfortable: 2-inch mattress padding provides superior support for extended sessions.


Hygienic: Its antimicrobial black PU leather upholstery helps maintain a clean and sanitary environment.

Protective: An edge bumper reinforces the table's perimeter, preventing wear and damage.

Robust: With a 500 lbs. weight capacity, this table supports many therapeutic applications.

Spacious: A large treatment area accommodates both patients and practitioners for hands-on therapy techniques.

Features



Black PU Leather Upholstery



Comfortable
2-inch Mattress Padding



Sanitation
Antimicrobial Mattress



Protective
Edge Bumper



Space-Saving
Wall-Mounted Foldable Design



Mounting Block and Chain



Folded

Depth of 6.3" When Folded



Stable
Three-Front-Leg Design

Foldable Bed and Legs




500 lbs.
Weight Capacity



Spacious

Large treatment area accommodates both patients and practitioners



Specifications

Materials	Rubber, MDF, plywood, PU leather
Color	Black upholstery and natural wood tone base
Weight Capacity	500 lbs.
Net Weight	154 lbs.
Overall Dimensions (Open)	84" L x 51" W x 17" H
Folded Height	62.5"
Surface Padding	2-inches Mattress
Edge Protection	Preinstalled Bumper
Leg Design	Three-Front-Leg Structure with Locking Mechanism
Storage Mechanism	Wall-Mounted, Mounting Block and Chain
Depth On Wall When Folded	6.3"