



PEDIATRIC GRANDE STEEL EXAM TABLE



[ADI996-09-BASE](#)

[ADI996-09-PMB-BLK](#)

[ADI996-09-PM-BLK](#)

[ADI996-09-PMB-BLU](#)

[ADI996-09-PM-BLU](#)



**Pediatric Grande Steel Exam Table,
Midnight Express**

When outfitting your pediatric exam room, ensuring safety and comfort for your smallest patients is paramount. The Pediatric Grande Steel Exam Table, designed with child-friendly features, offers the ideal solution for pediatric healthcare settings. Its robust steel body will continuously provide long-lasting durability and is effortless to clean between patients. The dual-sided cabinet ensures easy access to storage from either side of the table, while two Euro-style drawers with stainless steel slides offer ample space for essential items. An outfitted pull-out step stool with a safety tread provides secure access for young patients. At the same time, the cushioned, contoured mattress conforms to the children's bodies for ultimate comfort and safety. The antimicrobial material reduces the risk of cross-contamination, ensuring a safe environment for every patient.

Additionally, this table offers a pull-out leg rest for accommodating different patient sizes. The built-in paper dispenser, compatible with 21.8" paper rolls, also makes changing liners between patients seamless. With its thoughtful, baby-centric design and high functionality, the Pediatric Grande Series Steel Exam Table gives healthcare professionals the peace of mind they need while caring for their young patients.

FEATURES

Contoured: The curved edges and premium mattress design conform to the child's natural shape, ensuring a safe and gentle examination experience.

Safety: The soft, curved edges help prevent children from rolling off the table during the examination.

Dual-Sided Access: Storage from both sides of the table allows flexibility in placement and accessibility.

Soft Closing: Both the drawers and cabinet feature soft-closing mechanisms to minimize noise and enhance safety.

Step Stool: A pull-out step stool with safety tread allows young patients to access the table easily, and it pushes back smoothly when not in use.

Integrated: A built-in dispenser with storage for extra rolls ensures fast liner changes.

Antimicrobial Protection: The antimicrobial mattress helps maintain a hygienic space, reducing the risk of cross-contamination.

SPECIFICATIONS

Pediatric Grande Steel Exam Table	
Overall Size	50.98" L x 27.17" W x 35.72" H
Base Weight	168.65 lb.
Base Dimensions	45.4" L x 26.77" W x 29.13" H
Cabinet Dimensions	22.44" L x 17.32" H x 22.83" D
Upper Drawer Dimensions	Internal Dimensions: 4.7" H x 15.7" D - Overall Width: 22"
Lower Drawer Dimensions	Internal Dimensions: 5.5" H x 15.7" D - Overall Width: 22"
Material	Steel
Finish	Powder Coating
Length with Leg Rest Extended	14.1"
Weight Capacity	500 lb.
Compatible Paper Roll Size	Length: 21.8" Height: 2.28"
Step Stool Dimensions	13" L x 22" W x 8.3" H
Ground To Top of Step Stool	8.3"
Certification	California CARB II Compliant
Mattress Weight	27.60 lb.
Mattress	2-inch padded foam antimicrobial PU leather upholstery
Mattress Dimensions	50.98" L x 27.17" W x 4.72" H
Full Length with Pullout Shelf Mattress	67.3" L
Mattress Color Options	Black, and Midnight Express
Total Weight	196.25 lb.

FEATURES



Step Stool



Two Drawers



Dual-Sided Access



Pull-out Leg Rest

UPHOLSTERY COLOR OPTIONS



Black



Midnight Express

