



ARISTO TREATMENT TABLE with Tilted Mattress



ADI996-HBTDP-BLK



Aristo H-Braced Treatment Table with Tilted Mattress, Shelf, Drawers and Paper Roll Dispenser – Black

Indulge in comfort and efficiency with our Aristo Treatment Table with Tilted Mattress. Its adjustable backrest offers unparalleled versatility, allowing customization to six angles: 27, 36, 43, 50, 55, and 60 degrees, as well as adjusting to a completely flat surface. This ensures optimal positioning for various treatments and examinations. Its sturdy wooden base, coated with a natural clear PU paint, not only exudes modernity but also ensures long-lasting durability. Patients will experience sublime relaxation on the 2-inch-thick antimicrobial mattress with rounded corners, offering comfort and hygiene in equal measure. The antimicrobial upholstery will ensure sterilization between patients and can be easily sanitized. These exam tables are engineered to accommodate up to 400 lbs. and fortified by a support brace for added stability, instilling confidence in their reliability.

To allow for customization, optional accessories such as our full-length shelf to store medical supplies, two drawers that can be opened from both sides, and a paper roll dispenser can be added (available separately) that will fully optimize your exam table. You can further customize the therapy table by selecting from one of our upholstery colors to match your aesthetic: Black or Midnight Express. Either choice will complement your workspace seamlessly and will never become outdated. Choose our Aristo Treatment Table with Tilted Mattress that redefines excellence today, ensuring both you and your clients experience the epitome of comfort and style.

FEATURES

Sanitation: The antimicrobial mattress is designed to provide a safe and hygienic environment for patients by reducing the risk of cross-contamination.

Adjustable: Customize the angle of the backrest to six positions - 27, 36, 43, 50, 55, and 60 degrees, as well as entirely flat.

High-Quality: The treatment table is constructed with a sturdy wooden base coated with a natural clear PU paint, ensuring its durability.

Plush: Patients will experience relaxation on the 2-inch-thick mattress with rounded corners.

Capacity: Engineered to accommodate up to 400 lbs., ensuring stability and confidence in the reliability of the table.

Reinforced: The exam table is fortified with a support brace for added stability, enhancing its overall strength.

Customization: Optional accessories such as a full-length shelf, two drawers that can be opened from both sides, and a paper roll dispenser and cutter (available separately) can be added for personalized functionality.

ADI996-HBTDP-BLU



Aristo H-Braced Treatment Table with Tilted Mattress, Shelf, Drawers and Paper Roll Dispenser—Blue

FEATURES



Tilted Mattress



H-Brace Base



Two Drawers



Shelf



Paper Roll Dispenser



Paper Cutter

UPHOLSTERY COLOR OPTIONS



Black



Midnight Express

SPECIFICATIONS

Aristo H-Braced Treatment Table with Tilted Mattress		
Upholstery Colors	Midnight Express	Black
Materials	PU Leather, Plywood for Mattress, and Solid Wood and Laminate for Support Brace	
Overall Dimensions	30.75" H x 28.03" W x 72.05" L	
Mattress Dimensions	Adjustable Backrest: 24.80" L x 28.03" W Other Half of Mattress: 47.24" L x 28.03" W	
Backrest Adjustability	27°, 36°, 43°, 50°, 55°, & 60°	
Net Weight	64.81 lb.	
Weight Capacity	400 lb.	
Padding	2" Thick High-Density Foam	
Upholstery	Antimicrobial	
Optional Accessories Dimensions	Drawers: 17.44" H x 21.65" W x 17.71" D	
	Internal Drawer Depth: 9.45"	
	Shelf: 3.15" H x 17.72" W x 65.74" L	
	Paper Roll Dispenser: 2.36" H x 21.6" W x 1.57" D	

KIT SKUS FOR ADI996-HBT-SERIES:

Aristo H-Braced Treatment Table with Tilted Mattress	Upholstery Color	Shelf	Drawers	Paper Roll Dispenser
ADI996-HBT-BLK	■			
ADI996-HBTD-BLK	■	✓	✓	
ADI996-HBTDP-BLK	■	✓	✓	✓
ADI996-HBTP-BLK	■			✓
ADI996-HBTS-BLK	■	✓		
ADI996-HBTSP-BLK	■	✓		✓
ADI996-HBT-BLU	■			
ADI996-HBTD-BLU	■	✓	✓	
ADI996-HBTDP-BLU	■	✓	✓	✓
ADI996-HBTP-BLU	■			✓
ADI996-HBTS-BLU	■	✓		
ADI996-HBTSP-BLU	■	✓		✓

STOCK SKUS FOR ARISTO H-BRACED TREATMENT TABLE :

ADI996-HB	H-Brace Base
ADI996-HB-DR	Drawers
ADI996-HB-PR	Paper Roll Dispenser + Cutter
ADI996-HB-SHELF	Shelf
ADI996-HB-TM-BLK	Tilted Mattress - Black
ADI996-HB-TM-BLU	Tilted Mattress - Blue