

ADIME904-06



The Remote Display Digital Scale seamlessly combines precision and versatility for any professional environment, featuring a comprehensive range of functions designed to streamline use and ensure accuracy. Capable of supporting up to 550 lbs., this scale accommodates a broad range of users, making it ideal for clinical or commercial applications. It is equipped with a BMI function, allowing users to effortlessly track body mass index, while the zero/hold function ensures reliable measurements, even in fluctuating conditions. The scale features a large, clear LCD display for easy readings, with the ability to lock units of measure and toggle between pounds and kilograms. The inclusion of mounting hardware adds flexibility, enabling the remote display to be installed for optimal convenience or left as a standalone unit with the included bracket.

This unit operates efficiently on three 1.5V batteries or via an AC adapter, offering multiple power options for extended use. Precision is paramount in any setting, and this scale excels with its anti-slip platform surface and external calibration for unparalleled accuracy. Weighing just 8.31 lbs. and constructed from durable carbon steel, it ensures longevity while maintaining a sleek profile with compact dimensions of 12" x 12.5". Whether mounted or used in its portable form, this Remote Display Digital Scale delivers reliable performance and robust functionality, making it an invaluable tool for any professional setting.

FEATURES

Capacity: The scale supports a maximum weight of 550 lbs., accommodating a wide range of users.

BMI Function: Equipped with a BMI function, the scale allows users to easily track body mass index for health assessments.

Zero/Hold: The zero/hold function provides reliable, consistent measurements, even in fluctuating conditions.

Display: A large, clear LCD display ensures easy reading, while also offering the ability to lock units of measure and toggle between pounds and kilograms.

Mounting Flexibility: The included mounting hardware enables the remote display to be installed for optimal convenience or used as a standalone unit with the included bracket.

Power Options: The scale operates efficiently on three 1.5V batteries or via an AC adapter, providing flexible power options for extended use.

FEATURES





CAPACITY
THE SCALE SUPPORTS A
MAXIMUM WEIGHT OF 550 LBS.



DURABILITY

CONSTRUCTED FROM DURABLE
CARBON STEEL, THIS SCALE
OFFERS BOTH LONGEVITY AND
A SLEEK, COMPACT DESIGN.



ACCURACY
THE SCALE ENSURES PRECISE
MEASUREMENTS IN ANY SETTING.



SAFETY
AN ANTI-SLIP PLATFORM
SURFACE, TO PREVENT THE
PATIENT FROM SLIPPING

SPECIFICATIONS

SKU	ADIME904-06
Scale Dimensions	2" H x 12" W x 12.5" L
Platform Dimensions	12" x 12.4"
Base Weight	8.31 lbs.
Indicator Weight	0.55 lbs.
Cord Weight	0.22 lbs.
Adaptor Weight	0.22 lbs.
Weight of All Components	9.30 lbs.
Included	Mounting Hardware and Bracket
Weight Capacity	550 lbs.
Resolution	0.2 lb. (0.1 kg)
Power Supply	AC Adapter (Included) and 3 1.5 V AA Battery Operation (Batteries Not Included).
Functions	BMI, ON, OFF, UP, DOWN, KB/LB., ZERO, and HOLD
Cable Length	2.62' (0.8 meter)
Display	Backlit LCD
Shut Off	Manual or Auto Shut-Off after 2 minutes of inactivity.

