



Digital Chair Scale

ADIME904-07



Our **Digital Chair Scale** is designed to deliver uncompromising accuracy, comfort, and convenience in demanding healthcare environments. With a weight capacity of up to 550 lbs. and a UL certified power supply, it meets the highest standards of safety and reliability. Engineered with both patients and caregivers in mind, this scale features a lightweight frame with foldable foot pedals that automatically position when the patient is seated, ensuring ease of use and reducing physical strain. The compact, foldable design allows for effortless mobility and space saving storage, making it ideal for hospitals, clinics, dialysis centers, nursing homes, and mobile healthcare services.

A backlit, rechargeable LCD display offers clear visibility and extended operation without the hassle of frequent battery changes. The intuitive interface and weight hold function allow clinicians to record measurements with confidence and without rushing. For added accuracy, an integrated bubble level at the rear of the chair ensures the scale is properly positioned, even on uneven flooring, delivering consistent and reliable readings.

Features

Strength: Capable of supporting patients up to 550 lb, this scale provides the reliability needed in healthcare environments.

Precision: Built-in bubble level ensures the scale is positioned on a level surface, promoting accurate weight readings.

Comfort: A supportive seat and backrest ensure patients remain stable and at ease throughout the weighing process.

Portability: A lightweight frame and folding design allow staff to move and store the scale with minimal effort.

Clarity: The large LCD display offers clear weight readings with options for tare, unit conversion, and BMI calculation.

Convenience: A weight hold function locks the measurement, allowing staff to focus on patients before recording results.

Safety: Locking rear wheels provide secure positioning during use, preventing unintended movement.

Efficiency: Quick tare and zeroing functions simplify patient transfers and ensure accurate results every time.

Durability: Built with a sturdy frame and UL certified plug, the unit maintains dependable performance under daily clinical use.

Collapsible: A handle release mechanism allows the handle to be quickly folded or unfolded, simplifying storage and enabling rapid deployment.

Features



Dimensions



Maximum Capacity
550 lbs.



No
Tare

Off

Unit

Hold

Set

Display

Tare Function

Simplify patient transfers and ensure accurate results every time.

Unit Selection

Switches the unit between kg. and lbs.

Hold Function

Locks the measurement, allowing staff to focus on patients before recording results.

Set Function

Enter patient height and identifying information.

Specs

Materials	Powder Coated Steel Frame, Foam,PVC,Wood, Plastic
Net Weight	44 lbs.
Overall Dimensions	39.2" H x 24.4" W x 40.2" D
Seat Dimensions	21.5" H x 20" W x 19" D
Footrest Dimensions (Folded)	26" H x 35.4" D
Footrest Dimensions	7.68" L x 5.7" W
Power Supply	Rechargeable battery with UL certified plug
Display	Large LCD screen with weight, tare, unit conversion, BMI, and hold functions
Mobility	Two directional front wheels and two swivel rear wheels with locking brakes
Weight Capacity	550 lbs.