



Digital And Beam PHYSICIAN SCALES



[ADIME904-01](#)

[ADIME904-02](#)

[ADIME904-03](#)

[ADIME904-04](#)

[ADIME904-05](#)

Discover our comprehensive range of precision scales catering to varied needs, from our Digital Scale, which offers accurate digital readings, to our Physician Beam Scale, which provides traditional functionality. Our selection also includes specialized scales such as our premium Veterinary Scale designed to accommodate animals of all sizes, our Wheelchair Scale that ensures convenience and accuracy for patients with mobility needs, and our Medical Handrail Scale, which delivers a stable, secure way to achieve precise patient measurements. Whatever your specific needs are, we will surely have a scale that delivers and will withstand the test of time.


FEATURES

Platform: Our adequately sized platforms give patients plenty of room to stand comfortably while measurements are being conducted.

Durability: Their hardy metal body ensures longevity from continual wear and tear.

User-Friendly: Each are desinged to be straightfoward and allow you to get instant, accurate readings.

Sanitary: All of our scales are constructed to be easy to clean so that you can quickly eliminate pathogenic agents after each use.



**Stainless Steel
Veterinary Scale**



**Digital Physician Scale
With Height Rod**



**Physician Beam Scale
With Height Rod**

SPECIFICATIONS

Product Number	Type	Weight	Overall Dimensions	Weight Capacity	Platform Dimensions
ADIME904-01	Digital	19.84 lb	37.40" H x 11.81" W x 20.47" D	660 lb	15.75" x 11.81"
ADIME904-02	Beam	35.27 lb	58.2" H x 19.29" W x 14.76" D	440 lb	14.76" x 10.82"
ADIME904-03	Digital	143.3 lb	32.56" H x 39.36" W x 1.96" D 47.63" L	1100 lb	47.63" x 39.76"
ADIME904-04	Digital	44.09 lb	44.09" L x 20.47" W x 1.77" D	660 lb	44.09" x 20.47"
ADIME904-05	Digital	55.11 lb	51.6" H x 20.47" W x 24.4" D	660 lb	24.41" x 20.47"

Digital Physician Scale With Height Rod

ADIME904-01



- ✓ LCD Screen Shows Patient's Weight, Indicate BMI, and Can Switch Units Between Lb and Kg
- ✓ Works With Two Power Sources: AC Adapter (Included) or 6V Rechargeable Batteries

Physician Beam Scale With Height Rod

ADIME904-02



- ✓ Die Casted Dual Beams With Two Easy Manual Sliding Weights Can Be Read From Either Side
- ✓ Two Easy Gliding Wheels So That You Can Maneuver The Scale Wherever Needed



Medical Wheelchair Scale

ADIME904-03



- ✓ Built-In Ramp On Both Sides and Stainless-Steel Handrails Allow Easy Access
- ✓ LCD Screen Shows Patient's Weight, Indicate BMI, and Can Switch Units Between Lb and Kg

Stainless Steel Veterinary Scale

ADIME904-04



- ✓ Spacious Platform Design Reduces Stress and Facilitates Easy Loading of Animals Without Strain or Injury
- ✓ LCD Screen Shows Patient's Weight, Indicate BMI, and Can Switch Units Between Lb and Kg

Medical Handrail Scale

ADIME904-05



- ✓ Sturdy Handrail Provides Enhanced Support For Individuals of Varying Mobility Levels
- ✓ LCD Screen Shows Patient's Weight, Indicate BMI, and Can Switch Units Between Lb and Kg

