



FOLDING PARALLEL BARS (7 FEET)

ADIME905-01



Step into a world of limitless possibilities for your physical therapy practice with our remarkable 7-foot Folding Parallel Bars. Elevate your rehabilitation sessions to new heights as these bars blend convenience and strength, providing an exceptional experience for therapists and patients alike. Crafted with precision, these width-adjustable Floor Folding Parallel Bars boast double stiffening rails, ensuring unmatched stability and support throughout each exercise. With easy height adjustment through our screw knob and fourteen height options, achieving the perfect therapeutic setting becomes an effortless endeavor.

We've built our parallel bars with heavy-duty square telescoping steel uprights to endure the demands of intense physical therapy use, guaranteeing utmost safety and ensuring it can handle up to 350 lbs. The second built-in height knob eliminates unwanted movement, fostering a smooth and controlled session. Additionally, each handrail is thoughtfully equipped with end bumpers, safeguarding the equipment and the surrounding environment. When you need extra space or the parallel bars aren't being used, you can easily fold them to a compact size and store them away. Embrace a world of limitless therapeutic possibilities with our 7-foot Folding Parallel Bars, engineered explicitly for physical therapy use.

FEATURES

Customization: The width can be easily adjusted to suit user preferences and exercise requirements.

Durable: 1-1/2" one-piece stainless-steel handrails give patients a reliable grip.

Adjustable: Our bars feature screw knobs and fourteen height options allowing for effortless customization of the height setting.

Heavy-Duty: The square telescoping steel uprights are built to endure the demands of intense physical therapy use.

Capacity: Our Folding Parallel Bars can handle up to 350 lb. accommodating various patients.

Sturdy: The built-in height knob eliminates unwanted movement, fostering a smooth and controlled session.

Non-Marking: Equipped on each handrail are end bumpers to safeguard equipment and the surrounding environment.

Spacious: An ample 7-foot size provides plenty of room for accommodating various exercises, offering therapists and patients more freedom of movement and versatility.

SPECIFICATIONS

Folding Parallel Bars (7 Feet)	
Material	Stainless Steel, iron, Rubber Wood, Plastic
Overall Dimensions	28" – 41" H x 83.9" W x 16" – 23.8" D
Weight Capacity	350 lb.
Height Increments	Hole distance: 0.59": The Last Hole Distance: 0.86"
Number of Height Adjustment Holes	Fourteen
Weight	123 lbs.

FEATURES



CUSTOMIZATION

ADJUSTABLE HEIGHT AND WIDTH TO SUIT USER PREFERENCES AND EXERCISE




WIDTH-ADJUSTABLE

16"-23.8"



BUILT-IN HEIGHT KNOB

28"-41"
HOLE DISTANCE: 0.6"



COMPACT

EFFORTLESSLY FOLDS TO A COMPACT SIZE FOR A SPACE-SAVING STORAGE OPTION.

