


Aluminum Wall-mounted Height Measuring Rod



ADIME907-02

Step into precision measurement with our Aluminum Wall Mounted Height Measuring Rod, engineered for accuracy and ease. Its lightweight design ensures simple installation on any wall, providing two mounting options for seamless integration into your space. The body is made from high-quality aluminum, offering steadfast durability. Whether used in a clinic, gym, or school, this height rod's adaptability ensures versatility, while its double-graduated scale in both inches and centimeters offers universal convenience.

Experience the epitome of functionality with its folding headpiece and telescopic body. Formatted for practicality, its telescopic feature allows for effortless adjustment, accommodating heights from 36.73 inches to a maximum of 64.49 inches. With a measuring range spanning 23.62 to 78.74 inches, this rod caters to diverse needs with unparalleled accuracy. Its folding headpiece can effortlessly be folded away by nurses and doctors when not in use, minimizing clutter and maintaining a streamlined environment conducive to patient care. Elevate your measurement practices today with our sleek and reliable Aluminum Wall Mounted Height Measuring Rod, setting the standard for efficiency in any professional setting.

FEATURES

Mountable: Offered in two wall mounting options for seamless integration into any space.

Durable: Made from aluminum for long-lasting performance and robustness.

Double-Graduated: Features a scale in both inches and centimeters for convenience and ease of use.


Adjustable: Effortless telescopic body, accommodating heights from 36.73 inches to a maximum of 64.49 inches, providing versatility in measurement.

Accuracy: The measuring range spans from 23.62 to 78.74 inches, catering to diverse needs with unparalleled accuracy.


Functional: Enhances practicality by allowing the headpiece to be folded away when not in use, minimizing clutter, and optimizing space.



FEATURES



Foldable Headpiece



CENTIMETERS

INCHES

Double-Graduated



MOUNTING BRACKET

With Mounting Bracket



Drilling Template

DRILLING TEMPLATE

Without Mounting Bracket

SPECIFICATIONS

SKU	ALP4250-4
Material	Aluminum
Weight	1.65 lb.
Dimensions	36.81" H x 1.33" W x 2.08" D
Units of Length	Centimeters and Inches
Height Measuring Range	23.62 to 78.74 inches

