



Rehab Mobile Treatment Mirror Series







ADIME908-01


Rehab Mobile Adult Mirror

with Cuff Weight Rack

The Rehab Mobile Adult Mirror with Cuff Weight Rack provides a professional-grade solution for real-time patient feedback and resistance training. The two-sided, safety backed glass mirror allows multiple patients to monitor their form, gait, and alignment simultaneously. By integrating a 480 lb. weight capacity into a slim 18-inch profile, this unit maximizes clinic floor space while ensuring that all your essential tools, including 48 cuff weights, 22 dumbbells, and 6 rolls of tubing, remain organized and within arm's reach.



Product Features



Dumbbell Holder

Capacity
22 Dumbbells
Compatibility
Viva Essentials Hand Weights
#ADIME920-110



Cuff Weight Hook

Capacity
48 Cuff Weights



Locking Casters

Wheel Thickness
1.2"
Wheel Dimensions
7" L x 2.5" W x 4" H




Bumper Molding

Product Features

- Two-sided mirror with safety backing allows for simultaneous use by multiple patients for gait training and posture correction.
- Designed to keep most-used tools organized and accessible, with dedicated space for 48 cuff weights, 22 dumbbells and 6 rolls of exercise tubing.
- Maximize space with heavy-duty 480 lb. capacity that allows you to organize full resistance set.
- At 18 inches wide, this unit provides dual-purpose treatment station while taking up minimal floor space.
- Equipped with industrial grade locking casters, station can be easily repositioned and locked firmly into place.
- Unit is fully compatible with Viva Essentials ADIME920-110 Hand Weights to streamline clinical workflow.
- Natural wood finish provides modern, professional look that complements any physical therapy or sports medicine facility.
- High-capacity unit includes integrated bumper molding to safeguard facility walls from scuffs.

Technical Data



Specifications

Dimensions	72" H x 18" W x 31.5" D
Weight	180 lbs.
Material	Rubber, Plywood
Weight Capacity	480 lbs.
Holding Weight Capacity	196 lbs.
Storage Capacity	48 Cuff Weights 22 Dumbbells 6 Exercise Tubing Rolls

Optional Accessories



Viva Essentials Hand Weights (10 Pack)

ADIME920-110
Set Includes
Ten Dumbbells
Weight
1, 2, 3, 4, 5, 6, 7, 8, 9, 10 lbs.

Color



Grey



ADIME908-02

Rehab Mobile Treatment Mirror

In high-traffic physical therapy and sports medicine environments, the journey to recovery is built on the precision of every movement. The Rehab Mobile Treatment Mirror was engineered to be a silent partner in that process, providing the essential visual biofeedback patients need to master their form, gait, and alignment. It effectively reclaims valuable square footage, offering a dedicated station for posture correction and therapeutic exercise without the unnecessary footprint of integrated weight racks. By focusing on a clean, single-sided safety glass surface, it allows clinicians to create a professional-grade treatment bay even in the tightest rehab environments.



Product Features



Locking Casters

Wheel Thickness
1.2"
Wheel Dimensions
7" L x 2.5" W x 4" H



Safety-Backed Glass



Bumper Molding

Product Features

- Each mirror is constructed with specialized safety backing to prevent shattering in high-traffic environments.
- At only 18 inches deep, unit provides full-scale treatment station while preserving valuable floor space.
- Built from premium-grade plywood, unit features low center of gravity and rock-solid base. Heavy-duty construction eliminates "wobble" common in lighter mobile mirrors.
- Station is mounted on high-capacity, industrial-grade four locking casters.
- Base is wrapped in protective bumper molding to safeguard walls, door frames, and equipment during transport.
- Clean, grey wood-grain finish offers modern look that integrates into any clinical environment. Materials selected for durability and ease of cleaning.
- By combining high-mobility with two-sided visibility, unit serves as versatile anchor for therapeutic activities from gait training to post-op assessments.

Technical Data



Specifications

Dimensions	72" H x 26.5" W x 25.5" D
Weight	77 lbs.
Materials	Rubber, Plywood, and High-Quality Glass Mirror
Mirror	3" D
Wheel Thickness	1.2"
Wheel Dimensions	7" L x 2.5" W x 4" H

Color



Grey