



UNIVERSAL MOBILE CART with Mounting Plates



[ADIME912-C](#) [ADIME912-C-MP-01](#) [ADIME912-C-MP-02](#)
[ADIME912-C-MP-03](#) [ADIME912-C-MP-04](#) [ADIME912-C-MP-05](#)



Improve the efficiency of your medical facility with our Universal Mobile Cart and our five different mounting plates, which are tailored to integrate with various patient vital signs monitoring systems. Our mobile cart features a telescoping pole with a height range from 31.49" to 51.18". This versatile cart offers exceptional flexibility to meet your needs and can quickly be adapted. Additionally, it includes a sturdy wire basket and cord organizer to ensure orderliness and accessibility in busy clinical environments.

Designed to withstand the rigorous demands of medical environments, our Universal Mobile Cart is built with a robust stainless-steel body. It's supported by five 2" thick wheels, two of which are equipped with 3" brakes for enhanced safety. This ensures not only smooth maneuverability but also long-term durability, making it a reliable investment for your facility.

Our Universal Mobile Cart is compatible with monitors from four leading brands and offers a quick and secure installation process. Whether you're using a Welch Allyn Connex Spot Monitor, Edan Monitor, RVS-100 Vital Sign Monitor, or Covidien PM10N Monitor Series, our cart provides a stable and reliable platform for your vital signs monitoring equipment.

FEATURES

Compatibility: Our mobile cart seamlessly integrates with various patient vital signs monitoring systems, including the Welch Allyn Connex Spot Monitor, Edan Monitor, RVS-100 Vital Sign Monitor, and Covidien PM10N Series.

Adjustable: The mobile cart's telescoping pole offers a height range of 31.49" to 51.18", accommodating different needs and users.

Outfitted: The included wire basket enhances organization by providing storage space for accessories or supplies.

Incorporated: The cart's integrated cord organizer keeps cables tidy and prevents tangling.

Construction: The body of the cart is made from durable stainless steel, ensuring longevity and resistance to wear and tear.

Stability: Our mobile cart features five 2" thick wheels with a 3" diameter each offer optimal, smooth mobility.

Secure: The two lockable wheels will keep your cart in place during use, and ensure no unwanted movement occurs.

FEATURES



Telescoping Pole

OFFERS A HEIGHT RANGE OF 31.49" TO 51.18"

HEIGHT ADJUSTMENT KNOB



Wire Basket

PROVIDES STORAGE SPACE FOR ACCESSORIES OR SUPPLIES



Cord Organizer

KEEPS CABLES TIDY AND PREVENTS TANGLING



Caster Wheels

FIVE 2" THICK WHEELS WITH A 3" DIAMETER OFFER OPTIMAL, SMOOTH MOBILITY

TWO LOCKABLE WHEELS KEEP YOUR CART IN PLACE DURING USE



SPECIFICATIONS: CART

| Universal Mobile Cart Without Mounting Plates | |
|---|--|
| SKU | ADIME912-C |
| Weight | 22.04 lb. |
| Cart Dimensions | 31.4"H x 22.5" D |
| Height Range | 31.4" to 51.1" |
| Wheels | 5, 2" thick wheels |
| Compatible Mounting Plates | Welch Allyn Connex Spot Monitor, Edan Monitor, RVS-100 Vital Sign Monitor, and Covidien PM10N Series |





Compatibility

SEAMLESSLY INTEGRATES WITH:



SPECIFICATIONS: MOUNTING PLATES



| Universal Mobile Cart With Mounting Plates | | | | | |
|--|--|--|--|---|---|
| SKU | ADIME912-C-MP-01 | ADIME912-C-MP-02 | ADIME912-C-MP-03 | ADIME912-C-MP-04 | ADIME912-C-MP-05 |
| Item Name | Universal Mobile Cart With Welch Allyn Connex Mounting Plate | Universal Mobile Cart With Welch Allyn Connex Mounting Plate With Work Surface | Universal Mobile Cart with Edan Mounting Plate | Universal Mobile Cart with RVS-100 Vital Signs Monitor Mounting Plate | Universal Mobile Cart with Covidien PM10N Series Mounting Plate |
| Weight | 27.11 lb. | 27.11 lb. | 26.01 lb. | 26.01 lb. | 26.01 lb. |
| Compatible Mounting Plate | Welch Allyn Connex Spot Monitor | Welch Allyn Connex Spot Monitor | Edan Monitor | RVS-100 Vital Signs Monitor Mounting Plate | Covidien PM10N Series Monitor |
| Mounting Plate | Welch Allyn Connex | Welch Allyn Connex with Work Surface | Edan | RVS-100 Vital Signs | Covidien PM10N Series |
| Weight | 5.07 lb. | 5.07 lb. | 3.97 lb. | 3.97 lb. | 3.97 lb. |
| Dimensions | 8.27" L x 11.22" W x 11.61" D | 8.27" L x 11.22" W x 11.61" D | 9.84" L x 4.78" W x 9.84" D | 9.84" L x 4.78" W x 9.84" D | 9.45" L x 4.17" W x 9.45" D |