



Geri Chair 3-Position Recliner



SKU: ADIME915-01-BLU



Ideal for hospitals, nursing homes, assisted living, and rehabilitation centers, the bariatric **Geri Chair 3-Position Recliner** provides exceptional comfort and therapeutic benefits for patients requiring extended seating. With three reclining positions—upright, deep recline, and elevated leg rest, our chair encourages natural relaxation, improves circulation, and reduces discomfort. The 3-inch-thick seat cushion and 2.5-inch-thick back cushion offer proper ergonomic support, while the articulating headrest adjusts for personalized comfort. For effortless maneuverability, the chair features a cushioned rear push handle and locking casters, allowing caregivers to reposition the chair securely and efficiently.

We offer a versatile optional front tray to enhance patient convenience and engagement. The front activity tray serves as a stable surface for utilizing a laptop, having meals, reading or therapy sessions. When not in use, it conveniently tucks away on the backrest frame. The included left side tray further expands functionality, locking into place when needed and folding away when not used. In order to accommodate various patient sizes, our steel-framed medical recliner supports up to 500 lbs., ensures consistent stable and reliable use. Whether assisting in patient recovery or providing a comfortable seating solution for long-term care, this chair is essential to any medical facility.

Features

Supportive: Featuring a 3-inch-thick seat cushion and a 2.5-inch-thick back cushion, designed to promote proper posture and provide long-lasting comfort.

Hygienic: The antimicrobial upholstery resists bacteria and is easy to clean.

Durable: Its powder-coated steel frame provides exceptional strength and long-lasting reliability.

High-Capacity: It can hold up to 500 lbs., accommodating various patients including bariatric patients.

Convenient: The optional front tray with a stainless steel support bar serves as a stable surface for electronics, meals, and therapy activities.

Space-Saving: The included and optional trays can be folded and stowed away when not in use.

Maneuverable: The cushioned rear push handle allows caregivers to comfortably move the chair where needed.

Secure: The medical grade locking casters ensure stability when stationary and smooth mobility.

User-Friendly: Its reclining function is controlled by a lever, allowing users to find their preferred position effortlessly.

Concealed: Its outfitted footrest fully tucks under the recliner, preventing tripping hazards and ensuring a safer, clutter-free space.

Features

500 lbs.
Weight Capacity



Color
Midnight Express Upholstery



Reclining Lever



3-Position Reclining System

- Upright
- Deep Recline
- Elevated Leg Rest

Rear Medical Grade Swivel Wheels

Brakes with foot pedal activation



Reclining Function

Controlled by a button on the right side, allowing users to find their preferred position effortlessly.



Cushioned Rear Push Handle



Trays

Large molded trays serve as a stable surface for meals, reading, and therapy activities.



Optional

Front Tray

Can be tucked away on the backrest frame

ADIME915-01-T



Included

Folded Side Tray

Located on the armrest, further expands functionality, locking into place when needed

Specifications

Geri Chair 3-Position Recliner	
Materials	Steel, Stainless Steel, Antimicrobial Fabric, Plastic
Overall Product Weight	86 lbs.
Overall Height	46.8" H x 31.1" W x 38.6" D
Overall Length	67.3 inches (Fully Reclined)
Weight Capacity	500 lbs.
Reclining Positions	Upright, Deep Recline, Elevated Leg Rest
Side Tray Location	Left Armrest, Locks into Place
Side Tray Dimensions	21.2" L x 13" W
Mobility Features	Cushioned Rear Push Handle, Locking 5-inch Casters
Reclining Mechanism	Controlled By Lever On The Right Side
Cushion Thickness	Backrest (2.5"), Seat (3"), Armrest (2")
Thickness of Wheels	1.3"
Reclining Position	24°
Backrest Adjustment	69° – 24°
Leg Rest Adjustment	90° – 24°

Front Tray (Optional)	
Front Tray Weight	4.5 lbs.
Front Tray Dimensions	25.8" W x 14.7" D

