

Gastric Drainage Constant/Intermittent Aspirators

Models 6036, 6037



Gomco 6036 and 6037 are mobile aspirators which specialize in gastric drainage for hospitals and clinics. Both models generate up to -250 mm Hg for constant suction and up to -200 mm Hg for intermittent suction, with a 20 lpm and 3-4 lpm respective flow rates at open flow. The Gomco 6036 comes with reusable collection bottle, whereas the Gomco 6037 comes with disposable collection canisters.

Key Features

- Fully independent constant and intermittent modes; each mode has it's own collection system (collection unit and suction tubing)
- Adjustable suction timing allows the clinician to tailor treatment to fit individual patients
- Easy-to-store and transport compact design with a flat stainless steel top
- Convenient storage for 6' electrical cord with hospital-grade plug
- Disposable hydrophobic bacteria filter helps prevent fluid and aerosol contamination

Specifications

Vacuum Range:

Constant: -25 to -250 mm Hg
Intermittent: -25 to -200 mm Hg
Flow Rate:

Constant: 20 lpm at open flow
Intermittent: 3-4 lpm at open flow
Cycle Times: ON: 5-90 seconds, adjustable

OFF: 12 seconds, fixed
Sound Level:
Pump Type:
Electrical Cord:
Dimensions:

OFF: 12 seconds, fixed
Less than 52 dBA
Oil-less diaphragm
6' hospital grade
23" W x 17" D x 33" H

(109 x 43 x 84 cm) Weight: 62 lbs. (28 kg)

UL/ETL Listed: Yes Warranty: 1 year Regulatory Approvals: FDA,

HC (model 6036 only)

| Model | Voltage | Glass Bottle | Poly Bottle | Disposable Canister |
|------------|------------|--------------|-------------|----------------------------|
| 01-62-6036 | 115V/60 Hz | 1 - 600 ml | 1 - 2800 ml | |
| 01-63-6036 | 220V/50 Hz | 1 - 600 ml | 1 - 2800 ml | |
| 01-62-6037 | 115V/60 Hz | | 1 - 2800 ml | 1 - 1100 ml 1 - 2100 ml |

All models come with 2 pkgs. of disposable suction tubing and 2 pkgs. of disposable hydrophobic bacteria filters.

| 01-90-2000 | Disposable Suction Tubing Kit, one 18" & one 72" per pkg. | | |
|------------|---|--|--|
| 01-90-2454 | Glass Collection Bottle, with rubber cap, short & long tube assembly, 600 ml | | |
| 01-90-3100 | Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose barb (3/pkg.) | | |
| 01-90-3694 | Polycarbonate Collection Bottle, with metal cover and cap & float assembly, 2800 ml | | |
| 01-90-3695 | Disposable Collection Canister, stem inlet and hydrophobic bacteria filter, 1100 ml (12/case) | | |
| 01-90-3712 | Disposable Collection Canister, stem inlet and hydrophobic bacteria filter, 2100 ml (10/case) | | |





Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA 314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S. Toll Free: 800/444-3940 FAX: 800/477-7701 Canada & Mexico Toll Free: 800/446-0552 FAX: 800/246-6201 International Telephone: 314/268-1683 FAX: 314/771-5183

ISO 13485: 2003

All specifications are nominal and subject to change without notice.

Form 21-00-6036-6037 July 2013