



Premium Grade Non-Woven Gauze

Made with high grade premium blended materials, our gauze provides you with superior absorbency. These superior sponges are lint free and provide more comfort with cushioning for wounds.



Non-Sterile			Packs/ case	Sterile			Pkg/ box	Box/ Case
A7632	2" x 2" 4 ply	40	A2336	8" x 4" 4 ply, 2's	25	20		
A7633	3" x 3" 4 ply	20	A2337	4" x 4" 4 ply, 2's	25	24		
A7634	4" x 4" 4 ply	10	A2338	2" x 2" 4 ply, 2's	25	60		
200 sponges per pack			A2339	3" x 3" 4 ply, 2's	25	48		

General Purpose Medical Grade Non-Woven Gauze

Non-woven gauze sponges are highly absorbent, double that of cotton gauze sponges with the equivalent plies. Leave patients with a soft and soothing feel with these lint free sponges.



Non-Sterile			Packs/ case	Sterile			Pkg/ box	Box/ Case
A2101-CH	2" x 2" 4 ply	25	A2240	2" x 2" 4 ply, 2's	50	30		
A2103-CH	3" x 3" 4 ply	20	A2240-5	2" x 2" 4 ply, 5's		100		
A2100-CH	4" x 4" 4 ply	10	A3340	3" x 3" 4 ply, 2's	50	24		
A2105-CH	3" x 4" 4 ply	20	A4440	4" x 4" 4 ply, 2's	50	12		
A2106-CH	8" x 4" 4 ply	10	A8840	8" x 4" 4 ply, 2's	25	20		
A2107	4" x 4" 6 ply	10	A9210	4" x 4" 6 ply - Drain, 2's	25	12		
A2108	2" x 2" 6 ply	25	A9210-NP	4" x 4" 6 ply, 2's	25	12		
A2220	2" x 2" 6 ply - IV	25	A9220	2" x 2" 6 ply - IV, 2's	35	20		
A4410	4" x 4" 6 ply - Drain	10	A9220-NP	2" x 2" 6 ply, 2's	35	20		
200 sponges per pack			A4440-10	4" x 4" 4 ply, 10's	10	10		
			A4440-20	4" x 4" 4 ply, 10's		50		

Cellulose Filled Non-Woven Gauze

These extra-absorbent sponges are filled with cellulose fibres for an improved performance. They are ideal for procedure that have a high volume of fluid and require a high level of absorbency.

Non-Sterile			Packs/ case	Sterile			Pkg/ box	Box/ Case
A3301	3" x 3" 8 ply	40	A0234	3" x 4" 8 ply, 2's	25	24		
A4301	4" x 3" 8 ply	20	A0244	4" x 4" 8 ply, 2's	25	24		
A4401	4" x 4" 8 ply	20						
A4801	8" x 4" 8 ply	10						
100 sponges per pack								