

Diagnostic Station

Modular System



Part Number: DS-UDM-CHFHD-DXX

Features:

- Modular design allows for easy customization
- Ergonomic handle with automatic ON/OFF function
- Low power, long life LED options
- High efficiency power supply reduces energy costs
- Smart Bulb System with timer control extends bulb life and saves energy
- Space saving with smallest equipment footprint available
- Handles are treated with MicrobeCare™ antimicrobial to assist with infection control
- Rheostat on handle allows easy control of light intensity
- Handle's spiral cord expands to over 10 ft (3 m)
- Easy to clean and disinfect
- Comes standard with anti-theft locking instrument heads
- Latex and mercury free aneroid and cuff
- Multiple mounting possibilities to suit specific needs
- Compatible with majority of manufacturers' instrument heads

Amico Diagnostic Incorporated provides quality manufactured medical equipment. The diagnostic station with Smart Bulb system saves on energy, consumables and space, while our antimicrobial treated surfaces assist with infection control.



Part Number: DS-UA1-CHFHD-DAO

Mounting Options


Wall Board



- Pre-drilled wall board allows for simple installation of diagnostic station
- Wide variety of colors and wood laminates to meet all design needs
- Color can be matched to Amico Headwalls
- Cable management behind board

Part Number: DS-ULG-CHF-DAX

Specialty Wall Boards



- Helps to provide non-clinical ambience
- Pre-drilled board with cable management behind board
- Customized shapes available


Roll Stand



- Large 5 leg base for increased stability
- Large basket allows for easy storage of accessories
- Cord wrap for proper management of power cord
- Option to accommodate thermometer and aneroid

Part Number: DS-UR1-CHF-DXX

Rail Mount



- More flexible positioning over traditional wall board
- Equipment can be locked down or moved easily
- Can be mounted on any headwall, flatwall or direct rail option

Call for part number

Module Options

Base Station



- Automatic ON/OFF independent of handle cradle
- High efficiency power supply allows simultaneous instrument operation
- Integrated cord wrap eliminates excess wiring
- Timer shut-off preserves bulb life

7YEAR
warranty

Part Number: DS-UDM-2

Specula Dispenser



- Dispenses up to 4 specula sizes
- Compartment for deep or shallow storage
- Clear housing allows for quick inventory management
- Compatible with leading brand of specula

Part Number: DS-DCD-SPDISP

Wall Aneroid



- Accurate within 3 mmHg
- 110° swivel adjustment
- Large 6" (15 cm) recessed dial face for increased visibility
- Latex and mercury free
- Lifetime calibration warranty

Part Number: AM-DB-LF2118

Electronic Thermometry



Part Number: TP-X-OP4

- Predictive or infrared options available
- Fast read response in under 4 seconds
- 60 second pulse timer option
- Hygienic single use probe covers allow touchless loading and ejecting

Part Number: TT-D-IR

Instrument Head Options



Fiber Optic Otoscope



- Fiber optics project light without obstruction
- Swivelling lens allows instrument access
- Compatible with insufflation bulb
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-FL-L: LED Part Number: DH-W35-FH-L: Halogen

Operating Otoscope



- Fiber optic illumination
- Design allows unhindered access to ear
- Swivelling lens and adjustable speculum holder

Part Number: DH-W35-OL: LED Part Number: DH-W35-OH: Halogen

Coaxial Ophthalmoscope



- 28 adjustable diopter lenses (+40 to -25)
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-CL-L: LED Part Number: DH-W35-CH-L: Halogen

Precision Coaxial Ophthalmoscope



- 68 single-diopter steps (+38 to -30) for precise focus
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-PL-L: LED Part Number: DH-W35-PH-L: Halogen

Apertures



- Small, medium, large, slit, fixation

Filters



- Cobalt blue, polarization, red-free, unfiltered

Streak Retinoscope



- Single control sleeve for vergence and rotation
- Superior streak quality enables easier refractions
- Sealed design keeps internal components dust free
- Anti-theft locking feature

Part Number: DH-W35-RH-L: Halogen

Bent Illuminator



- Fiber optic light transmission
- Provides bright consistent illumination

Part Number: DH-W35-BH: Halogen Part Number: DH-W35-BL: LED

Bi-Valve Nasal Illuminator



- Fiber optic light transmission
- Adjustable speculum
- Autoclavable

Part Number: DH-W35-VL: LED Part Number: DH-W35-VH: Halogen

Nasal Illuminator



- Fiber optic light transmission
- Head detaches to convert to throat illuminator

Part Number: DH-W35-NH: Halogen Part Number: DH-W35-NL: LED

Larynx Illuminator



- Fiber optic light transmission
- Rotating wide angle mirror with pilot aperture
- Autoclavable

Part Number: DH-W35-LL: LED Part Number: DH-W35-LH: Halogen

Throat Illuminator



- Direct light transmission
- Converts diagnostic handle to pen light


Part Number: DH-W35-TL: LED Part Number: DH-W35-TH: Halogen

Diagnostic Station Part Matrix

DS-UXX-YYZZ-QRS

Mounting Options

Wall Board



DM - Direct Wall Mount

R1 - Roll Stand

A1 - Amico Horizontal Rail Mounts

V1 - Vertical Rail Mounts

F1 - Fairfield Rail Mounts

E1 - Din 25 mm (1") Rail Mounts

G1 - Hill-Rom Vert. Rail Mounts

Head Options



Fiber Optic Oscope
FL - LED
FH - Halogen



Coaxial Ophthalmoscope
CL - LED
CH - Halogen



Precision Coaxial Ophthalmoscope
PL - LED
PH - Halogen



Bent Illuminator
BL - LED
BH - Halogen



Streak Retinoscope
RH - Halogen



Operating Oscope
OL - LED
OH - Halogen



Throat Illuminator
TL - LED
TH - Halogen



Bi-Valve Nasal Illuminator
VL - LED
VH - Halogen



Nasal Illuminator
NL - LED
NH - Halogen



Larynx Illuminator
LL - LED
LH - Halogen

Thermometry Options

- O - Oral/Axillary Predictive
- I - Tympanic Infrared
- X - No Thermometer



Aneroid and Basket Options

- A - Aneroid
- L - L Bracket Aneroid
- X - No Aneroid



- D - Specula Dispenser
- X - No Dispenser



Diagnostic Accessories

Disposable Ear Specula



- Single use ear specula
- Available in a wide variety of sizes
- Operating otoscope specula also available

DC-SSD-1K-25M: 2.5 mm (0.11"), pack of 1,000
DC-SSD-1K-40M: 4.0 mm (0.16"), pack of 1,000

Reusable Ear Specula



- Autoclavable up to 134°C (273°F)
- Available in a wide variety of sizes

DC-SSR-01-25M: 2.5 mm (0.11")
DC-SSR-01-30M: 3.0 mm (0.14")
DC-SSR-01-40M: 4.0 mm (0.16")
DC-SSR-01-50M: 5.0 mm (0.20")

Insufflation Bulb



- Oscope insufflation bulb and tube
- Compatible with fiber optic otoscope

DA-05-IB

Predictive Therm. Probe Covers



- Disposable single-use
- Touchless loading and ejecting of covers

500500: pack of 500
502000: pack of 2000

Infrared Therm. Probe Covers



- Disposable single-use
- Touchless loading and ejecting of covers

303030: Pack of 2,112

Replacement Lamps



- Oscope & Illuminator
DC-3100-L01: LED pack of 1
DC-3100-H10: Halogen pack of 10
- Ophthalmoscope
DC-4900-L01: LED pack of 1
DC-4900-H10: Halogen pack of 10
- Retinoscope
DC-L23-H10: Halogen pack of 10

Two Piece Cuff




- Latex-free bladder
- Variety of sizes available
- Connector options to fit any NIBP device

Stethoscopes



- Stainless steel double-head
- Non-chill rim
- Adult and pediatric available
- Latex-free

Wire Baskets



- 19" x 5" x 3"
(48.26 cm x 12.7 cm x 7.62 cm)
- Various sizes available

Call for Part Number

Distributed By:

