

| | |
|----------------------|--|
| Name: | |
| School/Organization: | |
| Dealer/Rep: | |

| | | | |
|------------|--|--------|--|
| Phone: | | | |
| E-mail: | | | |
| Sales Rep: | | Quote: | |

AM SERIES

27" HI-LO TREATMENT TABLES, NO CASTERS

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3[®] bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base



AM-100

TABLE CONFIGURATIONS:

| | | |
|---|--|---|
| <input type="radio"/> AM-100 - 1 Section | <input type="radio"/> AM-1227 - 2 Section | <input type="radio"/> AM-353 - 3 Section |
| QTY | QTY | QTY |



TABLE FEATURES:

| | | |
|--|--|--|
|  <p>Flat Cushion</p> |  <p>Lower Extremity Angle Ranges from 0 to 75°</p> |  <p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p> |
| AM-100 | AM-1227 | AM-353 |

OPTIONS:

| |
|--|
| <input type="radio"/> AM-804 - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls) |
| <input type="radio"/> AM-91 - Add 1" Extra Soft Foam |
| <input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer) |
| <input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color) |

STANDARD UPHOLSTERY CHOICES:

| | |
|-------------------------------------|---------------------------------|
| <input type="radio"/> Imperial Blue | <input type="radio"/> Taupe |
| <input type="radio"/> Blue Ridge | <input type="radio"/> Burgundy |
| <input type="radio"/> Forest Green | <input type="radio"/> Dove Gray |
| <input type="radio"/> Black | <input type="radio"/> Greystone |
| <input type="radio"/> Cappuccino | <input type="radio"/> Merlot |
| <input type="radio"/> Patina | |

STANDARD BASE COLOR:

| |
|------------------------------------|
| <input type="radio"/> Office White |
|------------------------------------|

CUSTOM BASE COLORS:

| |
|---------------------------------------|
| <input type="radio"/> Black |
| <input type="radio"/> Steel Gray |
| <input type="radio"/> Metallic Silver |

| | |
|----------------------|--|
| Name: | |
| School/Organization: | |
| Dealer/Rep: | |

| | |
|------------|--------|
| Phone: | |
| E-mail: | |
| Sales Rep: | Quote: |

AM SERIES

27" HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base
- Dual front in-line casters with 3rd wheel drop down for easy maneuverability



AM-350

TABLE CONFIGURATIONS:

| | | | | |
|---|---|---|---|---|
| <input type="radio"/> AM-150 - 1 Section | <input type="radio"/> AM-200 - 2 Section | <input type="radio"/> AM-227 - 2 Section | <input type="radio"/> AM-250 - 2 Section with 3pc Head Section | <input type="radio"/> AM-350 - 3 Section |
| QTY | QTY | QTY | QTY | QTY |



| | | | | |
|---|---|--|--|---|
| <input type="radio"/> AM-300 - 3 Section with Elevating Midsection | <input type="radio"/> AM-550 - 5 Section with Tilt Down Armrests | <input type="radio"/> AM-555 - 5 Section with Adjustable Armrests | <input type="radio"/> AM-500 - 5 Section with Elevating Midsection & Tilt Down Armrests | <input type="radio"/> AM-505 - 5 Section with Elevating Midsection & Adjustable Armrests |
| QTY | QTY | QTY | QTY | QTY |



TABLE FEATURES:

| | | | | |
|--|--|--|--|---|
|  <p>Flat Cushion</p> |  <p>Head Section Angle Ranges from +45° to -90°</p> |  <p>Lower Extremity Angle Ranges from 0 to 78°</p> |  <p>Head Section Angle Ranges from +45° to -90° Armrests Angle Ranges 0° to -45°</p> |  <p>Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°</p> |
| AM-150 | AM-200 | AM-227 | AM-250 | AM-350 |
|  <p>Lower Extremity Angle Ranges from 0 to 78° Manual Elevating Midsection</p> |  <p>Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°</p> |  <p>Lower Extremity Angle Ranges from 0 to 78° Armrests Vertical Adjustment 8" down</p> |  <p>Armrests Angle Ranges 0° to -45° Manual Elevating Midsection</p> |  <p>Armrests Vertical Adjustment 8" down Manual Elevating Midsection</p> |
| AM-300 | AM-550 | AM-555 | AM-500 | AM-505 |

OPTIONS:

| |
|---|
| <input type="radio"/> AM-1024 - Add Base with 4 Total Lock Swivel Casters |
| <input type="radio"/> AM-1060 - Add Base with 4 Casters and 6" Clearance (Adds 4") |
| <input type="radio"/> AM-804 - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls) |
| <input type="radio"/> AM-18 - Face Cutout |
| <input type="radio"/> AM-17 - Face Plug (Must order AM-18 also if ordering AM-150) |
| <input type="radio"/> AM-91 - Add 1" Extra Soft Foam |
| <input type="radio"/> AM-807 - Rail Side Fold Down (Pair for 27"W Tables) |
| <input type="radio"/> AM-808 - Rail Side Fold Down (Pair for 30"W Tables) |
| <input type="radio"/> AM-875 - Shoulder Mobilization Accessory (27"W, 3 & 5 Section Models Only) |
| <input type="radio"/> AM-1025 - Adjustable Armrest Package (2 Section Only) |
| <input type="radio"/> SP30 - 30"W |
| <input type="radio"/> 19855 - Add-on Belt & Brackets (27" W) |
| <input type="radio"/> 19859 - Add-on Belt & Brackets (30" W) |
| <input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer) |
| <input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color) |
| <input type="radio"/> 12088 - Side Panel 9 x 33.25 Perforated |
| <input type="radio"/> 12075 - Side Panel 9 x 33.25 Flat |

STANDARD UPHOLSTERY CHOICES:

| | |
|-------------------------------------|---------------------------------|
| <input type="radio"/> Imperial Blue | <input type="radio"/> Taupe |
| <input type="radio"/> Blue Ridge | <input type="radio"/> Burgundy |
| <input type="radio"/> Forest Green | <input type="radio"/> Dove Gray |
| <input type="radio"/> Black | <input type="radio"/> Greystone |
| <input type="radio"/> Cappuccino | <input type="radio"/> Merlot |
| <input type="radio"/> Patina | |

STANDARD BASE COLOR:

| |
|------------------------------------|
| <input type="radio"/> Office White |
|------------------------------------|

CUSTOM BASE COLORS:

| |
|---------------------------------------|
| <input type="radio"/> Black |
| <input type="radio"/> Steel Gray |
| <input type="radio"/> Metallic Silver |