

AM SERIES

Hi-Lo Treatment Tables



AM-300 Treatment Table

- | | | | |
|---|---|---|---|
| <ul style="list-style-type: none">• 3-Section Treatment Table• 76"L x 27"W | <ul style="list-style-type: none">• Min. Height: 18"• Max Height: 37"• 1 1/2" Firm Density Foam | <ul style="list-style-type: none">• Head: 15" x 27" with contoured face opening• Center: 22" x 27"• Foot: 39" x 27" | <ul style="list-style-type: none">• 400lbs. Capacity• Power-assisted Caster System |
|---|---|---|---|

Clinicians appreciate the full height adjustment and section mobility options so the patient can be placed in the proper position for maximum results with minimum effort.

All tables have these features:

- Comfort and easy access getting on and off of the table
- 34oz. antibacterial vinyl
- 18 month warranty
- Heavy duty welded steel construction

AM-BA SERIES

BAR ACTIVATED HI-LO CONTROL



AM-BA300 Treatment Table

- | | | | |
|---|---|---|--|
| <ul style="list-style-type: none">• 3-Section Treatment Table• 27"W x 76"L | <ul style="list-style-type: none">• Min. Height: 18"• Max Height: 37"• 1 1/2" Firm Density Foam | <ul style="list-style-type: none">• Head: 15" x 27" with Contoured Face Opening• Center: 22" x 27"• Foot: 39" x 27" | <ul style="list-style-type: none">• 400 lb. Capacity• 4 Total Lock Swivel Casters |
|---|---|---|--|

No more searching for the control switch! This unique design makes the control easy for the clinician to use from anywhere around the table, yet out of the patient's way as they get on and off the table.