

Name:	
School/Organization:	
Dealer/Rep:	

Phone:	
E-mail:	
Sales Rep:	Quote:

MAT PLATFORM TABLES

WITH ELECTRIC HEIGHT ADJUSTMENT

Conforms to ANSI/UL Std. IEC 60601, Certified to CSA Std. C22.1, Std. 60601.1

ETL CLASSIFIED



Intertek
4009555

STANDARD FEATURES:

- Load capacity: 1000 lbs
- Electric elevation: 20-30" H
- ETL Certified
- Indestructible tubing to prevent kinks
- Light up display to indicate standby mode/power applied
- Audible warning with delay safety features when device lowers
- 34 oz. Heavy-duty vinyl with Permablok 3[®] bacterial protection
- 2" High density foam
- Single mounted hand control
- Rounded corners




AMTM-640

STANDARD CONFIGURATIONS:

<input type="radio"/> AMTM-640 - 48"W x 84"L x 20-30"H	<input type="radio"/> AMTM-641 - Adjustable Backrest (Raises to 80") 48"W x 84"L x 20-30"H	<input type="radio"/> AMTM-650 - 60"W x 84"L x 20-30"H	<input type="radio"/> AMTM-660 - 72"W x 96"L x 20-30"H
QTY:	QTY:	QTY:	QTY:

OPTIONS:

<input type="radio"/> AMTM-03611-FS - Substitute a Foot Switch for Single Mounted Hand Control
QTY:
<input type="radio"/> AMTM-03611-FS-2 - Add a Second Foot Switch (Not available with Hand Control)
QTY:
<input type="radio"/> AMTM-03619-BMS-2 - Add a Second Hand Control (Not available with Foot Switch)
QTY:
<input type="radio"/> AMSBO-047 - Side Bumpers Mat Table, 4' x 7' (Not available with AMTM-641)
QTY:
<input type="radio"/> AMSBO-057 - Side Bumpers Mat Table, 5' x 7'
QTY:
<input type="radio"/> AMSBO-068 - Side Bumpers Mat Table, 6' x 8'
QTY:



AMSBO-047

STANDARD BASE COLOR:

<input type="radio"/>	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

<input type="radio"/>	Black
<input type="radio"/>	Steel Gray
<input type="radio"/>	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

<input type="radio"/>	Imperial Blue	<input type="radio"/>	Taupe
<input type="radio"/>	Blue Ridge	<input type="radio"/>	Burgundy
<input type="radio"/>	Forest Green	<input type="radio"/>	Dove Gray
<input type="radio"/>	Black	<input type="radio"/>	Greystone
<input type="radio"/>	Cappuccino	<input type="radio"/>	Merlot
<input type="radio"/>	Patina		



Armedica[®]
BY PIVOTAL HEALTH SOLUTIONS

800.701.5122
armedicamfg.com