

AM-BA SERIES

BAR ACTIVATED HI-LO CONTROL



AM-BA300 Treatment Table

- 3-Section Treatment Table
- 27"W x 76"L
- Min. Height: 18"
- Max Height: 37"
- 1 1/2" Firm Density Foam
- Head: 15" x 27" with Contoured Face Opening
- Center: 22" x 27"
- Foot: 39" x 27"
- 400 lb. Capacity
- 4 Total Lock Swivel Casters

No more searching for the control switch! This unique design makes the control easy for the clinician to use from anywhere around the table, yet out of the patient's way as they get on and off the table.