

Disposable Light Handle and Camera Covers

Provides a sterile barrier between non-sterile surgical lights and cameras and the sterile surgical field.

 SURGICAL AND MEDICAL ESSENTIALS

 aspen surgical®

Specialty Light Handle and Camera Covers

SURGICAL AND MEDICAL ESSENTIALS

Family of Specialty handles and covers designed to directly interface with multiple vendors lights and cameras without the need of an adapter. These solutions provide the sterile interface between the non-sterile operating room camera and surgical light from the sterile field of the nurse and / or surgeon.

Features for all Specialty Solutions include:

- Helps to save time and money associated with reprocessing reusable handles & covers
- Helps ensure quicker Operating Room turnaround times
- Sterile, intended for single use, and not made with natural rubber latex

Key Features that are unique to STERIS® include:

- Our covers are a direct replacements for covers provided by STERIS
- The Aspen Surgical solution is more rigid enabling easier installation and provides more coverage of the light and camera buttons
- The Camera Cover solution has an ultra-durable, clear lens enabling distortion free, complete visualization of the surgical field

Key Features that are unique to Trumpf® include:

- Our handle is designed to have a similar look and feel as the Trumpf Medical reusable handle
- The Aspen Surgical solution offers a snap-on feature enabling quick and easy installation and removal
- Ultra-durable lens offers clarity equivalent to the reusable handle, providing rapid response of the laser and sensor
- Light intensity can be adjusted with the touch sensors while the disposable light handle is being used

Specialty Light Handle and Camera Covers

 SURGICAL AND MEDICAL ESSENTIALS



DISPOSABLE O.R. LIGHT HANDLE COVER

Disposable O.R. Light Handle Cover, comparable to STERIS® LB53 Blue Sterile 120/Box (Compatible with HarmonyAIR® E-Series, HarmonyAIR® M-Series, HarmonyAIR® G-Series, Harmony LED, Harmony vLED, Harmony LA)

Part no.	Color	Sterility	Qty
LT-53	Blue	Sterile	120/Box



DISPOSABLE O.R. CAMERA COVER

Disposable O.R. Camera Cover, comparable to STERIS® LB54 Blue Sterile 30/Box (Compatible with HarmonyAIR® E-Series, HarmonyAIR® M-Series, HarmonyAIR® G-Series, FREE5® Camera System, Other legacy STERIS® Harmony Surgical Lighting Systems)

Part no.	Color	Sterility	Qty
LT-54	Blue	Sterile	30/Box



DISPOSABLE LIGHT HANDLE COVER

For use with Trumpf Medical™ ALC Plus Surgical Light Models: TruLight® 500 and iLED® 3/5.

Part no.	Color	Sterility	Qty
LT-ALC01	Green	Sterile	18/Box

Standard Disposable Light Handle and Camera Covers

SURGICAL AND MEDICAL ESSENTIALS

Family of Standard handles and covers that provide the sterile interface between the non-sterile operating room surgical light and camera from the sterile field of the nurse and / or surgeon. Our Standard handles and covers require an adapter sold separately.

Features for all standard handles and covers include:

- Standard handles and covers are offered in both rigid and flexible configurations
- Creates uniformity throughout hospital by utilizing one light handle product for multiple surgical lights
- Helps to save time and money associated with reprocessing reusable handles
- Convenient Kit item, helps ensure quicker operating room turnaround times
- Sterile, intended for single use, and not made with natural rubber latex

Key Features that are unique to our patented Camera Covers include:

- The Camera Cover solution has an ultra-durable, clear lens enabling distortion free, complete visualization of the surgical field
- Patented designs for both rigid and camera covers deliver high levels of efficiency and safety

Standard Disposable Light Handle and Camera Covers

 SURGICAL AND MEDICAL ESSENTIALS



LT-C02 | O.R. CAMERA COVER

LT-CF02 | FLEXIBLE O.R. CAMERA COVER

Part no.	Description	Color	Sterility	Qty
LT-C02	O.R. Camera Cover	Green	Sterile	20/Box
LT-CF02	Flexible O.R. Camera Cover	Green	Sterile	40/Box



LT-R01 | RIGID HANDLE

LT-F01B | FLEXIBLE COVER

Part no.	Description	Color	Sterility	Qty
LT-R01	Rigid Light Handle	Green	Sterile	24/Box
LT-F01B	Flexible Light Handle Cover (1 each)	Green	Sterile	100/Box
LT-F02	Flexible Light Handle Cover (2 each)	Green	Sterile	80/Bow
LT-F03	Flexible Light Handle Cover (3 each)	Green	Sterile	80/Box



Please work with your local Aspen Surgical Representative for information regarding adapters that are required for the Aspen Surgical Standard Light and Camera handles and covers.



For more information, contact your Sales Representative or visit www.aspensurgical.com.