

Alliance Aluminum

The Alliance Aluminum is a lightweight table that offers exceptional strength, unbeatable client access, and positioning accuracy. This table is perfect for any of today's most popular treatments in the spa, medi-spa, and massage fields.




PATENTED INTEGRACHINGE™



UNILOCK™



ULTIMATE ACCESS END PANEL
& SHIATSU CABLES



EXCEPTIONAL
FRAME STRENGTH



Shown with Heron top

TABLE SPECIFICATIONS

Max. Operating Weight	550 lbs.
Weight	30-36 lbs.
Height Ranges	24"-34"
Widths	30"
Length	73"
Top	Backrest, 0-60°, 6 Positions
Padding	Semi-Firm
Standard	Access Panel
Accessories	QuickLock Face Rest Platform, Arm Hammock, and our full line of accessories
Fabric	TerraTouch

TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



Contact Oakworks
717.235.6807
www.massagetables.com
customerservice@oakworks.com