



Aristo H-Brace Line

ADI996-03-27-BLU ADI996-03S-30-BLU ADI996-04-27-BLU ADI996-04S-30-BLU
ADI996-03-27-BLK ADI996-03S-30-BLK ADI996-04-27-BLK ADI996-04S-30-BLK



ADI996-03-27-BLU
Aristo H-Brace Flat Treatment Table
with Adjustable Back



ADI996-03S-30-BLK
Aristo H-Brace Table 30\"/>



ADI996-04-27-BLU
27\"/>



ADI996-04S-27-BLU
27\"/>

Aristo H-Brace Flat Treatment Table

With Adjustable Back

- Adjustable backrest reclines from 0° to 67°.
- Fortified with a support brace for added stability.
- Weight capacity up to 500 lbs.
- Sturdy wooden base.
- 2" high-density foam and antimicrobial PVC vinyl upholstery.
- **Optional:** Double drawers, shelf cabinets, paper cutter and dispenser.

Upholstery Colors

Midnight Express
#ADI996-03-27-BLU

Black
#ADI996-03-27-BLK



Features



Solid Wooden Base



Adjustable Backrest



2-inch Thick Mattress



Antimicrobial Upholstery

Optional Accessories



Double Drawer
#ADI996-DD



Paper Holder
#ADI996-PA



Shelf Cabinet
#ADI996-SC



Paper Cutter
#ADI903-20PC

Dimensions



Specifications

Materials	PVC Vinyl, Foam, Plywood for Mattress, and Beech Wood
Overall Dimensions	31" H x 27" W x 72" L
Adjustable Backrest Dimensions	30" L x 27" W x 2.55 H
Other Half of Mattress Dimensions	42.5" L x 27" W x 2.55 H
Weight	61.51 lb.
Weight Capacity	500 lb.
Padding	2" thick high-density foam
Upholstery	Antimicrobial

Aristo H-Brace Table 30" Tilted Mattress

With Shelf

- Adjustable backrest reclines from 0° to 67°.
- Fortified with a support brace for added stability.
- Weight capacity up to 500 lbs.
- Sturdy wooden base.
- 2" high-density foam and antimicrobial PVC vinyl upholstery.
- **Included:** The open shelf provides convenient storage for essential medical supplies and equipment.
- **Optional:** Double drawers, shelf cabinets, paper cutter and dispenser.

Upholstery Colors

	Midnight Express #ADI996-03S-30-BLU		Black #ADI996-03S-30-BLK
--	---	---	------------------------------------



Features



Open Shelf



Adjustable Backrest



2-inch Thick Mattress



Antimicrobial Upholstery

Optional Accessories



Double Drawer
#ADI996-DD



Paper Holder
#ADI996-PA



Shelf Cabinet
#ADI996-SC



Paper Cutter
#ADI903-20PC

Dimensions



Specifications

Materials	PVC Vinyl, Foam, Plywood for Mattress, and Beech Wood
Overall Dimensions	31" H x 30" W x 72" L
Adjustable Backrest Dimensions	30" L x 30" W x 2.55 H
Other Half of Mattress Dimensions	42.5" L x 30" W x 2.55 H
Weight	101.31 lb.
Weight Capacity	500 lb.
Padding	2" thick high-density foam
Upholstery	Antimicrobial

27" Aristo H-Brace Flat Treatment Table

- Fortified with a support brace for added stability.
- Weight capacity up to 500 lbs.
- Sturdy wooden base.
- 2" high-density foam and antimicrobial PVC vinyl upholstery.
- **Optional:** Double drawers, shelf cabinets, paper cutter and dispenser.

Upholstery Colors

Midnight Express
#ADI996-04-27-BLU

Black
#ADI996-04-27-BLK



Features



Solid Wooden Base



Fortified with Support Brace



2-inch Thick Mattress



Antimicrobial Upholstery

Optional Accessories



Double Drawer
#ADI996-DD



Paper Holder
#ADI996-PA




Shelf Cabinet
#ADI996-SC



Paper Cutter
#ADI903-20PC

Dimensions



Specifications

Materials	PVC Vinyl, Foam, Plywood for Mattress, and Beech Wood
Overall Dimensions	31" H x 27" W x 72" L
Weight	61.51 lb.
Weight Capacity	500 lb.
Padding	2" thick high-density foam
Upholstery	Antimicrobial

27" Aristo H-Brace Flat Treatment Table

With Shelf


- Fortified with a support brace for added stability.
- Weight capacity up to 500 lbs.
- Sturdy wooden base.
- 2" high-density foam and antimicrobial PVC vinyl upholstery.
- **Included:** The open shelf provides convenient storage for essential medical supplies and equipment.
- **Optional:** Double drawers, shelf cabinets, paper cutter and dispenser.

Upholstery Colors

	Midnight Express #ADI996-04S-27-BLU		Black #ADI996-04S-27-BLK
--	---	---	------------------------------------



Features



Solid Wooden Base



Open Shelf



2-inch Thick Mattress



Antimicrobial Upholstery

Optional Accessories



Double Drawer
#ADI996-DD



Paper Holder
#ADI996-PA




Shelf Cabinet
#ADI996-SC



Paper Cutter
#ADI903-20PC

Dimensions



Specifications

Materials	PVC Vinyl, Foam, Plywood for Mattress, and Beech Wood
Overall Dimensions	31" H x 27" W x 72" L
Weight	101.31 lb.
Weight Capacity	500 lb.
Padding	2" thick high-density foam
Upholstery	Antimicrobial