

BD™ Universal Viral Transport System

Adding More Flexibility for Viral Testing

BD™ Universal Viral
Transport for Viral,
Chlamydial, Mycoplasmal and
Ureaplasma Specimens
3 mL Vial

Becton, Dickinson and Company
7 Loveland Circle
Sparks, Maryland 21152 USA
Made in Italy by Copan

05A146



REF 220220

50

IVD

CE

LOT

BD Universal Viral Transport
Medium for Viruses,
Chlamydiae, Mycoplasmas
and Ureaplasmas (3.0 mL)

LOT

BD Universal Viral Transport
Medium for Viruses,
Chlamydiae, Mycoplasmas
and Ureaplasmas (1.0 mL)

Flexible, Compatible and Uniquely Designed!

- Room temperature and refrigerated viability in ONE formulation
- Compatible transport medium for the BD Directigen™ EZ viral tests
- Designed with safety and convenience to meet your viral testing needs



BD

Helping all people
live healthy lives

BD™ Universal Viral Transport System



Designed to prevent spills.



Unique Formulation

- Sustains viability for viruses, chlamydiae, mycoplasmas and ureaplasmas in a single formulation.
- Offers flexibility in transport temperature ranges—room temperature or refrigerated.
- Suppresses bacterial and fungal contamination by incorporating antibiotics in the medium.
- Preserves viruses and chlamydiae for long-term frozen storage.
- Utilizes a phenol red pH indicator to ensure medium integrity at the time of specimen collection.

Efficient Design

- Available in 1 mL or 3 mL vials to allow for testing flexibility.
- Flat-bottomed conical tube stands on the benchtop which reduces the risk of biohazardous spills.
- Capture-cap on 3 mL vial eliminates the need for forceps to remove the swab and reduces the risk of contamination.
- Scored swabs eliminate the need for scissors.

Compatibility

- Acceptable transport medium for the BD Directigen EZ Flu A+B and BD Directigen EZ RSV test kits.
- Able to be used with flocked swabs which feature nylon strands that use capillary action to absorb and elute the entire sample into liquid media.

Cat. No.	Description	Qty/Pkg
220220	BD™ Universal Viral Transport, 3 mL vial	Pkg. of 50
220221	BD™ Universal Viral Transport Standard Kit: 3 mL vial with 2 regular polyester-tipped swabs	Pkg. of 50
220222	BD™ Universal Viral Transport Combo Kit: 3 mL vial with 1 regular and 1 minitip wire/plastic swab	Pkg. of 50
220239	BD™ Universal Viral Transport 2 regular polyester-tipped swabs	Pkg. of 100
220240	BD™ Universal Viral Transport Combo Polyester Swabs: 1 regular and 1 minitip wire/plastic swab	Pkg. of 100
220244	BD™ Universal Viral Transport, 1 mL vial	Pkg. of 50
220526	BD™ Universal Viral Transport Kit: 1 mL vial with a flexible minitip flocked swab	Pkg. of 50
220527	BD™ Universal Viral Transport Kit: 3 mL vial with a regular and flexible minitip flocked swab	Pkg. of 50
220528	BD™ Universal Viral Transport Kit: 3 mL vial with a regular flocked swab	Pkg. of 50
220529	BD™ Universal Viral Transport Kit: 3 mL vial with a minitip flocked swab	Pkg. of 50
220531	BD™ Universal Viral Transport Kit: 3 mL vial with a flexible minitip flocked swab	Pkg. of 50
256050	BD Directigen™ EZ Flu A+B Test Kit	30 Tests
256030	BD Directigen™ EZ RSV Test Kit	30 Tests



BD Diagnostics
7 Loveton Circle
Sparks, MD 21152-0999
800.638.8663
www.bd.com/ds