

Single-use surgical products

Unstrung tonsil sponges

Unstrung tonsil sponges provide excellent absorbency, and they can be used as prep sponges for surgery.



BD cat no.	Description	Packaging
23275-650	¾ in. (1.9 cm), unstrung, cotton	5/pack 50 packs/case
23275-660	1 ¼ in. (3.2 cm), unstrung, cotton	5/pack 50 packs/case

Double strung tonsil sponges

Double strung tonsil sponges provide excellent absorbency, and they are constructed with thread for ease of retrieval.




BD cat no.	Description	Packaging
23275-670	¾ in. (1.9 cm), double strung, cotton	5/pack 50 packs/case
23275-680	¾ in. (2.2 cm), double strung, cotton	5/pack 50 packs/case
23275-690	1 ¼ in. (3.2 cm), double strung, cotton	5/pack 50 packs/case

Patties

Rayon neuro patties


Rayon neuro patties are very pliable when wet, allowing them to conform to anatomy. They are constructed with thread for visibility. They are also x-ray detectable.



BD cat no.	Description	Packaging
23275-310	¼ in. x ¼ in. (.6 cm x .6 cm), strung, rayon	10/pack 20 packs/case
23275-330	¼ in. x 3 in. (.6 cm x 7.6 cm), strung, rayon	10/pack 20 packs/case
23275-340	¼ in. x 1 ½ in. (.6 cm x 3.8 cm), strung, rayon	10/pack 20 packs/case
23275-350	¾ in. x ¾ in. (1.9 cm x 1.9 cm), strung, rayon	10/pack 20 packs/case
23275-510	½ in. x ½ in. (1.3 cm x 1.3 cm), strung, rayon	10/pack 20 packs/case
23275-520	½ in. x 1 in. (1.3 cm x 2.5 cm), strung, rayon	10/pack 20 packs/case
23275-530	½ in. x 3 in. (1.3 cm x 7.6 cm), strung, rayon	10/pack 20 packs/case
23275-540	½ in. x 1 ½ in. (1.3 cm x 3.8 cm), strung, rayon	10/pack 20 packs/case
23275-550	1 in. x 1 in. (2.5 cm x 2.5 cm), strung, rayon	10/pack 20 packs/case
23275-560	1 in. x 3 in. (2.5 cm x 7.6 cm), strung, rayon	10/pack 20 packs/case
23275-570	½ in. x 2 in. (1.3 cm x 5 cm), strung, rayon	10/pack 20 packs/case
23275-580	3 in. x 3 in. (7.6 cm x 7.6 cm), strung, rayon	10/pack 20 packs/case

Cotton patties

Cotton patties are highly absorbent, and they conform to anatomy and maintain integrity when wet. They can be used to moisten or absorb fluid during neurosurgery, and may also be used to protect the brain during incision. They are also x-ray detectable.



BD cat no.	Description	Packaging
NL9051	½ in. x ½ in. (1.3 cm x 1.3 cm), double strung, cotton	10/pack 20 packs/case
NL9053	½ in. x 1 in. (1.3 cm x 2.5 cm), double strung, cotton	10/pack 20 packs/case
NL9055	½ in. x 1 ½ in. (1.3 cm x 3.8 cm), double strung, cotton	10/pack 20 packs/case
NL9058	½ in. x 3 in. (1.3 cm x 7.6 cm), double strung, cotton	10/pack 20 packs/case
NL9052	¾ in. x ¾ in. (1.9 cm x 1.9 cm), double strung, cotton	10/pack 20 packs/case
NL9054	1 in. x 1 in. (2.5 cm x 2.5 cm), double strung, cotton	10/pack 20 packs/case
NL9059	1 in. x 3 in. (2.5 cm x 7.6 cm), double strung, cotton	10/pack 20 packs/case
NL9060	3 in. x 3 in. (7.6 cm x 7.6 cm), double strung, cotton	10/pack 10 packs/case

Spears

Polyvinyl acetal (PVA) surgical spears



Polyvinyl acetal (PVA) surgical spears feature a high tensile strength and allow gentle manipulation of delicate tissues. They are designed for use in LASIK eye and cataract surgeries. They absorb up to five times more fluid than standard cellulose spears and expand instantly upon contact with fluids.

BD cat no.	Description	Packaging
23275-500	7 cm length	10/pack 20 packs/case

Trach/drain sponges

Trach/drain sponges are highly absorbent, and they can be used around any drainage area to protect it from chafing.



BD cat no.	Description	Packaging
23276-530	4 in. x 4 in. (10.2 cm x 10.2 cm), 12-ply gauze, cotton	2/pack 100 packs/case

Vascular loops

V. Mueller single-use vascular loops are designed for occluding, retracting and identifying arteries, veins, tendons and nerves in surgical procedures. Vascular loops are made of medical-grade, radiopaque silicone and are available in four colors for easy identification. They are supplied sterile and not made with natural rubber latex.



BD cat no.	Description	Packaging
Maxi thickness (0.10" W x 14.75" L x 0.04" D [2.5 mm x 375 mm x 1 mm])		
CH111	Red	2/pack 10 packs/box
CH112	Yellow	2/pack 10 packs/box
CH113	White	2/pack 10 packs/box
CH114	Blue	2/pack 10 packs/box
Mini thickness (0.06" W x 14.75" L x 0.035" D [1.5 mm x 375 mm x 0.9 mm])		
CH115	Red	2/pack 10 packs/box
CH116	Yellow	2/pack 10 packs/box
CH117	White	2/pack 10 packs/box
CH118	Blue	2/pack 10 packs/box

**All dissectors and sponges are x-ray detectable. They are provided sterile and not made with natural rubber latex.