

Single-use surgical products

Insufflation products and accessories

V. Mueller™ insufflation products include adapters, tubing and Veress needles, adaptable to all insufflator brands for tubing standardization. The rotating luer lock connector on the durable, non-collapsing tubing helps prevent twisting or disconnection from the insufflator. The tubing is available with bidirectional 0.1 or 0.3 micron hydrophobic filters to reduce the chance of contamination. It is sterile and not made with natural rubber latex.



| BD cat no. | Description | Packaging |
|--|---|-----------|
| Single-use insufflation tubing | | |
| LA101-RLL | 0.1 micron filter, 10' length | 20/box |
| LA103-RLL | 0.3 micron filter, 10' length | 20/box |
| Reusable insufflation tubing adapters | | |
| LAUNV-LL | Snowden-Pencer™, Circon, Smith & Nephew and Olympus® insufflators | 1/box |
| LASTK-LL | Stryker® insufflators | 1/box |
| LALNV-LL | Linvatec® insufflators | 1/box |
| LAWLF-LL | Richard Wolf™ insufflators | 1/box |
| LASTZ-LL | KARL STORZ™ insufflators | 1/box |
| Single-use Veress needle | | |
| 007.00 | 14 G, 120 mm | 12/box |
| 007.01 | 14 G, 150 mm | 12/box |

DeFogger™ anti-fog solution

Each bottle of DeFogger™ anti-fog solution comes with a nonabrasive, x-ray detectable foam pad that allows universal placement and will not damage the laparoscope. It is sterile and not made with natural rubber latex.



| BD cat no. | Description | Packaging |
|------------|----------------------|-----------|
| LAF0G-FT | 6 gm flip-top bottle | 20/box |

Dissectors and sponges

Round stick sponges

Round stick sponges provide excellent absorbency, and they can be used as prep sponges.



| BD cat no. | Description | Packaging |
|------------|--------------------------|-----------------------|
| 23275-740 | 1 in. (2.5 cm), cotton | 10/pack 25 packs/case |
| 23275-745 | 1 ½ in. (3.8 cm), cotton | 10/pack 25 packs/case |

Flat stick sponges

Flat stick sponges provide excellent absorbency, and they can be used as prep sponges for surgery or with a FOERSTER clamp for deep sponge sticking.



| BD cat no. | Description | Packaging |
|------------|--------------------------|-----------------------|
| 23275-750 | 1 ½ in. (3.8 cm), cotton | 10/pack 25 packs/case |

K-dissector sponges

K-dissector sponges can be used as longitudinal or folded sponges for dissection (they include a count holder for easy sponge counts).



| BD cat no. | Description | Packaging |
|------------|--|-----------------------|
| 23275-410 | ¼ in. x ¾ in. (.6 cm x 1.4 cm), cotton twill | 5/pack 100 packs/case |

Round hard dissectors

Round hard dissectors, also known as beehives, can be used for blunt dissection. They can be grasped with forceps through the small hole in the center (they include a count holder for easy sponge counts).



| BD cat no. | Description | Packaging |
|------------|-----------------------|-----------------------|
| 23275-440 | ¾ in. (.9 cm), cotton | 5/pack 100 packs/case |

Rosebud dissectors

Rosebud dissectors can be used for blunt dissection and are designed for use in delicate surgery (they include a count holder for easy sponge counts).



| BD cat no. | Description | Packaging |
|------------|-----------------------|-----------------------|
| 23275-460 | ¼ in. (.6 cm), cotton | 5/pack 100 packs/case |

Peanut sponges

Peanut sponges can be used for soft tissue dissection and sponging. They can be used as is or folded (they include a count holder for easy sponge counts).



| BD cat no. | Description | Packaging |
|------------|-----------------------|-----------------------|
| 23275-470 | ¾ in. (.9 cm), cotton | 5/pack 100 packs/case |

Cherry sponges

Cherry sponges provide excellent absorbency, and they can be used for general dissection (they include a count holder for easy sponge counts).



| BD cat no. | Description | Packaging |
|------------|------------------------|-----------------------|
| 23275-480 | ½ in. (1.3 cm), cotton | 5/pack 100 packs/case |

Cylindrical sponges

Cylindrical sponges are highly absorbent and designed for use in tonsil, adenoid and nasal surgeries.



| BD cat no. | Description | Packaging |
|------------|---|----------------------|
| 23276-250 | ½ in. x 1 ½ in. (1.3 cm x 3.8 cm), strung, cotton | 5/pack 50 packs/case |

**All dissectors and sponges are x-ray detectable