



# BD Recykleen™

## Sharps Collectors



Helping all people  
live healthy lives

Providing environmentally sound  
solutions for your disposal needs.



# The BD Recykleen™ Story

Being sensitive to environmental issues is important for people and businesses alike. BD recognizes this need and is searching for ways to be a better environmental steward while helping our customers. One result is the BD Recykleen

product line. The amount of recycled material depends upon the product. While BD strives to use as much recycled plastic as possible in each product, there are practical things to consider like:

- **Quality** – Consistency and reliability of recycled resin material
- **Performance** – Products made with recycled resin must meet the same stringent performance requirements as those from virgin plastic
- **Supply** – Availability of recycled resin sufficient to meet the market demand
- **Visual** – Some products need transparency, which cannot be provided with most recycled plastics

This product line will continue to grow and includes:



Collector not included.



**BD Recykleen™  
Cabinets**




















**BD RCRA  
Collectors**

**BD Recykleen™  
L and XL Collectors**

**BD Recykleen™  
Trolleys**

**BD Recykleen™  
Patient Room**

# BD Recykleen™ Collectors

Waste Stream	Type	Ref #	Product Size	Product Dimensions (H x W x D in.)	Qty Case	Compatible with	
Sharps and Infectious Waste	Patient Room	305056	5.4 qt BD Recykleen Pearl Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017, 305447, 305021	
		305057	5.4 qt BD Recykleen Red Counterbalanced Door	10.75 x 12 x 4.5	20		
		305053	3 gal BD Recykleen Pearl Counterbalanced Door	15.75 x 13.5 x 6	10	305097, 305099, 305420 305409, 305588, 305552	
		305054	3 gal BD Recykleen Red Counterbalanced Door	15.75 x 13.5 x 6	10		
							
	Multi-use Nestable	305059	8 qt BD Recykleen Funnel Top	9.75 x 6.75 x 6.75	24	305018, 305485, 305570	
							
	Large and Extra Large	305024	9 gal BD Recykleen Hinge Top	18.75 x 17.75 x 11.75	8	305090	
		305025	9 gal BD Recykleen Slide Top	18.75 x 17.75 x 11.75	8	305090, 305091	
		305023	17 gal BD Recykleen Hinge Top	23.5 x 20 x 14.75	5	305092	
305026		19 gal BD Recykleen Slide Top	26.25 x 20 x 14.75	5	305092, 305093		
		 					
Pharmaceutical	Patient Room	305058	5.4 qt BD Recykleen Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017	
		305055	3 gal BD Recykleen Counterbalanced Door	15.75 x 13.5 x 6	10	305097, 305099, 305420	
		    					
RCRA	Multi-use One Piece	305066	3 gal Plug Cap	12 x 7.5 x 10.5	12	305495	
		305075	3 gal Restrictive Funnel Cap	12 x 7.5 x 10.5	12	305497	
	Large and Extra Large	305068	9 gal Clear Hinge Top	18.5 x 17.75 x 11.75	8	305090	
		305069	9 gal Clear Slide Top	18.5 x 17.75 x 11.75	8	305090, 305091	
		305071	17 gal Clear Hinge Top	23.5 x 19.5 x 14.25	5	305092	
		305072	19 gal Clear Slide Top	26.25 x 19.5 x 14.25	5	305092, 305093	
		      					
Cabinets and Trolleys	Cabinets	305096	5.4 qt BD Recykleen Brown Cabinet	14 x 6 x 15	1	305056, 305057, 305517, 305551, 305543, 305426, 305427, 305428	
		305098	5.4 qt BD Recykleen Green Cabinet	14 x 6 x 15	1		
		305097	2 and 3 gal BD Recykleen Brown Cabinet	15 x 7 x 19	1	305053, 305054, 305422, 305435, 305423, 305436	
		305099	2 and 3 gal BD Recykleen Green Cabinet	15 x 7 x 19	1		
			 				
	Trolleys	305090	9 gal Basic BD Recykleen Trolley	23.5 x 16.3 x 28.6	1	All 9 gal	
		305091	9 gal Foot-operated BD Recykleen Trolley	23.5 x 16.3 x 28.6	1	All 9 gal slide top	
		305092	17 and 19 gal Basic BD Recykleen Trolley	24.1 x 18 x 27.5	1	All 17 and 19 gal	
		305093	19 gal Foot-operated BD Recykleen Trolley	24.3 x 17.7 x 35.2	1	All 19 gal slide top	
							



20%

Product made from a minimum of 20% recycled materials

# BD RECYKLEEN™ COLLECTORS PATIENT ROOM & PHARMACEUTICAL

## Key product features include:

- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
- Counterbalanced door facilitates hands-free disposal and is designed to reduce the potential of sharps injuries
- Designed to improve tamper resistance by preventing hand access to the contents



305054

305057



PATIENT ROOM



305053

305056




305055

305058



PATIENT ROOM



PHARMACEUTICAL  
PATIENT ROOM



25%

Product made from a minimum of 25% recycled materials

# BD RECYKLEEN™ COLLECTORS TROLLEYS



305091



305090



305092



Collectors not included.

305093



305091

## Key product features include:

- A design that reduces the risk of collectors tipping over
- Wheels that enable the collector to be conveniently placed where it will be most used
- Locking wheels, ensuring the product is mobile only when you want it to be
- The ability to open and close the collector for hands-free disposal (foot-operated trolleys only)
- A two-year warranty providing product replacement in case of manufacturer defects



TROLLEY



50%

Product made from a minimum of 50% recycled materials

# BD RECYKLEEN™ COLLECTORS L & XL AND MULTI-USE NESTABLES

## Large and Extra Large

### Key product features include:

- Temporary and final closure features
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to provide security against tampering
- Double wall handles for safe transportation
- Clear tops to ensure visibility of collector fill level
- Easy accessibility for disposal of large and bulky waste
- Dual Opening Hinge Top designed so that the small or the large side can remain open and the other locked based on the size of access opening desired




305026

305024



SLIDE TOP



DUAL OPEN

LARGE AND  
EXTRA LARGE

## Multi-Use Nestables

### Key product features include:

- Clear top designed for better visibility of contents inside the collector
- Funnel cap designed to restrict unwanted access to contents inside the collector can also be used as a guide to direct used sharps inside the collectors
- Needle removal ports designed to remove needle when needed
- Side handles on each side of the collectors to better handle collector when full. Handles folds into the top and stays engaged during use.
- Where storage space is limited, new top and base design allows for stackability after collector is closed for final disposal



305059



Stackable



STACKABLE

MULTI-USE NESTABLE




Product made from a minimum of 70% recycled materials

# BD RECYKLEEN™ COLLECTORS RCRA HAZARDOUS WASTE COLLECTORS

## Key product features of non-DOT Collectors include:

- Clear tops with multiple opening styles to accommodate different types of waste and environments
- Dual Opening Hinge Top designed so that one side can remain open and the other locked based on department needs
- Dual Lock Cap (3 and 5 gallon) designed to provide temporary closure for safe intermediate transportation and secure final closure with two locks
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transportation
- Temporary closure feature and handles designed for safe transportation
- Large and extra large collectors are fully compatible with BD trolleys



RCRA



Product made from 100% recycled materials

# BD RECYKLEEN™ CABINETS

## Key product features include:

- A universal lock providing security and preventing unwanted access to collectors
- A one-piece transparent door facilitating easy cleaning
- The ability to see the collector fill level indicator
- A two-year warranty providing product replacement in case of manufacturer defects



CABINET



305096

# BD ecoFinity™ Life Cycle Solution

Powered by an alliance with Waste Management™

## Introducing the only closed-loop solution for hospital sharps.

The **BD ecoFinity™ Life Cycle Solution** is an innovative sharps recycling program that can recover and recycle **up to 70% of your facility's used sharps materials.\*** Waste management provides safe, convenient on-site collection of medical sharps, plus treatment and recovery of these materials to reduce your environment footprint. BD uses these recycled materials to create clean new containers – closing the loop. Now, single-use can mean both safe and sustainable. **Waste Diverted. No Compromise.**

Visit [bd.com/ecoFinity](http://bd.com/ecoFinity) to learn more.



\*BD data on file. Claim based on pilot studies completed by BD. Services are provided by WM Healthcare Solutions, Inc., a Waste Management company, except in the state of New York where services are provided by Waste Management of New York.



THINK GREEN.®



BD Medical  
[www.bd.com/sharps](http://www.bd.com/sharps)

WM Waste Management logo is a registered service mark of Waste Management, Inc. BD, BD Logo and BD Recykleen are trademarks of Becton, Dickinson and Company. © 2011 BD and WM Healthcare Solutions 0796 Rev 2 (8/11)

