

# BD Disposal Solutions Product Catalog



Helping all people  
live healthy lives

### BD ENVIRONMENTAL STEWARDSHIP

BD's approach to corporate citizenship is embedded in our culture, and is intrinsic to our purpose, *"Helping all people live healthy lives."* Responsible corporate citizenship is integrally tied to our global business strategy. Sustained performance and a dedication to improvement are at BD's core. Recognizing that human health and a healthy environment are inseparable, BD is committed to responsible environmental stewardship in our business conduct.

We are dedicated to the design, manufacture and distribution of safe, effective medical technology products while minimizing environmental impacts.

We strive to conserve energy, water and natural resources and minimize the amount of hazardous waste from our manufacturing and administrative operations. We also promote recycling and resource recovery in our manufacturing processes.

### BD OCEANSIDE MANUFACTURING INNOVATIONS

At BD's manufacturing facility in Oceanside, CA, mechanisms are in place to help minimize the environmental impact of our operations.

Using a lean manufacturing approach towards environmental management, this facility has prioritized to the following areas for manufacturing innovations:

- **Reduce** energy consumption by > 170,000 KW / year
- **Reuse** more than 400,000 lbs of post-industrial plastic resin per year in our products
- **Recycle** > 10,000 lbs / year of paper, plastic, and aluminum generated at the facility
- **Reduce** carbon dioxide emissions by > 85,000 lbs CO<sub>2</sub> / year
- **Reduce** waste oil disposal by > 700 gallons / year


### BD DISPOSAL SOLUTIONS: A COMMITMENT TO QUALITY

Sharps collectors are manufactured using polypropylene or polyethylene resins and colorants that are nontoxic and comply with the Toxics in Packaging Clearing House (TPCH) heavy metal limits. No latex materials are used in these products. The collectors are safe for incineration. Under combustion with adequate access of air, the resin burns to water vapor, carbon dioxide and carbon ash.

BD Sharps Collectors meet or exceed the puncture force requirements defined by the American Society for Testing and Materials (ASTM) and British Standards. Sharps collectors are classified as medical devices and are manufactured according to the FDA Quality System Regulation and the International Quality Systems Standard ISO 9001.

BD Sharps Collectors are listed with the FDA and meet pre-market registration requirements in accordance with 21CFR807. BD Sharps Collectors also meet the Occupational Safety and Health Administration (OSHA) Bloodborne Pathogens Standards 29 CFR 1910.1030.

# Table of Contents



Patient / Exam Room .....	4
Extra Large Sharps .....	5
Multi-Use Nestable .....	6
Multi-Use One Piece .....	7
Phlebotomy .....	7
Chemotherapy .....	8
Pharmaceutical .....	9
RCRA Hazardous Waste Collectors .....	10
BD Recykleen™ Sharps Collectors .....	11
Accessories .....	12
Trolleys .....	13

## Patient / Exam Room

These collectors feature a variety of sizes and colors, counterbalanced doors that are tamper resistant, visible fill levels and handles for safe transportation. The horizontal drop promotes safety and maximizes disposal volume. They provide security in the patient/public accessible areas when used with a locking accessory (see page 12).

### High Sensitivity Counterbalanced Door

- Feature:** Hands-free disposal of sharps with counterbalanced door.
- Advantage:** Hands stay out of the way to fully deposit sharps inside the collector; from 1/2 cc syringe to your largest sharps.
- Benefit:** Designed to reduce potential of sharps injuries by your staff during disposal while counterbalanced door returns to open position without any human intervention.



305017 Cabinet  
with 305551



305423

### Tamper Resistance

- Feature:** Counterbalanced door design prevents unwanted access.
- Advantage:** Eliminates possibility of hands accessing the collector.
- Benefit:** Provides tamper resistance to prevent hand access.

### Handle

- Feature:** Large, ergonomically-designed handle on the front.
- Advantage:** Designed for easy and safe removal.
- Benefit:** One-hand movement for transportation.

### Locking Cabinets with Large Viewing Window

- Feature:** Cabinet with large window for patient room applications.
- Advantage:** Provides extra security in the patient room.
- Benefit:** Designed to prevent unwanted removal of patient room collector; features clear view of fill level.



305443

## Patient / Exam Room

Reference#	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305517	5.4 qt Red Collector / Counterbalanced Door	10.75 x 12 x 4.5	20
305551	5.4 qt Clear Collector / Counterbalanced Door	10.75 x 12 x 4.5	20
305543	5.4 qt Pearl Collector / Counterbalanced Door	10.75 x 12 x 4.5	20
305426	5.4 qt Red Collector / Counterbalanced Door	10.75 x 12 x 4.5	12
305427	5.4 qt Clear Collector / Counterbalanced Door	10.75 x 12 x 4.5	12
305428	5.4 qt Pearl Collector / Counterbalanced Door	10.75 x 12 x 4.5	12
305422	2 gal Pearl Collector / Counterbalanced Door	11.5 x 13.5 x 6	10
305435	2 gal Red Collector / Counterbalanced Door	11.5 x 13.5 x 6	10
305423	3 gal Pearl Collector / Counterbalanced Door	15.75 x 13.5 x 6	10
305436	3 gal Red Collector / Counterbalanced Door	15.75 x 13.5 x 6	10
305443	5.4 qt Red Side Entry	10.75 x 12 x 4.5	20
305444	5.4 qt Pearl Side Entry	10.75 x 12 x 4.5	20
305425	5.4 qt Pearl Side Entry	10.75 x 12 x 4.5	12

See Accessories on page 12

## Extra Large Sharps

With large access openings, including temporary and final closure, these collectors feature dual opening hinge and sliding tops with VAT (visual • audible • tactile) perimeter locks, temporary closure, handles and clear tops to view fill level. They are sturdy, stackable when full and designed for easy handling and disposal.

### Large Opening Slide Top (9/19 gal) sizes

**Feature:** 5.25" x 7.25" (9 gal), 6" x 9" (19 gal) opening.

**Advantage:** Disposal of large amount of sharps waste.

**Benefit:** Provides easy accessibility to the collector when sharps waste is large and cumbersome.

### VAT Perimeter Locks visual • audible • tactile

**Feature:** 18 (9 gal) or 20 (17 & 19 gal) VAT perimeter locks to fully secure collector during assembly.

**Advantage:** Provides added security against tampering.

**Benefit:** VAT perimeter locks add security against tampering of sharps waste contents.



Dual Opening Hinge Top

### Dual Opening Hinge Top

**Feature:** Dual opening hinge top.

**Advantage:** Decide the side to remain opened based on your department needs. Locked side remains permanently locked with the help of the VAT perimeter locks.

**Benefit:** Accessibility on both sides of collector based on sharps waste generator needs and added security against penetration with VAT perimeter locks.

## Extra Large Sharps

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305615	9 gal Clear Hinge Top	18.75 x 17.75 x 11.75	8
305616	9 gal Clear Slide Top	18.75 x 17.75 x 11.75	8
305601*	9 gal Clear Hinge Top with gasket	18.75 x 17.75 x 11.75	8
305602*	9 gal Clear Slide Top with gasket	18.75 x 17.75 x 11.75	8
305610	17 gal Clear Hinge Top	23.5 x 20 x 14.75	5
305665*	17 gal Clear Hinge Top with gasket	23.5 x 20 x 14.75	5
305609	19 gal Clear Slide Top	26.25 x 20 x 14.75	5
305666*	19 gal Clear Slide Top with gasket	26.25 x 20 x 14.75	5

\* Meets or exceeds **US DOT** packaging requirements 49CFR173.24, 49CFR173.24a, and 49CFR173.134.

See Trolleys on page 13



305567 Trolley  
with 305616

## Multi-Use Nestable

These collectors feature a variety of sizes, clear lids with multiple opening styles, temporary closure and handles designed for safe transporting. The funnel top is designed to provide extra security by restricting unwanted access.

### Tamper Resistance – Funnel

**Feature:** Funnel designed to prevent unwanted access.

**Advantage:** Designed to reduce the possibility of hands accessing the contents.

**Benefit:** Designed to promote safe disposal of sharps waste.



305480



Regular Funnel Clear Top



Large Funnel Clear Top



Open Clear Top (8 qt)



Small Open Clear Top



Open Clear Top (6 gal)



305457

305460



305485 Bracket with 305646

### Open Top

**Feature:** 2" diameter opening — 8 qt  
4.25" diameter opening — 14 qt  
2.5" diameter opening — 14 qt small opening  
4.25" diameter opening — 6 gal

**Advantage:** Large opening to accommodate large disposal needs in non-patient/public areas.

**Benefit:** Disposal of large variety of sharps waste.



305417 with 305570 Stabilizer

## Multi-Use Nestable

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305417	8 qt Open Clear Top	11.25 x 10 x 6.75	24
305460	8 qt Regular Funnel Clear Top	11.25 x 10 x 6.75	24
305463	8 qt Pearl Regular Funnel Clear Top	11.25 x 10 x 6.75	24
305456	14 qt Open Clear Top	11.5 x 12.5 x 8.5	20
305464	14 qt Regular Funnel Clear Top	11.5 x 12.5 x 8.5	20
305480	14 qt Large Funnel Clear Top	11.5 x 12.5 x 8.5	20
305646	14 qt Small Open Clear Top	11.5 x 12.5 x 8.5	20
305457	6 gal Open Clear Top	17.5 x 12.5 x 8.5	12
305465	6 gal Large Funnel Clear Top	17.5 x 12.5 x 8.5	12

See Accessories on page 12

## Multi-Use One Piece



305477

305490

These collectors are preassembled and designed to meet the needs of a variety of clinical settings. They feature temporary closure and final closure, handles and multiple opening styles.

### Hinge Cap 6.9 qt and 8.2 qt

- Feature:** Hinge cap for the 6.9 and 8.2 collectors.
- Advantage:** Provides temporary closure for in-between usage and intermediate transportation.
- Benefit:** Designed to provide temporary closure for safe intermediate transportation and secure closure for final disposal.



305489

305488

305557

## Multi-Use One Piece

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case
305557	1.4 qt Regular Funnel	7.75 x 3.75 x 3.75	36
305488	3.3 qt Regular Funnel	5.75 x 8.5 x 5	24
305489	6.9 qt Regular Funnel	11.5 x 8.75 x 5.5	12
300467*	6.9 qt Regular Funnel	11.5 x 8.75 x 5.5	12
305490	8.2 qt Regular Funnel	13.5 x 8.75 x 5.5	12
300470*	8.2 qt Regular Funnel	13.5 x 8.75 x 5.5	12
305491	5 gal Open Top	18 x 10.5 x 7.5	8
300473*	5 gal Open Top	18 x 10.5 x 7.5	8
305477	5 gal Large Funnel	18 x 10.5 x 7.5	8
305577	5 gal Screw Cap	18 x 10.5 x 7.5	8

\* Cannot be autoclaved and not vented

See Accessories on page 12

### Handle

- Feature:** Ergonomically designed handle.
- Advantage:** Designed for comfort and safety.
- Benefit:** Designed for ease of movement and transportation.

### Large Opening 5 gal

- Feature:** 4.5" diameter opening.
- Advantage:** Large opening to accommodate large disposal needs in non-patient/public areas.
- Benefit:** Disposal of large variety of sharps waste.

## Phlebotomy

These collectors feature clear tops, temporary closures, needle removal ports and a trim design for trays and bags.



305635

305487

### Lid Design – Needle Removal

- Feature:** Built-in needle removal port.
- Advantage:** Designed to easily remove needle from syringes and blood collection devices.
- Benefit:** Helps prevent contact with sharps for safe disposal.



305970 Bracket  
with 305487

## Phlebotomy

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case
305635	1 qt Small Open Top	6.5 x 4.75 x 4	60
305487	1.5 qt Removal Port	9 x 4.75 x 4	36



305635

305487

See Accessories on page 12

# Chemotherapy

These collectors feature a variety of lid designs and sizes for disposal of trace chemotherapy waste. Trace waste is defined by CDC/NIOSH as “those that contain less than 3% by weight of the original quantity of hazardous drugs.”<sup>1</sup> The collectors are designed to be leak resistant and easily identified by color and label. As an added feature, each chemotherapy collector comes with a liquid-absorbent material pre-attached to the inside bottom of the collector.

1 CDC/NIOSH Alert Publication No. 2004-165

## Large Opening Slide Top (9/19 gal) sizes

- Feature:** 5.25" x 7.25" (9 gal), 6" x 9" (19 gal) opening.
- Advantage:** Disposal of large amount of chemotherapy waste.
- Benefit:** Provides easy accessibility to the collector when chemotherapy waste is large and cumbersome.



305611 Trolley with 305614



## Hinge and Slide Top with Inner Gasket Lining

- Feature:** Hinge and slide top with inner gasket lining.
- Advantage:** Minimizes leak of unwanted liquid if accidentally tipped over.
- Benefit:** Hinge and slide top provides one-hand operation to open/close to gain access to collector. The inner gasket lining potentially reduces leaks and adds security during final closure.



Dual Opening Hinge Top

## Chemotherapy

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305492	3 gal Plug Cap	12 x 7.5 x 10.5	12
305493	5 gal Plug Cap	18 x 10.5 x 7.5	8
305578	5 gal One Piece Screw Cap	18 x 10.5 x 7.5	8
305603	9 gal Clear Hinge Top	18.75 x 17.75 x 11.75	8
305604	9 gal Clear Slide Top	18.75 x 17.75 x 11.75	8
305614	17 gal Clear Hinge Top	23.5 x 20 x 14.75	5
305613	19 gal Clear Slide Top	26.25 x 20 x 14.75	5



Low temperature incinerator

# Pharmaceutical

These collectors feature a variety of lid designs and sizes to collect and dispose unused, expired, or older pharmaceutical waste. They are designed to be leak resistant and easily identified by color and label. As an added feature, each collector comes with a liquid-absorbent material pre-attached to the inside bottom of the collector.



305028



305329

## VAT Perimeter Locks visual • audible • tactile

**Feature:** 18 (9 gal) or 20 (17 & 19 gal) VAT perimeter locks to fully secure collector during assembly.

**Advantage:** Provides added security against tampering.

**Benefit:** VAT perimeter locks add security against tampering of pharmaceutical waste contents.



305420 Cabinet with 305622

## High Sensitivity Counterbalanced Door

**Feature:** Hands-free disposal of pharmaceutical waste.

**Advantage:** Hands stay out of the way to fully deposit waste inside the collector.

**Benefit:** Designed to reduce potential of injuries by your staff during disposal while counterbalanced door returns to open position without any human intervention.



305622



305633

## Handle

**Feature:** Large, ergonomically-designed handle.

**Advantage:** Designed for easy and safe removal.

**Benefit:** One-hand movement for transportation.

## Dual Opening Hinge Top

**Feature:** Dual opening hinge top.

**Advantage:** Decide the side to remain opened based on your department needs. Locked side remains permanently locked with the help of the VAT perimeter locks.

**Benefit:** Accessibility on both sides of collector based on pharmacy waste generator needs and added security against penetration with VAT perimeter locks.



Dual Opening Hinge Top

## Pharmaceutical

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305330	5.4 qt Counterbalanced Door	10.75 x 12 x 4.5	20
305622	3 gal Counterbalanced Door	15.75 x 13.5 x 6	10
305633	3 gal Plug Cap	12 x 7.5 x 10.5	12
305634	9 gal Hinge Top	18.75 x 17.75 x 11.75	8
305329	9 gal Slide Top	18.75 x 17.75 x 11.75	8
305028	17 gal Hinge Top	23.5 x 20 x 14.75	5



Low temperature incinerator

These collectors are made from a minimum of 70% recycled material and are intended for Resource Conservation and Recovery Act (RCRA) solid hazardous waste. As an added feature, each RCRA collector comes with a liquid-absorbent material pre-attached to the inside bottom of the collector.



## Two-Lock Cap — 3 gal and 5 gal

**Feature:** Two-lock hinge system.

**Advantage:** Provides temporary closure for in-between usage and intermediate transportation.

**Benefit:** Designed to provide temporary closure for safe intermediate transportation and secure final closure with two locks.

Two-lock cap examples



## Custom-fit Top Foam Gasket — 9, 17 and 19 gal

**Feature:** Proprietary custom fit process to form foam gasket around the inside of the hinge top.

**Advantage:** Allows adaptability to hinge top shape regardless of top size.

**Benefit:** Minimizes leaks and adds security during final closure.



Dual Opening Hinge Top

## RCRA Hazardous Waste Collectors

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305066	3 gal One Piece	12 x 7.5 x 10.5	12
305067	5 gal One Piece	18 x 10.5 x 7.5	8
305068	9 gal Clear Hinge Top	18.5 x 17.75 x 11.75	8
305069	9 gal Clear Slide Top	18.5 x 17.75 x 11.75	8
305071	17 gal Clear Hinge Top	23.5 x 19.5 x 14.25	5
305072	19 gal Clear Slide Top	26.25 x 19.5 x 14.25	5



## BD introduces BD Recykleen™ Sharps Collectors.

This is a new sharps collector line made with a minimum of 50% recycled material. These collectors feature dual opening hinge top with VAT perimeter locks, temporary closure, handles and clear top to view fill levels.



Dual Opening Hinge Top

### Dual Opening Hinge Top with VAT Perimeter Locks visual • audible • tactile

- Feature:** Dual opening hinge top with 18 (9 gal) or 20 (17 gal) VAT perimeter locks.
- Advantage:** Decide the side to remain opened based on your department needs. Locked side remains permanently locked with the help of the VAT perimeter locks.
- Benefit:** Accessibility on both sides of collector based on sharps waste generator needs and added security against penetration with VAT perimeter locks.

## BD Recykleen™ Sharps Collectors

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case
305024	9 gal Hinge Top	18.75 x 17.75 x 11.75	8
305023	17 gal Hinge Top	23.5 x 20 x 14.75	5

# Accessories

## Cabinets, Brackets and Stabilizers

Category	Reference Number	Description	Dimensions (H x W x D inches)	Qty/Case	Use With	Compatible with Products
Patient/ Exam Room	305017	Cabinet†	13.5 X 13 X 5.5	1	5.4 qt Counterbalance Door	305517, 305551, 305543, 305426, 305427, 305428
	305447	Locking Bracket	9.5 x 9.5 x 4.25	10	All 5.4 qt	305517, 305551, 305543, 305426, 305427, 305428, 305443, 305444, 305425
	305021	Stabilizer	3 x 13 x 6.5	5	All 5.4 qt	305517, 305551, 305543, 305426, 305427, 305428, 305443, 305444, 305425
	305420	Cabinet	18.5 x 14.5 x 6.75	1	2/3 gal Counterbalance Door	305422, 305435, 305423, 305436
	305409	Locking Bracket	4.25 x 14 x 4.75	10	2/3 gal Counterbalance Door	305422, 305435, 305423, 305436
	305588	Wire Bracket	6 x 12.5 x 6	6	2/3 gal Counterbalance Door	305422, 305435, 305423, 305436
	305552	Stabilizer	3 x 13.5 x 7.5	5	2/3 gal Counterbalance Door	305422, 305435, 305423, 305436
	305448	Glove Box	7 x 12 x 4.25	6	Attaches to cabinets and 5.4 qt locking bracket	Attached to 305017, 305420, 305447
Multi-Use Nestable	305485	Bracket	6 x 3.5 x 1.25	12	Multi-Use Nestable 8 qt, 14 qt & 6 gal	305417, 305460, 305463, 305456, 305464, 305480, 305646, 305457, 305465
	305589	Wire Bracket	9.5 x 10.25 x 7.75	4	Multi-Use Nestable 8 qt	305417, 305460, 305463
	305570	Stabilizer	2 x 10 x 7	5	Multi-Use Nestable 8 qt	305417, 305460, 305463
	305494	Padlock	1.75 x 1 x .25	12	Bracket, Wire Bracket	Use with 305485, 305589, 305570 to make locking
Multi-Use One Piece	305497	Quick Release Strap	1.5 x 8 x 1	12	Multi-Use One Piece 3.3 qt, 6.9 qt, 8.2 qt, 3 gal (Plug Cap), 5 gal	305488, 305489, 300467, 305490, 300470, 305491, 300473, 305477, 305577
	305496	Security Bracket	1.5 x 8 x 8	12	Multi-Use One Piece 6.9 qt, 8.2 qt	305489, 300467, 305490, 300470
	305495	XL Security Bracket	2 x 10.5 x 8	12	5 gal Multi-Use One Piece	305491, 300473, 305477, 305577
	305494	Padlock	1.75 x 1 x .25	12	Security Bracket, XL Security Bracket	Use with 305496, 305495 to make locking
Phlebotomy	305970	Cup Bracket	3 x 3.5 x 3	12	1 qt, 1.5 qt	305635, 305487
Chemotherapy	305497	Quick Release Strap	1.5 x 8 x 1	12	3 gal (Plug Cap), 5 gal	305492, 305493, 305578
	305495	XL Security Bracket	2 x 10.5 x 8	12	3 gal (Plug Cap), 5 gal	305492, 305493, 305578
	305494	Padlock	1.75 x 1 x .25	12	XL Security Bracket	Use with 305495
Pharmaceutical	305017	Cabinet	13.5 X 13 X 5.5	1	5.4 qt Counterbalance Door	305330
	305447	Locking Bracket	9.5 x 9.5 x 4.25	10	All 5.4 qt	305330
	305021	Stabilizer	3 x 13 x 6.5	5	All 5.4 qt	305330
	305420	Cabinet	18.5 x 14.5 x 6.75	1	2/3 gal Counterbalance Door	305622
	305409	Locking Bracket	4.25 x 14 x 4.75	10	2/3 gal Counterbalance Door	305622
	305588	Wire Bracket	6 x 12.5 x 6	6	2/3 gal Counterbalance Door	305622
	305552	Stabilizer	3 x 13.5 x 7.5	5	2/3 gal Counterbalance Door	305622
	305448	Glove Box	7 x 12 x 4.25	6	Can be attached to 305017, 305420 & 305447	Can be attached to 305017, 305420, 305447
	305497	Quick Release Strap	1.5 x 8 x 1	12	3 gal (Plug Cap)	305633
	305495	XL Security Bracket	2 x 10.5 x 8	12	3 gal (Plug Cap)	305633
305494	Padlock	1.75 x 1 x .25	12	XL Security Bracket	Use with 305495	
RCRA	305497	Quick Release Strap	1.5 x 8 x 1	12	3 gal (Plug Cap), 5 gal	305066, 305067
	305495	XL Security Bracket	2 x 10.5 x 8	12	3 gal (Plug Cap), 5 gal	305066, 305067
	305494	Padlock	1.75 x 1 x .25	12	XL Security Bracket, Quick Release Strap	Use with 305495

† Cannot be used with Side Entry 5.4 qt

# Trolleys

For the extra large collectors, standard trolleys and foot-operated trolleys are available. These trolleys provide stability and mobility, with handles and 3-inch locking wheels. The foot-operated trolley has the ability to open and close the collector via foot-operated pedal for hands-free disposal.

305567 Trolley with 305329



305556 Trolley with 305024



305612 Trolley with 305609



305567 Trolley with 305069


305611 Trolley with 305614



## Trolleys

Category	Reference Number	Description	Dimensions (H x W x D inches)	Qty/Case	Use With	Compatible with Products
Extra Large	305556	Standard Trolley	24 x 15.5 x 22	1	9 gal Hinge/Slide Top	305615, 305616, 305601, 305602
	305567	Foot-Operated Trolley	26.5 x 15.5 x 22.25	1	9 gal Slide Top	305616, 305602
	305611	Standard Trolley	26.5 x 17 x 23	1	17 gal Hinge Top & 19 gal Slide Top	305610, 305609, 305665, 305666
	305612	Foot-Operated Trolley	34 x 17 x 23	1	19 gal Slide Top	305609, 305666
Chemotherapy	305556	Standard Trolley	24 x 15.5 x 22	1	9 gal Hinge/Slide Top	305603, 305604
	305567	Foot-Operated Trolley	26.5 x 15.5 x 22.25	1	9 gal Slide Top	305604
	305611	Standard Trolley	26.5 x 17 x 23	1	17 gal Hinge Top & 19 gal Slide Top	305614, 305613
	305612	Foot-Operated Trolley	34 x 17 x 23	1	19 gal Slide Top	305613
Pharmaceutical	305556	Standard Trolley	24 x 15.5 x 22	1	9 gal Hinge/Slide Top	305634, 305329
	305567	Foot-Operated Trolley	26.5 x 15.5 x 22.25	1	9 gal Slide Top	305329
	305611	Standard Trolley	26.5 x 17 x 23	1	17 gal Hinge Top	305028
RCRA	305556	Standard Trolley	24 x 15.5 x 22	1	9 gal Hinge/Slide Top	305068, 305069
	305567	Foot-Operated Trolley	26.5 x 15.5 x 22.25	1	9 gal Slide Top	305069
	305611	Standard Trolley	26.5 x 17 x 23	1	17 gal Hinge Top & 19 gal Slide Top	305071, 305072
	305612	Foot-Operated Trolley	34 x 17 x 23	1	19 gal Slide Top	305072
BD Recycleen™	305556	Standard Trolley	24 x 15.5 x 22	1	9 gal Hinge/Slide Top	305024
	305611	Standard Trolley	26.5 x 17 x 23	1	17 gal Hinge Top	305023

# WASTE STREAM



Segregation of dedicated pharmaceutical waste versus other regulated medical waste is an important safety and environmental activity that starts with you. Selecting the appropriate collector for your waste segregation needs help minimize the impact on the environment and provide proper waste disposal.

Contact Customer Service  
at **1-888-237-2762**  
or your BD Sales Representative  
for additional information



**BD Medical**  
4665 North Avenue  
Oceanside, CA 92056  
[www.bd.com/sharps](http://www.bd.com/sharps)