

BD Disposal Solutions

Product Catalog

Recognizing that human health and
a healthy environment are inseparable



Sharps

Pharmaceutical

Chemotherapy

RCRA

Accessories



Helping all people
live healthy lives

The BD Recykleen™ Story

Being sensitive to environmental issues is very important for all people and businesses. BD recognizes this need and is searching for ways to be a better environmental steward while helping our customers. One result is the BD Recykleen™ product line.

These products are made from post-industrial or post-consumer recycled plastic. The amount of recycled material depends upon the product. While BD strives to use as much recycled plastic as possible in each product, there are practical things to consider like:

- **Supply** – Availability of recycled resin sufficient to meet the market demand
- **Quality** – Consistency and reliability of recycled resin material
- **Performance** – Products made with recycled resin must meet the same stringent performance requirements as those from virgin plastic
- **Visual** – Some products need transparency, which cannot be provided with most recycled plastics

This product line will continue to grow and includes:



Collector not included.



**BD Recykleen™
Cabinets**

**BD RCRA
Collectors**

**BD Recykleen™
L and XL Collectors**





**BD Recykleen™
Trolleys**

**BD Recykleen™
Patient Room Collectors**

Table of Contents

Waste Stream Segregation	4
How to Choose the Right Collector	5
Patient Room Sharps Collectors	6
Large and Extra-Large Sharps Collectors	7
Multi-use Nestable Sharps Collectors	8
Multi-use One Piece Sharps Collectors	9
Phlebotomy Sharps Collectors.....	9
Pharmaceutical Sharps Collectors.....	10
Chemotherapy Sharps Collectors.....	11
RCRA Hazardous Waste Collectors—DOT	12
RCRA Hazardous Waste Collectors—Non-DOT	13
Cabinets and Accessories.....	14
Trolleys	15

Look for these symbols throughout the catalog to easily identify key features and benefits:




Waste Stream Segregation

BD Disposal Solutions offers collectors to meet the disposal needs of the four categories of waste generated in medical settings. There are two main categories of waste: regulated medical waste and hazardous waste. Regulated medical waste includes Sharps / Infectious Waste, Pharmaceutical (non-hazardous) and Chemotherapy waste.

The hazardous waste category is for the collection of RCRA (Resource Conservation and Recovery Act by the U.S. EPA) waste.

This table indicates general guidelines for proper waste stream segregation.



How to Choose the Right Collector

Waste Stream	Type	Capacity	Available in BD Recykleen™	Access Opening Size (inches)	Sample Products that Fit in Collector
Sharps	Patient Room	5.4 qt with Counterbalanced Door	●	2.0 x 8.0	Blood collection tubes and collection assembly, syringes up to 60 mL with needle attached
		2 gal with Counterbalanced Door		2.0 x 8.0	
		3 gal with Counterbalanced Door	●	2.0 x 8.0	
	Phlebotomy	1 qt		2.5 x 1.2 Open (305635)	1 mL - 10 mL syringes, blood collection tubes and assembly, slides, pipettes
		1.5 qt		1.5 x 1.5 Open (305487)	
	Multi-use Nestable	8 qt		2.5 x 1.75 Dia Open (305343 only); 1.4 x 1.0 Dia Petals (305344, 305059)	1 mL - 20 mL syringes with needle, blood collection tubes, slides, pipettes
		14 qt		2.5 Dia Open (305646); 4.22 Dia Petals (305456); 1.8 x 1.1 Dia Petals (305464); 1.0 Dia Petals (305480)	1 mL - 60 mL syringes with needle, blood collection assembly, blood collection tubes, slides, pipettes
		6 gal		4.22 Dia Open (305457); 1.0 Dia Petals (305465)	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
	Multi-use One Piece	1.4 qt		1.8 Dia Petals	Small syringes, needles, pipettes, slides
		3.3 qt		1.8 Dia Petals	
		6.9 qt		1.8 Dia Petals	
		8.2 qt		1.8 Dia Petals	
		5 gal		4.2 Dia Open	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
	Large and Extra Large	9 gal Hinge Top	●	9.0 x 10.0	All items mentioned above plus small laparoscopic instruments, phlebotomy procedure bottles, dialysis catheters
		9 gal Slide Top	●	5.0 x 7.5	
		17 gal Hinge Top	●	10.5 x 12.5	Medical devices, spinal tap devices, laparoscopy instrumentation
		19 gal Slide Top	●	6.0 x 8.5	
Non-Hazardous Pharmaceutical	Patient Room	5.4 qt with Counterbalanced Door	●	2.0 x 8.0	All syringe sizes with needle attached, small IV bags 50 mL - 250 mL, medication vials, ampules
		3 gal with Counterbalanced Door	●	2.0 x 8.0	All syringe sizes with needle attached, small IV bags 50 mL - 250 mL
	Multi-use One Piece	3 gal		4.2 Dia Open (305633)	All syringe sizes and bags up to 250 mL
	Large and Extra Large	9 gal Hinge Top		9.0 x 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
		9 gal Slide Top		5.0 x 7.5	
		17 gal Hinge Top		10.5 x 12.5	
		19 gal Slide Top		6.0 x 8.5	
Trace Chemotherapy	Multi-use One Piece	3 gal		4.22 Dia Open (305076)	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
		5 gal		4.2 Dia Open (305493)	
		5 gal Screw Cap		4.5 Dia Open (305578)	
	Large and Extra Large	9 gal Hinge Top		9.0 x 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
		9 gal Slide Top		5.0 x 7.5	
		17 gal Hinge Top		10.5 x 12.5	
		19 gal Slide Top		6.0 x 8.5	
Hazardous RCRA	Multi-use One Piece	3 gal	●	4.22 Dia Open (305066) 1.0 Dia Petal (305075)	Smaller size of solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		5 gal		4.2 Dia Open (305067)	
	Large and Extra Large	9 gal Hinge Top	●	9.0 x 10.0	Solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		9 gal Slide Top	●	5.0 x 7.5	
		17 gal Hinge Top	●	10.5 x 12.5	
		19 gal Slide Top	●	6.0 x 8.5	
	PGII	13 gal		6.0 Dia Open	Solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		16 gal		6.0 Dia Open	

Patient Room Sharps Collectors

BD Patient Room Collectors are available in BD Recykleen™ and BD standard brands and come in a variety of sizes and colors. All BD Patient Room Collectors have visible fill-level indicators and handles for safe transportation.

Key product features include:

- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
- Counterbalanced door facilitates hands-free disposal and is designed to reduce the potential of sharps injuries
- Designed to improve tamper resistance by preventing hand access to the contents
- Available in BD Recykleen™ brand



305517



305053



20%

Product made from a minimum of 20% recycled materials

Patient / Exam Room


Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
BD Recykleen™ Sharps Collectors				
305056	5.4 qt BD Recykleen™ Counterbalanced Door Pearl	12 x 12 x 4.8	20	305096, 305098, 305447, 305021
305057	5.4 qt BD Recykleen™ Counterbalanced Door Red	12 x 12 x 4.8	20	
305053	3 gal BD Recykleen™ Counterbalanced Door Pearl	16.6 x 10.7 x 6	10	305097, 305099, 305409, 305588, 305552
305054	3 gal BD Recykleen™ Counterbalanced Door Red	16.6 x 10.7 x 6	10	
305517	5.4 qt Counterbalanced Door Red	12 x 12 x 4.8	20	305096, 305098, 305447, 305021
305551	5.4 qt Counterbalanced Door Clear	12 x 12 x 4.8	20	
305543	5.4 qt Counterbalanced Door Pearl	12 x 12 x 4.8	20	
305426	5.4 qt Counterbalanced Door Red	12 x 12 x 4.8	12	
305427	5.4 qt Counterbalanced Door Clear	12 x 12 x 4.8	12	
305428	5.4 qt Counterbalanced Door Pearl	12 x 12 x 4.8	12	
305422	2 gal Counterbalanced Door Pearl	12.45 x 10.72 x 6	10	
305435	2 gal Counterbalanced Door Red	12.45 x 10.72 x 6	10	305097, 305099, 305409, 305588, 305552
305423	3 gal Counterbalanced Door Pearl	16.6 x 10.7 x 6	10	
305436	3 gal Counterbalanced Door Red	16.6 x 10.7 x 6	10	305447, 305021
305443	5.4 qt Counterbalanced Side Entry Red	11.7 x 11.6 x 4.5	20	
305444	5.4 qt Counterbalanced Side Entry Pearl	11.7 x 11.6 x 4.5	20	
305425	5.4 qt Counterbalanced Side Entry Pearl	11.7 x 11.6 x 4.5	12	

Large and Extra-Large Sharps Collectors

BD Large and Extra Large Collectors are available in BD Recykleen™ and BD standard brands. These collectors are fully compatible with trolleys to provide stability and mobility (see page 15 for details).



305023



Key product features include:

- Temporary and final closure features
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to provide security against tampering
- Double-wall handles for safe transportation
- Clear tops to ensure visibility of collector fill level
- Easy accessibility for disposal of large and bulky waste
- Dual-opening hinge top designed so that the small or the large side can remain open and the other locked based on the size of access opening desired



305616

Large and Extra-Large Sharps

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
BD Recykleen™ Sharps Collectors				
305024	9 gal BD Recykleen™ Hinge Top	18.5 x 17.75 x 11.75	8	305090
305025	9 gal BD Recykleen™ Slide Top	18.5 x 17.75 x 11.75	8	305090, 305091
305023	17 gal BD Recykleen™ Hinge Top	23.5 x 20 x 14.75	5	305092
305026	19 gal BD Recykleen™ Slide Top	26.25 x 20 x 14.75	5	305092, 305093
305615	9 gal Hinge Top	18.5 x 17.75 x 11.75	8	305090
305601*	9 gal Hinge Top with Gasket	18.5 x 17.75 x 11.75	8	
305616	9 gal Slide Top	18.5 x 17.75 x 11.75	8	305090, 305091
305602*	9 gal Slide Top with Gasket	18.5 x 17.75 x 11.75	8	
305610	17 gal Hinge Top	23.5 x 20 x 14.75	5	305092
305665*	17 gal Hinge Top with Gasket	23.5 x 20 x 14.75	5	
305609	19 gal Slide Top	26.25 x 20 x 14.75	5	305092, 305093
305666*	19 gal Slide Top with Gasket	26.25 x 20 x 14.75	5	

* Meets or exceeds **US DOT** packaging requirements 49 CFR 173.24, 49 CFR 173.24a, and 49 CFR 173.134.

See Trolleys on page 15

Multi-use Nestable Sharps Collectors

BD Multi-use Nestable Collectors are available in three sizes and can be used with accessories for added security and stability.

Key product features include:

- Clear tops with multiple opening styles to accommodate different sizes of waste and environments
- Handles designed for safe transport
- Funnel-style lids designed to prevent accessibility to the collector contents and minimize spillage if the collector tips over
- Open-top lids feature larger diameter openings to accommodate disposal needs in non-patient / public areas



305457

8 qt



Regular Funnel Clear Top
Hinge Cap with Petals

8 qt



Open Clear Top
Hinge Cap

14 qt and 6 gal



Large Funnel Clear Top
Hinge Cap with Petals

14 qt



Open Clear Top
with Tethered Cap



305343 with
305570 Stabilizer



305646 with
305485 Bracket



Multi-use Nestable

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305343	8 qt Hinge Cap	10.3 x 11.3 x 6.8	24	305485, 305018, 305570
305344	8 qt Hinge Cap with Petals	10.3 x 11.3 x 6.8	24	
305646	14 qt Tethered Cap	11.5 x 12.8 x 8.8	20	305485
305456	14 qt Hinge Cap	11.5 x 12.8 x 8.8	20	
305464	14 qt Hinge Cap with Needle Ports	11.5 x 12.8 x 8.8	20	
305480	14 qt Hinge Cap with Petals	11.5 x 12.8 x 8.8	20	
305457	6 gal Hinge Cap	17.5 x 12.8 x 8.8	20	
305465	6 gal Hinge Cap with Petals	17.5 x 12.8 x 8.8	12	

Brackets and Stabilizers

305485	Bracket	6 x 3.5 x 1.25	12	8 qt, 14 qt, 6 gal
305018	Wire Bracket	9.5 x 10.25 x 7.75	6	All 8 qt
305570	Stabilizer	2 x 10 x 7	5	All 8 qt
305494	Security Lock	1.75 x 1 x .25	12	305485, 305018, 305570

See Accessories on page 14

Multi-use One Piece Sharps Collectors

BD Multi-use One Piece Collectors are available in five sizes and can be used with accessories for added security.

Key product features include:

- Handles designed for safe movement and transportation
- Hinge caps designed to provide temporary and final closure features



305557



305488

Multi-use One Piece

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305557	1.4 qt Hinge Cap with Petals	7.5 x 3.6 x 3.6	36	
305488	3.3 qt Hinge Cap with Petals	5 x 9.1 x 5.3	24	305497
305489	6.9 qt Hinge Cap with Petals	11.5 x 9.4 x 5.3	12	
300467*	6.9 qt Hinge Cap with Petals	11.5 x 8.75 x 5.5	12	305496
305490	8.2 qt Hinge Cap with Petals	13.4 x 9.4 x 5.3	12	305497
300470*	8.2 qt Hinge Cap with Petals	13.4 x 9.4 x 5.3	12	
305491	5 gal Hinge Cap	18 x 10.5 x 7.5	8	
300473*	5 gal Hinge Cap	18 x 10.5 x 7.5	8	
305477*	5 gal Hinge Cap with Petals	18 x 10.5 x 7.5	8	305495
305577	5 gal Screw Cap	18 x 10.5 x 7.5	8	305497
305100	5 gal Tethered Cap	18 x 10.5 x 7.5	8	
305078	6 gal Tethered Cap	17.5 x 12.5 x 8.5	12	305485

Brackets and Stabilizers

305497	Quick Release Strap	1.5 x 8 x 1	12	3.3 qt, 6.9 qt, 8.2 qt, 5 gal
305496	Security Bracket	1.5 x 8 x 8	12	6.9 qt, 8.2 qt
305495	XL Security Bracket	2 x 10.5 x 8	12	5 gal
305494	Padlock	1.75 x 1 x .25	12	305495, 305496

* Cannot be autoclaved and is not vented

Phlebotomy Sharps Collectors

BD Phlebotomy Collectors are available in two sizes and can be used with brackets for added stability.

Key product features include:

- Built-in standard, safety and pen needle removal ports designed to prevent contact with sharps during disposal
- Clear tops featuring a clear view of the fill level
- Temporary and final closures for safety and security



305635



305487 with
305970 Bracket

Phlebotomy

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305635	1 qt with Needle Ports	7 x 4.9 x 3.9	60	305970
305487	1.5 qt Needle Port	9.4 x 4.9 x 3.9	36	305970

Brackets

305970	Cup Bracket	3 x 3.5 x 3	12	305635, 305487
--------	-------------	-------------	----	----------------

Pharmaceutical Sharps Collectors

BD Pharmaceutical Collectors feature a variety of lid designs and sizes to collect and dispose of unused, expired or other pharmaceutical waste (see Waste Stream Segregation on page 4).

Key product features include:

- Gaskets or liquid-absorbing materials to prevent leakage
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transportation
- Dual-opening hinge top designed so that the small or the large side can remain open and the other locked based on the size of access opening desired
- Large and extra-large collectors are fully compatible with BD Trolleys (see page 15)
- 5.4 quart and 3 gallon Pharmaceutical Collectors with counterbalanced door:
 - Facilitate hands-free disposal and are designed to reduce the potential for sharps injuries
 - Designed to improve tamper resistance by preventing hand access to the contents
 - Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
 - Available in BD Recykleen™ brand



305055





305633



305028



305329



Pharmaceutical

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
BD Recykleen™ Sharps Collectors				
305058	5.4 qt BD Recykleen™ Counterbalanced Door	12 x 12 x 4.8	20	305096, 305098, 305021, 305447
305055	3 gal BD Recykleen™ Counterbalanced Door	16.6 x 10.7 x 6	10	305097, 305099, 305409, 305552, 305588
305330	5.4 qt Counterbalanced Door	12 x 12 x 4.8	20	305096, 305098, 305021, 305447
305622	3 gal Counterbalanced Door	16.6 x 10.7 x 6	10	305097, 305099, 305409, 305552, 305588
305633	3 gal Hinge Cap	12 x 10.5 x 7.5	12	305495, 305497
305634	9 gal Hinge Top with Gasket	18.5 x 17.75 x 11.75	8	305090
305329	9 gal Slide Top with Gasket	18.5 x 17.75 x 11.75	8	305090, 305091
305028	17 gal Hinge Top with Gasket	23.5 x 20 x 14.75	5	305092
305328	19 gal Slide Top with Gasket	26.25 x 20 x 14.75	5	305092, 305093

See Accessories on page 14 and Trolleys on page 15

Chemotherapy Sharps Collectors

BD Chemotherapy Collectors feature a variety of lid designs and sizes for disposal of trace chemotherapy waste (see Waste Stream Segregation on page 4). Trace waste is defined by CDC / NIOSH as "those that contain less than 3% by weight of the original quantity of hazardous drugs." *

* CDC/NIOSH Alert Publication No. 2004-165

Key product features include:



- Gaskets or liquid-absorbing materials to prevent leakage
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Temporary closure feature and handles designed for safe transportation
- Large and extra-large collectors are fully compatible with BD Trolleys (see page 15)
- Clear lids with multiple opening styles to accommodate different types of waste and environments
 - Dual-opening hinge top designed so that one side can remain open and the other locked based on the size of access opening desired



305613



305076



Chemotherapy

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305076	3 gal Hinge Cap	12 x 10.5 x 7.5	12	305495, 305497
305493	5 gal Hinge Cap	18 x 10.5 x 7.5	8	
305578	5 gal Screw Cap One Piece	18 x 10.5 x 7.5	8	
305114	5 gal with Tethered Cap	18 x 10.5 x 7.5	8	
305603	9 gal Hinge Top with Gasket	18.5 x 17.75 x 11.75	8	305090
305604	9 gal Slide Top with Gasket	18.5 x 17.75 x 11.75	8	305090, 305091
305614	17 gal Hinge Top with Gasket	23.5 x 20 x 14.75	5	305092
305613	19 gal Slide Top with Gasket	26.25 x 20 x 14.75	5	305092, 305093

See Accessories on page 14 and Trolleys on page 15

RCRA Hazardous Waste Collectors — DOT

BD RCRA Collectors are intended for Resource Conservation and Recovery Act (RCRA) solid hazardous waste and trace liquids associated with these wastes (see Waste Stream Segregation on page 4).


BD offers two versions of RCRA collectors, **DOT PGII / UN-rated** and **non-DOT**, to meet a variety of needs. DOT / UN-rated RCRA collectors are designed to be transported without being overpacked by the waste haulers. As most waste haulers will overpack smaller collectors whether they are DOT rated or not, BD offers these in 13 and 16 gallon sizes.

What is UN?

The UN or *United Nations Recommendations for the Transport of Dangerous Goods Model Regulations* are suggested regulations for controlling the transport of hazardous materials, including medical waste. The regulations define Packing Groups and test requirements to qualify for a UN-compliant package marking.

A UN-compliant package must have a marking that describes its maximum content weight and the Packing Group that it passed.

It looks like:



The regulations describe packaging, labeling and documentation required for transporting hazardous materials, contain a table with codes for each specific type of hazardous material and have several codes for pharmaceutical waste.

What is DOT?

The DOT, or *U.S. Department of Transportation Regulations*, is a government agency regulating fast, safe and efficient transportation of dangerous goods.

What is PGII?

PG (Packing Group) is a Department of Transportation term that groups hazardous materials according to the danger presented by the materials. There are three packing group classifications.

I - great II - medium III - minor

Key product features of RCRA Hazardous Waste Collectors include:

- Black base with black tops to distinguish the waste stream
- Leak resistant gasketed tops with 7" port opening



305074

305073

Leak Resistant Gasketed Top with 7" port opening



RCRA Hazardous Waste Collectors DOT PGII / UN-rated

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case
305073	13 gal DOT Port Top	16 x 21.5 x 11	6
305074	16 gal DOT Port Top	16 x 24.75 x 11	6

RCRA Hazardous Waste Collectors — Non-DOT

Non-DOT Collectors feature a variety of lid designs and sizes and come with a liquid-absorbent material inside the collector for trace waste.*

* CDC/NIOSH Alert Publication No. 2004-165



305069





305066



305075

Key product features of non-DOT Collectors include:

- Made from a minimum of 70% recycled plastic
- Clear tops with multiple opening styles to accommodate different types of waste and environments
 - Dual-opening hinge top designed so that one side can remain open and the other locked based on department needs
 - Dual-lock cap (3 and 5 gallon) designed to provide temporary closure for safe intermediate transportation and secure final closure with two locks
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transportation
- Temporary closure feature
- Large and extra-large collectors are fully compatible with BD Trolleys (see page 15)



All RCRA collectors must be incinerated for final disposal.

RCRA Hazardous Waste Collectors Non-DOT

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305075	3 gal BD Recykleen™ Hinge Cap with Petals	12 x 10.5 x 7.5	12	305495, 305497
305068	9 gal BD Recykleen™ Hinge Top with Gasket	18.5 x 17.75 x 11.75	8	305495, 305497
305069	9 gal BD Recykleen™ Slide Top with Gasket	18.5 x 17.75 x 11.75	8	305495, 305497
305071	17 gal BD Recykleen™ Hinge Top with Gasket	23.5 x 20 x 14.75	8	305090
305072	19 gal BD Recykleen™ Slide Top with Gasket	18.5 x 17.75 x 11.75	8	305090, 305091
305066	3 gal Hinge Cap	12 x 10.5 x 7.5	12	305495, 305497
305067	5 gal Hinge Cap	18 x 10.5 x 7.5	5	305092, 305093

See Accessories on page 14 and Trolleys on page 15

Cabinets and Accessories

BD Recykleen™ Cabinets are available in two sizes: small for use with 5.4 quart collectors and large for use with 2 and 3 gallon collectors. Each size is available in both green and brown semi-transparent colors.



Key product features include:

- A universal lock providing security and preventing unwanted access to collectors
- A one-piece transparent door facilitating easy cleaning
- The ability to clearly see the collector fill-level indication
- A two-year warranty providing product replacement in case of manufacturer defects



BD Recykleen™ Cabinets

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305096	5.4 qt BD Recykleen™ Sharps Collector Wall Cabinet Brown	14 x 13.3 x 5.5	1	305056, 305057, 305517, 305551, 305543, 305426, 305427, 305428, 305448
305098	5.4 qt BD Recykleen™ Sharps Collector Wall Cabinet Green	14 x 13.3 x 5.5	1	
305097	2 and 3 gal BD Recykleen™ Sharps Collector Wall Cabinet Brown	19 x 14.5 x 7.5	1	305053, 305054, 305422, 305435, 305423, 305436, 305448
305099	2 and 3 gal BD Recykleen™ Sharps Collector Wall Cabinet Green	19 x 14.5 x 7.5	1	

Brackets

Brackets are designed to secure the collector to a wall or cart. Some brackets include additional security via a locking mechanism.

Stabilizers

Stabilizers are designed to stabilize collectors on a countertop or other flat surface.

Brackets and Stabilizers

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305447	Locking Bracket	9.5 x 9.5 x 4.25	10	All 5.4 qt
305409	Locking Bracket	4.25 x 14 x 4.75	10	2 and 3 gal Counterbalanced Door
305588	Wire Bracket	4.5 x 12 x 5.3	6	2 and 3 gal Counterbalanced Door
305021	Stabilizer	3 x 13 x 6.5	5	All 5.4 qt
305552	Stabilizer	2.75 x 13.5 x 7.3	5	2 and 3 gal Counterbalanced Door
305448	Glove Box Cabinet Bracket	6 x 12 x 4.4	6	Attaches to Cabinets and 5.4 qt Locking Bracket
305495	Wire Security Bracket	6 x 3.5 x 1.2	12	Multi-use One Piece Collectors
305497	Quick Release Strap	1.5 x 8 x 1	12	Multi-use One Piece Collectors

Trolleys

Our full line of BD Recykleen™ Trolleys are available in basic and foot-operated versions for use with large and extra-large collectors. Trolleys provide stability and mobility for large and extra-large collectors.



Key product features include:

- A design that helps minimize the risk of collectors tipping over
- Wheels that enable the collector to be conveniently placed where it will be most used
- Locking wheels, ensuring the product is mobile only when you want it to be
- The ability to open and close the collector for hands-free disposal (foot-operated trolleys only)
- A two-year warranty providing product replacement in case of manufacturer defects



BD Recykleen™ Trolleys

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305090	9 gal BD Recykleen™ Basic Trolley	23.5 x 16.3 x 28.6	1	All 9 gal
305091	9 gal BD Recykleen™ Foot-operated Trolley	23.5 x 16.3 x 28.6	1	All 9 gal Slide Top
305092	17 and 19 gal BD Recykleen™ Basic Trolley	24.1 x 18 x 27.5	1	All 17 and 19 gal
305093	19 gal BD Recykleen™ Foot-operated Trolley	24.3 x 17.7 x 35.2	1	All 19 gal Slide Top

BD Environmental Stewardship

BD's approach to corporate citizenship is embedded in our culture, and is intrinsic to our purpose, *Helping all people live healthy lives.* Responsible corporate citizenship is integrally tied to our global business strategy. Sustained performance and a dedication to improvement are at BD's core. Recognizing that human health and a healthy environment are inseparable, BD is committed to responsible environmental stewardship in our business conduct.




BD Disposal Solutions: A Commitment to Quality

BD Medical Waste Collectors are manufactured using polypropylene or polyethylene resins and colorants that are nontoxic and designed to comply with the Toxics in Packaging Clearing House (TPCH) heavy metal limits. The collectors are safe for incineration. Under combustion with adequate access of air, the resin burns to water vapor, carbon dioxide and carbon ash.

BD Sharps Collectors meet or exceed the puncture force requirements defined by the American Society for Testing and Materials (ASTM) and British Standards. Sharps collectors are classified as medical devices and are manufactured according to the FDA Quality System Regulation and the International Quality Systems Standard ISO 9001.

BD Sharps Collectors are listed with the FDA and meet pre-market registration requirements in accordance with 21 CFR 807. BD Sharps Collectors also meet the Occupational Safety and Health Administration (OSHA) Bloodborne Pathogens Standards 29 CFR 1910.1030.



Contact Customer Service at **1-888-237-2762**
or your BD Sales Representative for additional information.



BD Medical
1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/sharps