



The BD Epilor™ Plastic Loss of Resistance (LOR) Syringe

Designed to help achieve greater success in detecting the epidural space.

- **Optimized stopper design and stopper material** aids in improving sensitivity and reduces the chances of air or saline from leaking around the stopper which can lead to false detection of the loss of resistance (LOR).
- **Stabilized plunger rod design and a positive plunger rod stop** provide stability and prevent “plunger wobble” that could cause leaking around the stopper and hinder detecting loss of resistance.
- **Other important features include:** Latex Free construction, enhanced graphics and a clear barrel, textured plunger rod and thumb pad resulting in an ergonomic design for ease of use.
- The BD Epilor plastic LOR syringe is available in a convenient single 7ml size.

New Catalog Number	Description	Packaging
405291	7ml BD Epilor™ Plastic LOR syringe with Luer-Lok™ Tip	10 per Box 50 per Case
405292	7ml BD Epilor™ Plastic LOR syringe with Luer-Slip Tip	10 per Box 50 per Case

To order contact the BD Customer Service Center at 1.888.237.2762. The BD Epilor Plastic LOR syringe is also available in Baxter/BD Alliance regional anesthesia standard and custom trays. See *reverse side*.



BD Medical Systems
1 Becton Drive - MC211
Franklin Lakes, NJ 07417-1884

BAXTER & BD REGIONAL ANESTHESIA TRAYS WITH BD EPILOR™ PLASTIC LOSS OF RESISTANCE SYRINGE

Teaming up to offer you choice in your Epidural tray components, Baxter and BD are pleased to provide you with the BD Epilor Plastic Loss of Resistance Syringe. The BD Epilor Plastic Loss of Resistance Syringe is designed:

- To enhance the detection of loss of resistance.
- With non-latex components, enhanced graphics, a clear barrel, textured thumb pad and an ergonomic design for greater ease of use.



The BD Epilor Plastic Loss of Resistance Syringe is available in the following Baxter/BD Epidural Configurations:

Code Number	Procedural Epidural Needle	20ga Nylon Epidural Catheter, closed-tip	19ga Nylon Epidural Catheter, closed-tip	BD Epilor plastic LOR syringe, 7ml, LS	3ml plastic syringe	20ml plastic syringe	19ga x 1-1/2" Blunt aspirating filter needle	22ga x 1-1/2" hypodermic needle	25ga x 1" hypodermic needle	1% Lidocaine HCL, 5ml	0.9% Sodium Chloride Injection, USP, 10ml	1.5% Lidocaine HCL Injection w/Epidephrine 1:200,000, 5ml (for test dose)	0.22 micron filter syringe adaptor/ protector set	Fenestrated Drape (C=Clear, P=Paper)
1T5683	18ga x 3-1/2" clear-hub Weiss	X		X	X	X	X	X	X	X	X	X	X	P
1T5672	17ga x 3-1/2" clear-hub Weiss		X	X	X	X	X	X	X	X	X	X	X	P
1T5684	18ga x 3-1/2" clear-hub Weiss	X		X	X	X	X	X	X	X	X	X	X	C
1T5673	17ga x 3-1/2" clear-hub Weiss		X	X	X	X	X	X	X	X	X	X	X	C

Notes: LL denotes BD Luer-Lok™ connection
LS denotes a Luer-Slip™ connection

To Order a Baxter/BD Regional Anesthesia Tray with the BD Epilor Plastic Loss of Resistance syringe, please contact the Baxter Medication Delivery Center for Service at 1.888.229.0001

**Baxter and BD offer a full range of Regional Anesthesia
Trays and choice quality components you expect.**

Baxter

Baxter Healthcare Corporation
Medication Delivery
Deerfield, IL 60015

BD

BD Medical Systems
1 Becton Drive - MC211
Franklin Lakes, NJ 07417-1884

In Canada: **BD Medical Systems**
2771 Bristol Circle
Oakville, Ontario, L6H 6R5

www.bd.com/surgical

eGeneralMedical+com

The Best Equipment & Supplies for Students and Healthcare Professionals



eGeneralMedical.com is a company owned and operated by a practicing Physician Assistant. Specializing in medical, health and wellness products for healthcare professionals and individuals, we are a growing distributor of medical equipment, supplies, disposables, pharmaceuticals, medical apparel, and supplements.



Stethoscopes
& Accessories



Diagnostic
Sets



Logo Embroidered
Labcoats



P.A. Study Guides
& Review Books



PDA Bundles
specifically designed for
the practicing P.A.

At eGeneralMedical.com, we sell only high quality product line that caters to physician patients and to the consumer. We are located in Raleigh, North Carolina, the largest city in the Research Triangle Region, one of the leading areas for medical and pharmaceutical research in the world.

We also provide services and incentives for healthcare professionals, education centers and community organizations including: PIPEC (Purchasing Incentive Program for Education & Community), Medical Rewards, Honorarium Redemptions, Coupon Order Discounts. In addition to our newsletters, eGeneralMedical is further developing a resource center which delivers healthcare information, wellness tips and the latest medical news to professionals and consumers.

Top Quality Health Products at Low Prices!