

# Catheter placement kits and accessories

## PleurX catheter and starter kit 50-7700 and 50-9900

The PleurX catheter and starter kit contains the items needed for catheter placement as well as discharge vacuum bottles with patient education and training materials to help ensure a smooth transition to a home or alternative care center. The PleurX catheter kit contains the components needed for catheter placement and initial drainage. The starter kit ensures patients have the drainage supplies they need upon discharge and includes training materials to help patients and caregivers confidently drain at home.

### PleurX pleural (50-7700) or peritoneal (50-9900) kits contain:

#### 1. Preparation components

- ChloraPrep® 10.5 mL applicator with Hi-Lite Orange® tint (*qty. 2*)
- Fenestrated drape
- Filter Straw filter
- Lidocaine HCL USP, 1%, 5 mL (*qty. 3 ampoules*)
- SafetyGlide™ needle, 22 G x 1½" (*3.8 cm*)
- SafetyGlide needle, 25 G x 1" (*2.5 cm*)
- Syringe, 10 mL

#### 2. Placement components

- PleurX catheter (*15.5 Fr*)
  - 66 cm (*pleural*)
  - 71 cm (*peritoneal*)
- Safety scalpel
- Guidewire introducer with needle (*.038"*)
- Syringe, 10 mL
- J-tip guidewire
  - 45 cm (*pleural*)
  - 60 cm (*peritoneal*)
- Forceps
- Dilator, 8 Fr
- Dilator, 12 Fr
- Peel-away introducer, 16 Fr
  - 13.7 cm dilator length, 10 cm sheath length (*pleural*)
  - 17.8 cm dilator length, 14 cm sheath length (*peritoneal*)
- Tunneler
- Point-Lok® sharps safety device

#### 3. Closing components

- Tweezer
- Silk suture, 2-0 straight needle
- Silk suture, 3-0 curved needle

#### 4. Drainage components


- Lockable drainage line with access tip
- Needle, 17 G x 1" (*2.5 cm*)
- 5-in-1 drainage line adapter
- Valve cap
- Alcohol pad

#### 5. Dressing components

- Gauze pads, 4" x 4" (*10.2 cm x 10.2 cm, qty. 6*)
- Foam catheter pad
- Self-adhesive dressing



PleurX catheter kits are also available separately with cat. no. 50-7000B (*pleural*) and 50-9000B (*peritoneal*).



PleurX catheter (15.5 Fr)

### Starter kit components

- PleurX drainage kits, including 1,000 mL bottles and procedure packs (qty. 4)
- Convenient carrying box
- Patient information kit

### Patient information kit components

- Patient education DVD
- Introductory letter with ordering information
- Frequently asked questions (FAQ) brochure
- Instructions for Use booklet
- Reference wall chart
- Emergency information card
- Drainage log



Starter kits are also available separately with 500 mL bottles (cat. no. 50-0075) and 1,000 mL bottles (cat. no. 50-0071).