

Drainage kits and accessories

PleurX drainage kit 50-7500B (500 mL) and 50-7510 (1,000 mL)

The PleurX drainage kit includes a procedure pack that provides the components needed each time the patient drains fluid.

Each drainage kit (50-7500B and 50-7510) contains:

- Plastic vacuum bottle with attached drainage line

Procedure pack components

- Blue wrapping
- Alcohol pads (qty. 3)
- Pair of gloves
- Valve cap
- Blue emergency slide clamp
- Gauze pads, 4" x 4" (10.2 cm x 10.2 cm, qty. 4)
- Foam catheter pad
- Self-adhesive dressing

